

Check out the great selection and great pricing on electrosurgical cables compatible with a wide range of generator / instrument combinations.

Disposable

Sold Box/10. Sterile. Single-Use. 10ft length. Off-white.



Monopolar: NAQD1001*



NAQD1004*



NAQD2201*



Bipolar: NAQD2103*



NAQD2407*

- NAQD1001 Monopolar Bovie® Plug (male) to Monopolar Pin (female, straight)
- NAQD1004 Monopolar Bovie® Plug (male) to Monopolar Pine (female, 90 degree)
- NAQD2201 Monopolar Wolf® Generator (female) to Wolf Kleppinger (female)
- NAQD2013 Bipolar Valleylab® Banana Pin (male) to Storz/Stryker Kleppinger (female)
- NAQD2407 Bipolar Valleylab® Banana Pin (male) to Twin Pin (female, 45 degree)

* Style available in reusable or disposable

**Don't see
what you
need here?
Call us!**

MONOPOLAR + BIPOLAR REUSABLE + DISPOSABLE

Capital Medical Resources offers electrosurgical cables for generators including: Bovie®, Wolf®, Valleylab®, Banana plug, Twin pin, and more, to connect to variety of instrument connectors such as Stryker®, Storz®, Wolf®, Aesculap®, Conmed®, etc. Let us cross reference what you need!

Reusable:

Sold Each. Re-usable. Non-Sterile. 10ft Length.



Monopolar: NAQI1001*



NAQI1003



NAQI1004*



NAQW2201*



Bipolar: NAQI12301*



NAQI12401



NAQI2402



NAQI2407*



NAQI2900



NAQS2101



NAQS2103



NAQW2204

- NAQI1001 Monopolar Bovie® Plug (male) to Monopolar Pin (female, straight)
- NAQI1003 Monopolar Valleylab® Banana Pin (male) to Valleylab® Banana Pin (male)
- NAQI1004 Monopolar Bovie® Plug (male) to Monopolar Pin (female, 90 Degree)
- NAQW2201 Monopolar Wolf® Generator (female) to Wolf® Kleppinger (female)
- NAQI2301 Bipolar Valleylab® Banana Pin (male) to Wolf® Kleppinger (female)
- NAQI2401 Bipolar Valleylab® Banana Pin (male) to Twin Pin (female, Straight)
- NAQI2402 Bipolar Fischer® Connector (OMS type) to Twin Pin (female, Straight)
- NAQI2407 Bipolar Valleylab® Banana Pin (male) to Twin Pin (female, Straight)
- NAQI2900 Bipolar Valleylab® Banana Pin (male) to Metz/Mayo Scissors (female)
- NAQS2101 Bipolar Valleylab® Banana Pin (male) to Storz® Fischer Connector
- NAQS2103 Bipolar Valleylab® Banana Pin (male) to Storz® Kleppinger (female)
- NAQW2204 Bipolar Wolf® Generator (female) to Twin Pin (female, Straight)
- NAQW2202 Tripolar Valleylab® Banana Pin (female) to Wolf® Generator (female)

© All Registered Trademarks are owned by their respective companies

Tripolar: NAQW2204

