

ELECTROSURGICAL CABLES

Check out the great selection and great pricing on electrosurgical cables compatible with a wide range of generator / instrument combinations.



Disposable Monopolar



Reusable Bipolar

**Don't see
what you
need here?
Call us!**

Monopolar	Disposable	600-100	Box/10	Bovie® Plug (male) to Monopolar Pin (female, 90°), Single Use
Monopolar	Disposable	600-101	Box/10	Bovie® Plug (male) to Monopolar Pine (straight), Single Use
Monopolar	Reusable	600-200	Each	Bovie® Plug (male) to Monopolar Pin (female 90°)
Monopolar	Reusable	600-201	Each	Bovie® Plug (male) to Monopolar Pine (straight)
Bipolar	Reusable	600-303	Each	Banana Pins (male) to Bladed Connector (female), Silicone
Bipolar	Reusable	600-304	Each	Banana Pins (male) to Bipolar (female), Silicone
Bipolar	Reusable	600-305	Each	Banana Pins (male) to Bipolar (female), Thermoplastic Rubber

MONOPOLAR + BIPOLAR REUSABLE + DISPOSABLE

Capital Medical Resources offers electrosurgical cables for generators including: Bovie®, Wolf®, Valleylab®, Banana plug, Twin pin, and more, to connect to variety of instrument connectors such as Stryker®, Storz®, Wolf®, Aesculap®, Conmed®,