

# CASE CART STATUS INDICATORS

The ideal solution for easy and affordable identification and communication of the status of case carts and other medical equipment.

## STATUS INDICATOR DIAL


The 5-position locking dial can be rotated to identify status including: **CLEAN** | **BIO-HAZARD** | **COMPLETE** | **INCOMPLETE** | **BLANK DRY ERASE** for custom identification. Includes mounting tape for installation. An optional stainless steel mounting plate is available for permanent welded installation on metal surfaces.

## STATUS INDICATOR SLIDERS

The 2-position locking plate slides into position to indicate **CLEAN** | **BIOHAZARD** status. Sliders come with mounting tape for quick and easy installation on a variety of surfaces and equipment. Regular and mini size available!

### Ordering info:

- 300700 8" Rotating Dial
- 300701 SS Mounting Plate (optional)
- 300702 3" x 6-1/2" Slider
- 300703 1" x 3" Mini Slider



Status indicators are designed and tested to withstand use on cart washers and in autoclaves at 250 degrees for 30-minute cycles without disfigurement or damage. 2-sided mounting tape is also durable in these conditions and can be easily replaced if damaged or for repositioning.

## IDENTIFICATION + **BETTER** COMMUNICATION

*Sterile Processing Professionals can save time, increase efficiency, and deliver safe and effective equipment by using status indicators to improve communication.*