

MOBILE & MEDICAL MOUNTING SOLUTIONS FOR EVERY APPLICATION

- + SUPPORTARMS & FIXINGS
- + MONITOR & DISPLAY MOUNTS
- + LIGHT & HEAVY DUTY CARTS
- + WALL CHANNEL & COLUMN SYSTEMS
- + WORKSTATIONS



MOBILE SOLUTIONS FOR EVERY APPLICATION

CART LIGHT DUTY

A roll stand for light applications of up to 14 kg/ 30.9 lbs. All carts are tested and certified as safe before supply. The robust construction makes it a favourite amount nursing staff. Smooth manoeuvrability, ease of cleaning combined with a 2-year warranty make this roll stand out.

DESIGNED FOR


- + Daily repetitive use
- + Reliability
- + Current regulation (EU) 2017/745 (MDR) & IEC 60601-1 safety standards

FEATURES

- + Easy to clean
- + No pooling areas for liquids
- + Easy-glide movement
- + Lockable castors for stability
- + Small footprint

ACCEPTS

- + All standard fixings
- + External power supply holsters
- + Baskets
- + Many other accessories




MONITOR PRESENTATIONS

Which monitor would you like to mount?
CIM med offers solutions for a wide range of medical device manufacturers.



ACCESSORIES

A huge variety of accessories is available which you will find in our accessories catalogue.
Please note that not every accessory is suitable for every cart.



MOBILE SOLUTIONS FOR EVERY APPLICATION

CART HEAVY DUTY

The Heavy Duty Cart can be configured according to your needs. An extensive range of accessories is available to match your requirements. All individual solutions are tested and certified before delivery and are conform with EN 60601-1, 3rd Edition.

DESIGNED TO

- + Meet your requirements**
Literally hundreds of configurations available.
- + Save space**
Minimizing required space while maximizing access to your equipment.
- + Promote highest standards**
Regulation (EU) 2017/745 (MDR) & IEC 60601-1

FEATURES

- + Modularity**
Allows individual configurations realized efficiently and cost-effectively.
- + Variability**
The variable concept permits upgrading at any time.
- + Functionality**
Maximum modularity and ideal ergonomic design combined with refined functionality.

ACCEPTS

- + Drawers and shelves**
Single and double drawer consoles or simply shelves to place your equipment.
- + Support arms**
CIM height adjustable and non-height adjustable arms to mount monitors, devices and other medical equipment.
- + Accessories**
A huge variety of accessories is available: cable hooks, baskets, IV poles.




MONITOR PRESENTATIONS

Which monitor would you like to mount?
CIM med offers solutions for a wide range of medical device manufacturers.



ACCESSORIES

A huge variety of accessories is available which you will find in our accessories catalogue.
Please note that not every accessory is suitable for every cart.





Flush mount



Support arm pivot



Extended pivot arm light duty



Extended pivot arm heavy duty

Monitor presentations

- VESA 75/100 monitor adaptation
- Quick release VESA
- 5" plunge plate
- Manufacturer-specific mounting interface
- Keyboard tray
- Ergonomic W-bracket
- Fold-up Z-bracket

Fixings



4V

Adapter for vertical wall channel



2F

Fixed rail clamp



3

Vertical rail clamp



3DR
3TX

Dual vertical rail clamp



1
1B
1F
1C
1V

Fixed / universal pole clamp



13

Adapter for Mindray medical supply units



17

Adapter for Amico pendants & headwalls



9

Adapter for Löwenstein Medical anaesthesia machines



4K

Adapter for anaesthesia machines with existing side channel



Height adjustable arm



Extended height adjustable arm

Monitor presentations

- VESA 75/100 monitor adaptation
- Quick release VESA
- 5" plunge plate
- Manufacturer-specific mounting interface
- Keyboard tray
- Ergonomic W-bracket
- Fold-up Z-bracket

Fixings

 <p>4V</p> <p>Adapter for vertical wall channel</p>	 <p>2F</p> <p>Fixed rail clamp</p>	 <p>3</p> <p>Vertical rail clamp</p>	 <p>3DR 3TX</p> <p>Dual vertical rail clamp</p>	 <p>1 1B 1F 1C 1V</p> <p>Fixed / universal pole clamp</p>
 <p>13</p> <p>Adapter for Mindray medical supply units</p>	 <p>17</p> <p>Adapter for Amico pendants & headwalls</p>	 <p>9</p> <p>Adapter for Löwenstein Medical anaesthesia machines</p>	 <p>4K</p> <p>Adapter for anaesthesia machines with existing side channel</p>	



Available in pivot arm*

Workstations

CIM Cable Integrated Mounting System



Height adjustable arm for flat panel and keyboard / work surface



Extended height adjustable arm for flat panel and keyboard / work surface

Fixings



4V

Adapter for vertical wall channel



2F

Fixed rail clamp



3

Vertical rail clamp



3DR
3TX

Dual vertical rail clamp



1
1B
1F
1C
1V

Fixed / universal pole clamp



13

Adapter for Mindray medical supply units



17

Adapter for Amico pendants & headwalls



9

Adapter for Löwenstein Medical anaesthesia machines



4K

Adapter for anaesthesia machines with existing side channel

Combi systems with fold-up Z-bracket

CIM Cable Integrated Mounting System



Height adjustable arm with fold-up Z-bracket



Extended height adjustable arm with fold-up Z-bracket

Fixings



4V

Adapter for vertical wall channel



2F

Fixed rail clamp



3

Vertical rail clamp



3DR
3TX

Dual vertical rail clamp



1
1B
1F
1C
1V

Fixed / universal pole clamp



13

Adapter for Mindray medical supply units



17

Adapter for Amico pendants & headwalls



9

Adapter for Löwenstein Medical anaesthesia machines



4K

Adapter for anaesthesia machines with existing side channel