



## SARS-CoV-2 Antibody Seroconversion Panel

### MRNDx PANEL ID: MRNCOV-501, Donor ID: 776777

**All Plasma was collected from a 20 year old Hispanic female donor. The donor tested PCR positive for SARS-CoV-2 on 11/23/2020.**

Symptoms of Fever, Cough, Tiredness, Night Sweats, Loss of Smell, Loss of Taste, and Vomiting. Patient was not Hospitalized.

Panel ID	Abbott Architect SARS-CoV-2 IgM (Cut-off $\geq$ 1.0)	Abbott Architect SARS-CoV-2 IgG (Cut-off $\geq$ 1.4)	Chembio DPP SARS-CoV-2 IgM/IgG Rapid IgM (Cut-off $\geq$ 24)	Chembio DPP SARS-CoV-2 IgM/IgG Rapid IgG (Cut-off $\geq$ 24)	Ortho Eci SARS-CoV-2 IgG (R&D) (Cut-off $\geq$ 1.0)	DATE OF LAST COVID SYMPTOMS	DAYS SINCE LAST SYMPTOMS	DOD
MRNCOV-501-01	0.08	0.01	2	2	0.01	11/23/2020	2	11/25/2020
MRNCOV-501-02	0.19	0.02	2	2	0.01	11/23/2020	5	11/28/2020
MRNCOV-501-03	5.21	2.54	76	22	0.59	11/23/2020	10	12/03/2020
MRNCOV-501-04	10.32	6.49	132	171	4.95	11/23/2020	14	12/07/2020
MRNCOV-501-05	12.6	7.34	166	155	6.83	11/23/2020	18	12/11/2020
MRNCOV-501-06	12	7.45	181	185	8.49	11/23/2020	22	12/15/2020
MRNCOV-501-07	12.09	7.37	171	223	9.25	11/23/2020	26	12/19/2020
MRNCOV-501-08	10.45	6.99	156	255	9.92	11/23/2020	30	12/23/2020
MRNCOV-501-09	10.76	5.45	149	286	9.48	11/23/2020	33	12/26/2020
MRNCOV-501-10	10.35	7.19	116	213	9.6	11/23/2020	37	12/30/2020
MRNCOV-501-11	6.4	5.85	118	256	0.98	11/23/2020	46	01/08/2021

\*Abbott Architect SARS-CoV-2 IgM - Kit Lot#25047FN00 - EXP:06/06/2021, Abbott Architect SARS-CoV-2 IgG - Kit Lot#25159FN00 - EXP:06/29/2021\*

\*CE - Chembio DPP SARS-CoV-2 IgM/IgG Rapid Test - Kit Lot#2CV120420/A - EXP:10/28/2022, Ortho Eci SARS-CoV-2 IgG (R&D) - Kit Lot#0335 - EXP:06/06/2021\*

**All plasma was collected at a U.S. Licensed Plasmapheresis Center.**

#### CAUTION: HUMAN BLOOD

All materials listed above are considered to be potentially infectious and should be handled with universal precautions.

FOR RESEARCH USE ONLY. Not for use in Diagnostic Procedures.