



## SARS-CoV-2 Antibody Seroconversion Panel

### MRN PANEL ID: MRNCOV-505, Donor ID: 775154

**All Plasma was collected from a 45 year old Hispanic male donor. The donor tested PCR positive for SARS-CoV-2 on 11/06/2020.**

Symptoms of Loss of Smell, Loss of Taste, and Headache. Patient was not Hospitalized.

Panel ID	Abbott Architect SARS-CoV-2 IgM (Cut-off $\geq$ 1.0)	Abbott Architect SARS-CoV-2 IgG (Cut-off $\geq$ 1.4)	Chembio DPP SARS-CoV-2 IgM/IgG Rapid IgM (Cut-off $\geq$ 24)	Chembio DPP SARS-CoV-2 IgM/IgG Rapid IgG (Cut-off $\geq$ 24)	Ortho Eci SARS-CoV-2 IgG (R&D) (Cut-off $\geq$ 1.0)	DATE OF LAST COVID SYMPTOMS	DAYS SINCE LAST SYMPTOMS	DOD
MRNCOV-505-01	3.39	2.94	44	22	0.05	ACTIVE	ACTIVE	11/11/2020
MRNCOV-505-02	6.43	4.49	90	52	0.83	ACTIVE	ACTIVE	11/17/2020
MRNCOV-505-03	6.11	5.24	124	73	2.97	ACTIVE	ACTIVE	11/26/2020
MRNCOV-505-04	5.65	5.15	121	95	4.26	ACTIVE	ACTIVE	11/30/2020

\*Abbott Architect SARS-CoV-2 IgM - Kit Lot#25047FN00 - EXP:06/06/2021, Abbott Architect SARS-CoV-2 IgG - Kit Lot#25159FN00 - EXP:06/29/2021\*

\*CE - Chembio DPP SARS-CoV-2 IgM/IgG Rapid Test - Kit Lot#2CV120420/A - EXP:10/28/2022, Ortho Eci SARS-CoV-2 IgG (R&D) - Kit Lot#0335 - EXP:06/06/2021\*

**All plasma was collected at a U.S. Licensed Plasmapheresis Center.**

**The Plasma was tested and found Negative for Anti-HCV and Anti-HIV 1/2 by FDA approved tests.**

#### CAUTION: HUMAN BLOOD

All materials listed above are considered to be potentially infectious and should be handled with universal precautions.

FOR RESEARCH USE ONLY. Not for use in Diagnostic Procedures.