



## SARS-CoV-2 Antibody Seroconversion Panel

### MRNDx PANEL ID: MRNCOV-508, Donor ID: 776963

**All Plasma was collected from a 37 year old Hispanic female donor. The donor tested PCR positive for SARS-CoV-2 on 11/11/2020.**  
Symptoms of Cough, Tiredness, Diarrhea, Body Aches, Chills, Runny Nose, Loss of Smell, Sore Throat, Loss of Taste, and Headache. Patient was not Hospitalized.

Panel ID	Abbott Architect SARS-CoV-2 IgM (Cut-off ≥ 1.0)	Abbott Architect SARS-CoV-2 IgG (Cut-off ≥ 1.4)	Chembio DPP SARS-CoV-2 IgM/IgG Rapid IgM (Cut-off ≥ 24)	Chembio DPP SARS-CoV-2 IgM/IgG Rapid IgG (Cut-off ≥ 24)	Ortho Eci SARS-CoV-2 IgG (R&D) (Cut-off ≥ 1.0)	DATE OF LAST COVID SYMPTOMS	DAYS SINCE LAST SYMPTOMS	DOD
MRNCOV-508-01	0.34	<b>3.39</b>	9	23	0.99	11/30/2020	9	12/09/2020
MRNCOV-508-02	0.27	<b>4.11</b>	15	<b>69</b>	<b>1.97</b>	11/30/2020	12	12/12/2020
MRNCOV-508-03	0.23	<b>4.41</b>	8	<b>45</b>	<b>2.7</b>	11/30/2020	15	12/15/2020
MRNCOV-508-04	0.2	<b>4.15</b>	1	10	<b>2.84</b>	11/30/2020	19	12/19/2020
MRNCOV-508-05	0.15	<b>3.72</b>	9	23	<b>2.5</b>	11/30/2020	22	12/22/2020
MRNCOV-508-06	0.15	<b>2.93</b>	7	17	<b>1.75</b>	11/30/2020	26	12/26/2020
MRNCOV-508-07	0.12	<b>2.78</b>	5	<b>61</b>	<b>1.93</b>	11/30/2020	18	12/18/2020
MRNCOV-508-08	0.19	<b>1.74</b>	6	19	<b>2.14</b>	11/30/2020	78	02/16/2021

\*Abbott Architect SARS-CoV-2 IgM - Kit Lot#25047FN00 - EXP:06/06/2021, Abbott Architect SARS-CoV-2 IgG - Kit Lot#25159FN00 - EXP:06/29/2021\*

\*CE - Chembio DPP SARS-CoV-2 IgM/IgG Rapid Test - Kit Lot#2CV120420/A - EXP:10/28/2022, Ortho Eci SARS-CoV-2 IgG (R&D) - Kit Lot#0335 - EXP:06/06/2021\*

**All plasma was collected at a U.S. Licensed Plasmapheresis Center.**  
**The Plasma was tested and found Negative for Anti-HCV and Anti-HIV 1/2 by FDA approved tests.**

#### CAUTION: HUMAN BLOOD

All materials listed above are considered to be potentially infectious and should be handled with universal precautions.  
**FOR RESEARCH USE ONLY. Not for use in Diagnostic Procedures.**