

# THERMA ZONE<sup>®</sup>

THERMAL THERAPY

for  
**sports medicine**

As a professional, you should leverage the latest technology to treat patients and maintain a safe and efficient training room. Unfortunately, traditional equipment does not meet your needs.

Ice machines must run continuously and require heavy maintenance. Hydroculators can take up to 3 hours to heat.

**ThermaZone**, on the other hand, does not require ice, operates on two ounces of distilled water, and is virtually maintenance-free.

**ThermaZone** has a broad temperature range from 34°F – 125°F, 20 comfort levels, and 5 pre-set timer modes.

**ThermaZone** is the convenient, consistent, heating and cooling therapy you and your patients deserve.

Training room equipment does not need to be big to make an impact.



## ThermaZone<sup>®</sup>

iceless. programmable. opioid-free.

**ThermaZone** is a powerful and compact device that uses thermoelectric technology to provide heating and cooling therapy where and when you need it, without the use of ice.

**ThermaZone** satisfies a clear unmet clinical need – targeted, localized, non-opioid therapy with consistent temperature management.



FEATURING



Temperature controlled water flowing through welded microchannels

Our pads feature **Motion Advantage™** technology, which leverages precise fluid movement to deliver uniform temperature distribution across the surface of the pad for maximum coverage.

ON/OFF control panel with mode and temperature control buttons

5 preset timer options

Lighted comfort level window display


Lightweight, easy-to-carry handle



Adapter plug for AC/DC power supply or auto power cord

Connector ports for easy connection to pads

Unit Size: 11" L x 4.5" W x 4.5" H  
Unit Weight: 3 lbs. 7 oz.




## Compelling Features

- Non-opioid therapy
- Heat or cool therapy in one device
- 3 minutes to set up requires only 2 oz. of distilled water
- Rapid achievement of temperature – 20 different comfort levels
- 5 preset timer options:
  - 10, 20, or 30 minutes
  - 30 minutes on and 30 minutes pause, or
  - 30 minutes on and 60 minutes pause



**The Ice Age Is Over!**  
No more ice, no more hydroculators.




**H<sub>2</sub>O = 2 ounces**  
Uses only 1/4 of a typical glass of water



**Motion Advantage Pads**  
Targeted therapy designed by surgeons.



**Smaller by Design**  
Lighter and smaller for easy portability.



**Broad Temp. Range**  
20 controlled temp settings from 34°–125°F.



**Programmable Therapy**  
5 preset timer modes; heating and cooling therapy all in one device



**Consistent Temp.**  
20 temp levels; maintains temp as long as you need it.

### ThermaZone® Unit



Includes Unit, Power Supply, and Tubing Assembly

ThermaZone Therapy System

003-99 Unit with Medical Grade Power Supply

### ThermaZone® Pads (Reusable)

Individual pads work in conjunction with the **ThermaZone** Thermal Therapy technology and are interchangeable with the unit. Pads are universal and fit right or left side.



003-10 Front & Side Head Pad



003-12 Eye Pad



003-11 Occipital (Back of Head) Pad



003-17 Knee | Elbow Pad



003-47 Range of Motion (ROM) Knee Brace\*  
\*Insurance Reimbursable  
HPCS code L1833



003-15 Shoulder Pad Fits up to a 44" Chest




003-18 Back, Abdomen, Hip Pad




003-22 Universal Pad



003-19 Ankle | 90° Elbow Pad



003-20 Large Shoulder Pad Fits up to a 54" Chest



### Innovative Medical Equipment, LLC

23950 Commerce Park, Suite 100 | Beachwood, OH 44122 USA  
(877) 901-ZONE (9663) | [www.therma-zone.com](http://www.therma-zone.com)



#M006 v5 05/2021