


GE Signa HDx, 16x,16 ch with Coils



- Error Logs
- Diagnostics
- Image Quality
- Calibration
- Configuration
- Utilities
- Replacement
- PM
- Home

Hospital Information
Hospital Name : ""
System Id : 2
Time Zone : America/Chicago CST6CDT US
x86_64
Suite Id : MUK
Language : English
Magnet Serial Number : R3730
Host IP : 10.20.4.5
FE Laptop IP : 192.168.4.4
Gradient Amplifier : 8916
Field Strength : 15000

Release Information
Hostname : lx-mr
Build number for MROS is GEHCMR_1.1-8_cselx1.1_OS
Build number for MrpApps is HD16.0_V03_1638.a
Build number for OS is cselx_1.1
Build number for OnlineHelp is Help_HD16.0_V02_HDxt_English
Build number for PostSdc is HD16.0_V03_1638.a
Build number for TERRA is IOS_20110510_HDx2
Build number for cclass is HD16.0_V03_1638.a
Build number for doc is 15.0-HD16.0_DOC_HDxt-1.5T-5440001-3-rev8.0
Build number for driverSupport is HD16.0_V03_1638.a
Build number for insite is HD16.0_V03_1638.a
Build number for install is HD16.0_V03_1638.a
Build number for mx_apps is HD16.0_V03_1638.a
Build number for mx_plat is HD16.0_V03_1638.a
Build number for pc is HD16.0_V03_1638.a
MR Software release: HD16.0_1638a
Install in Spec Key : Installed

System Health Information
 The system health check is not available under non-proprietary mode

System Information
[0:0:0:0] disk SEAGATE ST3300657SS HPS1 /dev/sda
[0:0:1:0] disk SEAGATE ST3300657SS HPS1 /dev/sdb
[0:0:2:0] disk SEAGATE ST3300657SS HPS1 /dev/sdc
[1:0:0:0] cd/dvd hp DVD-RAM GH80N RF04 /dev/sr0

C Shell Unmount CD

Update

Guided Install MR Software release: HD16.0_1630a

ation Help

Hardware Configure

Field Strength	1.5T	
Package Name	ECHOSPEED SR120	
Gradient Type	8916(FFD)	
Resonance Module	BRM	
Magnet Serial Number	R3730	
Table Limit	31750	
Line Frequency	60Hz	
RF Amp Type	1.5T SRFD2 RF amp	
Governing Body	Fda2	
ISO Vector Z	8981 CX / LCC Magnet Wide Open Encl.	8909
Magnet Enclosure	WideOpen (light IQ)	
PAC Type	PAC 2B	
RF-Systems Cabinet(RFS)	no	
Spectro RFamp Type	None	
Shim supply type	None	More info
Lpca Configuration Type	Legacy	More info

Configure Previous Next Home

Page

Option Key

Select	Key Name	Key Type	Expiration Date
<input type="checkbox"/>	2D Fat Sat Fiesta	permanent	0
<input type="checkbox"/>	MR-Echo	Not Installed	-----
<input type="checkbox"/>	8 Channel System	Not Installed	-----
<input type="checkbox"/>	16 Channel System	permanent	0
<input type="checkbox"/>	Lava	permanent	0
<input type="checkbox"/>	FiberTrak	Not Installed	-----
<input type="checkbox"/>	FLAIR3D	permanent	0
<input type="checkbox"/>	Multi-Phase (variable delays)	permanent	0
<input type="checkbox"/>	Dual Array	Not Installed	-----
<input type="checkbox"/>	Timeout Application	Not Installed	-----
<input type="checkbox"/>	Blood Flow and Volume Measure...	permanent	0

Install/Uninstall Option Keys

Manual

License #:

Uninstall

eLicense CD

Install

Install Temporary Service Key

Host ID

148444788

Service License #:

Install

Configure

Previous

Next

Home

COIL, 3 CH PHASED ARRAY SHOULDER

COIL, QUAD KNEE

COIL, 4 CH PHASED ARRAY SHOULDER

COIL, 4 CH PHASED ARRAY SMALL EXTREMITY

COIL, 8 CH CTL

COIL, 8 CH NEURO VASCULAR ARRAY

COIL, 1.5T QUAD SPLIT HEAD ASSEMBLY

COIL, BODY ARRAY, POSTERIOR

COIL, HD 8 CH FOOT/ANKLE

HIGH RESOLUTION WRIST COIL

HD 8 CHANNEL HI RES BRAIN ARRAY COIL