


body-jet® evo

Fast, Versatile and Intuitive Body Contouring

body-jet® evo offers premium water-jet assisted liposuction for fat transfer and body contouring in-office and clinic. With the help of a fan-shaped, pulsating water spray, fat is gently separated and aspirated simultaneously in a closed system.

The innovative body-jet® technology is designed for **the fastest, most gentle and effective removal of small or large fat volumes** using the water-jet assisted liposuction (WAL) method. A fan-shaped, pulsating water-jet combines with a tissue-sparing, defined vacuum allowing for the **harvest of intact and highly viable fat cells in a closed system.**

- **No more tumescence technique**
- **No more change of cannulas for spray or aspiration**
- **No more centrifugation required**




Using the gentle force of water-jet, the adipose tissue is detached resulting in better preserved fat cells and surrounding tissue. This results in **higher patient safety, faster recovery and less pain.** Learn more today about the efficient, versatile and intuitive body-jet® evo.



Kwiatkowska K, Krapohl BD, Tanzella U, Ueberreiter K. Long-term clinical results and quality of life in patients undergoing autologous fat transplantation for breast augmentation using the BEAULI™ protocol. GMS Interdiscip Plast Reconstr Surg DGPW. 2019;8:Doc10. DOI: 10.3205/iprs000136, URN: urn:nbn:de:0183-iprs0001367.




body-jet® Harvesting Technology:



Phase 1 | The fan-shaped spray gently loosens fat for in-office harvesting through a closed loop system.



Phase 2 | Simultaneously separate and suction with the two-in-one cannula technology for more efficient aspirating.



body-jet® evo CUTTING-EDGE TECHNOLOGY


Precision high volume adipose tissue harvesting and body contouring.

Thanks to a novel technology for the body-jet® evo, the **intensity and effectiveness** of the water-jet can be adjusted to the users needs, **ensuring a virtually pain-free treatment**. The **surrounding connective tissue, nerves and vascular structures remain undamaged** resulting in higher patient safety and faster healing with less pain and bruising.

The results of water-jet assisted lipofilling using the body-jet® technology are **predictable**, as proven by a multitude of published clinical studies, and feature a user-friendly interface for simple setting adjustments.



Reach out today for more information or to connect with a local representative.



| body-jet® evo | |
|---------------|------------------------------------|
| DIMENSIONS | 54 cm (W) x 60 cm (D) x 178 cm (H) |
| WEIGHT | 65 kg |
| FLOW | Approximately 60 ml/min |
| VACUUM | Approximately -500 mbar |
| POWER SUPPLY | 230 VAC, 6 A, 50 Hz |
| PULSE EFFECTS | Short, Medium, Long |

Premium accessories for a wide field of applications.




LipoCollector®



FillerCollector®



Applicators



Cannulas