

It's What Platelets do that Matters

CAREstream America's Liberty MTX System converts PRP to PRFM through a very simple closed device. While other systems are also outstanding the below demonstrates how Liberty MTX is a better solution for your practice providing you with superior outcomes:

1. PRFM's physical and performance characteristics are substantially different than those of PRP. Our system seamlessly combines PRP with calcium to extend the platelet release curve from 4 hours to 2 weeks. PRFM remains in the injection site setting up a fibrin matrix as scaffolding. There are also options to use PRFM as a membrane of varying densities.
2. Our centrifuge is yours at no-charge. Kit cost is extremely reasonable; with generous volume discounts and/or no-initial-cost consignment.
3. Liberty MTX requires an 18ml blood draw; yielding 9ml of PRP.
4. Liberty MTX requires a 7-minute spin with no extra steps or preparation.
5. The Liberty MTX gel separator creates an automated process. No steps, settings, etc. are required.
6. Liberty MTX removes all red blood cells and most neutrophils. The new cells that you're trying to regenerate will not be exposed to a destructive environment and injection site pain is greatly minimized. Inflammation is part of the healing process, but not this part of it. Liberty MTX creates an ideal regenerative/anabolic environment. We separate components physically, scientifically and seamlessly; not just visually.
7. The minimal preparation process of Liberty MTX means that fewer platelets will be damaged; producing a more viable yield.
8. You will produce better outcomes with Liberty MTX as well as a more pleasant patient and practitioner experience.

Cory Steiner

(800) 366-8051

eortho.com

(215) 791-0674 cell

eortho@aol.com

(215) 321-2329 home