



Tropocells[®] PRP

Autologous platelet rich plasma & cell therapy technologies

Tropocells PRP Volume by Treatment Site

Treatment Site	PRP Volume	Volume of PPP to Remove	# Tropocells Kits Needed
Finger and Toe Joint	1 ml	5 ml	One 11 ml Kit
Wrist, Ankle, Ligament, Tendon, Facet Joint, Bursae	2 ml	4 ml	One 11 ml Kit
Shoulder, Hip, SI Joint, Knee	3 ml	3 ml	One 11 ml Kit

Notes:

- In general, one 11 ml Tropocells PRP Kit will produce a quantity of PRP for a single treatment site and one 22 ml kit will produce enough for two treatment sites.
- Multiple Treatment Site Examples:
 - 2 Knees - total of 6 ml of PRP needed (one 22 ml kit)
 - 1 Knee + 1 Shoulder - total of 6 ml of PRP needed (one 22 ml kit)
 - 1 Knee + 1 patellar tendon - total of 5 ml of PRP needed (one 22 ml kit)



CORY A. STEINER
EASTERN ORTHOPAEDICS
OFFICE: (800) 366-8051
www.eortho.com