

# Alternating Pressure Air Mattress

## Features

- Alternating pressure, two groups of air chambers alternate pressure on a regular basis to distribute body weight and prevent bedsores;
- Detachable cell, each air bladder can be removed and used independently, easy to replace;
- Adjustable pressure, 10 pressure levels make it easily adjust a tailored pressure according to your preference;
- Low noise, "whisper quiet" pump produces consistent air flow without interference with sleep or rest;
- Breathable medical PVC material allows heat to dissipate for the patient comfort, is easy to clean and durably designed, safely supporting up to 330 Lbs.;

### Air Mattress

Model No.	KRP/RCD-B-01	KRP/RCD-B-02
Dimension (L*W*D inflated)	190*85*8cm	190*90*7.5cm
Weight capacity	150kg (330 lbs.)	135kg (300 lbs.)
No. of air cells	24 pcs	22 pcs
Cell material	PVC coated nylon fabric (0.35mm)	PVC coated nylon fabric (0.32mm)
Pressure mode	Alternating pressure	Alternating pressure
Cover material	PU	—
Warranty	12 months	12 months

### Pump Unit

Model No.	KRP/RCD-B-01	KRP/RCD-B-02
Input voltage	110V/220V 50-60HZ	220V 50-60HZ
Pump output pressure	50-120mmHg	50-120mmHg
Cycle times	12mins	10mins
Air Output	5-7L/min	5-7L/min
Noise	30 dB	45 dB
Dimension	23.5*12*9.5cm	25*12*10cm
Construction	ABS	ABS
Warranty	18 months	12 months



Made up of PVC-coated nylon material, durable and not easy to age, easy for clean and maintenance.



Alternating



Low Noise



Breathable

