



# Heat and Cold Therapy System



Cold Therapy



Heat Therapy



Control Temperature




Semi-Conductor



No Ice Needed

This is a cutting edge and professional heat and cold therapy recovery system. The system integrates cold therapy and heat therapy in one device by using semiconductor technology. The equipment is easy to use with safety protections. It brings convenient and flexible treatment options to the physical therapists, sports medicine experts, athletic trainers and professional athletes.



## Features

- Integrate cold and heat therapy solution for different stages of rehabilitation which help to accelerate recovery;
- Be free from ice and refrigeration to stay cool, allowing for longer and more-sustained cold therapy;
- Support with high and low temperature protection, various detections and automatic drainage and shutdown to ensure an easy and a safe operation;
- Provide treatment presets for a quick start, and the control can be in your hands with custom temperature and treatment time;





#prestigecare



New



 4°C~18°C

 35°C~41°C



Input Voltage	Power	Water Reservoir Capacity	Net Weight	Dimension
110-240V AC 50/60Hz	420W	Max. 500ml	16.5 Kgs	415*225*470mm

#prestigecare