

PRESTIGE

WITH OVER 30 YEARS OF EXPERIENCE






KNEE BRACE SUPPORT

3x

MUSCLE POWER
SPORT PERFORMANCE
ENDURANCE
STABILITY
ARTHRITIS CARE
BREATHABLE



- ✔ Original 100% Prestige Product
- ✔ Trusted By Health Professionals
- ✔ Over 30 Years of Establishment
- ✔ High Quality Construction
- ✔ Ensured Recovery




TESTED APPROVED | CERTIFIED AUTHORIZED

KNEE BRACE- BW1



Supply Chain

INDICATIONS

Recommended for acute mild to moderate anterior knee pain, support for patella instability, lateral subluxation and maltracking.

CHARACTERISTICS

- Patented dual durometer (silicone) buttress system provides consistent corrective force on the patella
- Bifurcated strap dynamically pulls the patella into proper alignment
- Semi-flexible removable nylon hinges included
- 21cm long support giving maximum power to pull patella into alignment
- Fabric is hypoallergenic, breathable, and lightweight
- 6 month warranty

ORDERING INFORMATION

Part Number			
Left	Right	Thigh Circumference	Size
11-0261-1	11-0260-1	33cm-39cm	XS
11-0261-2	11-0260-2	39cm-47cm	S
11-0261-3	11-0260-3	47cm-53cm	M
11-0261-4	11-0260-4	53cm-60cm	L
11-0261-5	11-0260-5	60cm-67cm	XL
11-0261-6	11-0260-6	67cm-75cm	XXL
11-0261-7	11-0260-7	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella

#prestigecare