

AQA Required Practical:

MAKING A SOLUBLE SALT

Learn how to make a pure, dry soluble salt in the lab!

EQUIPMENT NEEDED:



Beaker & Bunsen Burner



Acid



Insoluble Base



Filter Paper & Funnel



METHOD:

1. ADD ACID TO BEAKER



2. HEAT & STIR



3. ADD BASE UNTIL REACTION COMPLETE



4. FILTER THE MIXTURE



5. EVAPORATE TO CRYSTALLISE



COLLECT YOUR PURE, DRY SALT CRYSTALS!

