

Surgistool®



Three-spindle base provides superior control and pedal access

Dual-density foam seat and adjustable backrest offer maximum support

Optional padded armrests allow versatile positioning and stability

Multiple seat configurations offer personalized comfort

The Ideal Surgeon's Stool

The Surgistool's customizable seat configurations and fully-adjustable backrest and armrests provide optimal support to meet your specific positioning preference. Superior pedal access and a stable yet mobile design combine to deliver the best possible clinical performance.




Teardrop Seat



Wedge Seat



Saddle Seat (Optional)



Teardrop Seat



Wedge Seat



Saddle Seat (Optional)

Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

Specifications

• Model Number	830	
• Weight Capacity	300 lbs	(136 kg)
• Height Range		
High	28.5"	(72 cm)
Low	21"	(53 cm)
• Height Range With Riser Kit		
High	30.5"	(77 cm)
Low	23"	(58 cm)
• Seat Cushion Style	Teardrop	
	Wedge	
	Saddle	

Stryker reserves the right to change specifications without notice.

Optional Features

- 2" (5 cm) riser kit
- Saddle seat
- Multi-position armrest