

# Eye Surgery Stretcher and Surgistool® Chair



## Eye Surgery Stretcher

### Access. Simple. Functionality.

**Stryker's Eye Surgery Stretcher** offers maximum headend clearance, giving surgeons ample leg room to maneuver while performing delicate eye procedures. The reliable work surface and mobile design combine to help with surgical throughput.



### Enhanced head end clearance

Maximizes the surgeon's access for either a temporal or a superior approach to surgery.



**Optional contoured or flathead piece** comfortably supports patient during procedures.

**Retractable fifth wheel** helps improve mobility.

**Metal tooth and groove braking system** prevents swivel and roll.

**Option of two- or three-position fold-down siderails.**

**Multiple bumper color options** allow you to customize the stretcher for identification purposes.



Navy    Gray    Red

## Surgistool Chair

### The ideal surgeon's stool

**Surgistool Chair** offers customizable seat configurations with fully adjustable backrest and armrests provide optimal support to meet the surgeon's specific positioning preferences. Superior pedal access and a stable, yet mobile, design helps to facilitate clinical performance.



**Three spindle base** provides superior control and pedal access.



**Optional padded armrests** allows versatile positioning and stability.

**Multiple seat configurations** offer personalized comfort.



Teardrop seat



Wedge seat



Saddle seat  
(optional)

# Eye Surgery Stretcher and Surgistool Chair



## Standard features: Eye Surgery Stretcher

- 3" (8 cm) Enhanced Comfort Mattress
- 8" (20 cm) Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake and steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

## Standard features: Surgistool Chair

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Surgistool, Stryker. All other trademarks are trademarks of their respective owners or holder.

Copyright © 2017 Stryker  
Mkt Lit-63 01262005 Rev H.2

3800 E. Centre Avenue  
Portage, MI 49002 USA  
t: 269 329 2100  
toll free: 800 327 0770

[patienthandling.stryker.com](http://patienthandling.stryker.com)

## Specifications

|                           | Eye Surgery Stretcher            | Surgistool Chair                               |
|---------------------------|----------------------------------|--|
| <b>Model number</b>       | 1089                             | 830  |
| <b>Overall length</b>     | 90" (228.6 cm)                   |  |
| <b>Overall width</b>      | 31.5" (80 cm)                    |  |
| <b>Weight capacity</b>    | 500 lb. (225 kg)                 | 300 lb. (136 kg)                               |
| <b>Height range</b>       |                                  |  |
| High                      | 34" (86.4 cm) (to litter top)    | 28.5" (72 cm) / 30.5" (77 cm) (with riser kit) |
| Low                       | 22.25" (56.5 cm) (to litter top) | 21" (53 cm) / 23" (58 cm) (with riser kit)     |
| <b>Seat cushion style</b> |                                  | Teardrop / wedge / saddle                      |
| <b>Litter positioning</b> |                                  |  |
| Backrest                  | 0°– 90°                          |  |
| Knee gatch                | 0°– 30°                          |  |
| Trend./Reverse Trend.     | ± 18°                            |  |
| <b>Patient surface</b>    | 26" x 87" (66 x 221 cm)          |  |
| <b>Siderails</b>          | 13" x 55" (33 x 139.5 cm)        |  |
| <b>Caster diameter</b>    | 8" (20 cm)                       |  |

Stryker reserves the right to change specifications without notice.

## Optional features: Eye Surgery Stretcher

- 4" (10 cm) Ultra Comfort Mattress
- Pre/post-op shoulder pads assembly (crank backrest only)
- Air delivery system/drape support
- Crank backrest
- Crank knee gatch
- Defibrillator tray/foot extender
- Temporal or superior wristrest
- Footboard/chartholder
- Three-sided hydraulic controls
- Three-sided brake and steer pedals
- Integrated patient transfer/arm board
- IV caddy
- Restraint straps
- Serving/instrument tray
- Siderail pads
- Surgery accessory rail
- Three-position siderails
- Upright oxygen bottle holder
- Oxygen bottle retaining ring
- Serving/instrument tray holder
- Head pads - Ultra Comfort
  - 3" (8 cm), flat
  - 3" (8 cm), concave
- IV poles
  - Permanent folding two- or three-stage at foot end
  - Removable
- Extended warranty available

## Optional features: Surgistool Chair

- 2" (5 cm) riser kit
- Saddle seat
- Multi-position armrest