

# Surgistool®

Medical



**Three-spindle base provides superior control and pedal access**

**Dual-density foam seat and adjustable backrest offer maximum support**

**Optional padded armrests allow versatile positioning and stability**

**Multiple seat configurations offer personalized comfort**

### The Ideal Surgeon's Stool

The Surgistool's customizable seat configurations and fully-adjustable backrest and armrests provide optimal support to meet your specific positioning preference. Superior pedal access and a stable yet mobile design combine to deliver the best possible clinical performance.




Teardrop Seat



Wedge Seat



Saddle Seat (Optional)



Teardrop Seat



Wedge Seat



Saddle Seat (Optional)

## Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

[www.stryker.com](http://www.stryker.com)

## Standard Features

- Adjustable contoured backrest
- Dual-density foam seat cushions (except saddle seat option)
- Foot-controlled, hydraulic, three-spindle base
- Front-end locking system
- Thermoformed plastic hood
- Teardrop or wedge seat

## Specifications

• <b>Model Number</b>	<b>830</b>	
• <b>Weight Capacity</b>	300 lbs	(136 kg)
• <b>Height Range</b>		
High	28.5"	(72 cm)
Low	21"	(53 cm)
• <b>Height Range With Riser Kit</b>		
High	30.5"	(77 cm)
Low	23"	(58 cm)
• <b>Seat Cushion Style</b>	Teardrop	
	Wedge	
	Saddle	

*Stryker reserves the right to change specifications without notice.*

## Optional Features

- 2" (5 cm) riser kit
- Saddle seat
- Multi-position armrest