

The next generation of warming cabinets

D-Series Data Logging Blanket & Fluid Warming Cabinets



The D-Series (Data Logging) Blanket & Fluid Warming Cabinets provide independent, digitally controlled heating chambers that offer actual temperature and set point displays. Data is recorded every 30 minutes and stored for up to two years. Choose from single, dual, or triple chamber units.

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at **sales@macmedical.com**.

STANDARD FEATURES

- "Easy Log" integrated data logger with removable USB – *now standard!*
- Keyed door lock for single & upper chamber units – *now standard!*
- Adjustable feet – *now standard!*
- Automatically records data every 30 minutes
- Internal storage up to 2 years
- Displays actual temperature and set points at all times
- Independent digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Keyed temperature lock out
- Push button digital operation
- Large digital display(s)
- Fully insulated to provide uniform heating
- Audible & visual over-temp alarm
- Remote alarm output

OPTIONAL FEATURES

- Full view glass door
- Left hand door hinge
- Mobile base or stand
- 220/240V power supply
- Celsius
- Two (2) recessed unit options
- Seismic braces
- Roll out basket(s) with cool touch handles
- Pass through chamber
- Sloped top
- Keyed door lock for intermediate & lower chamber units
- Electronic keypad lock
- Custom sizes available



FDA
Listed

OSHPD
Pre-Approved



MADE
IN THE USA

How to Order D-Series


Model #	Overall Size	Chambers	Door Type	Door Hinge	Base Style	Other Options	
SWC151822	15.5"D x 18"W x 22"H	Single	G = Glass Door	LH = Left Hinge	† 2B = 2" Base	220 = 220/240V Power Supply	
SWC182424	20.5"D x 24"W x 24.5"H	Single	Blank = Stainless Steel Door (Standard)	Blank = Right Hinge (Standard)	† 4B = 4" Base	C = Celsius	
SWC242424	26.5"D x 24"W x 24.5"H	Single			† 6B = 6" Base	R1 = Recessed Unit with Insulation Wrap, No Top Panel, No Side Panels, and No Trim Kit	
SWC183024	20.5"D x 30"W x 24.5"H	Single			NB = No Base	R2 = Recessed Unit with Top Panel, Side Panels, and Trim Kit	
SWC243024	26.5"D x 30"W x 24.5"H	Single			MB = Mobile Base		
SWC183036	20.5"D x 30"W x 36"H	Single			SB = Seismic Braces	WB = Roll Out Basket	
SWC243036	26.5"D x 30"W x 36"H	Single					
SWC182464	20.5"D x 24"W x 64.75"H	Single					P = Pass Through Chamber
SWC242464	26.5"D x 24"W x 64.75"H	Single					
SWC183064	20.5"D x 30"W x 64.75"H	Single					
SWC243064	26.5"D x 30"W x 64.75"H	Single					
SWC182474	20.5"D x 24"W x 74.5"H	Single					
SWC183074	20.5"D x 30"W x 74.5"H	Single					
SWC242474	26.5"D x 24"W x 74.5"H	Single					
SWC243074	26.5"D x 30"W x 74.5"H	Single					
DWC182464T*	20.5"D x 24"W x 64.75"H	Dual			† ST = Sloped Top		
DWC242464T*	26.5"D x 24"W x 64.75"H	Dual			†† IDL = Intermediate Chamber Door Lock		
DWC183064T*	20.5"D x 30"W x 64.75"H	Dual			††† LDL = Lower Chamber Door Lock		
DWC243064T*	26.5"D x 30"W x 64.75"H	Dual			†††† EL = Electronic Keypad Lock		
DWC182474T*	20.5"D x 24"W x 74.5"H	Dual					
DWC242474T*	26.5"D x 24"W x 74.5"H	Dual					
DWC183074T*	20.5"D x 30"W x 74.5"H	Dual					
DWC243074T*	26.5"D x 30"W x 74.5"H	Dual					
DWC183074E**	20.5"D x 30"W x 74.5"H	Dual					
DWC243074E**	26.5"D x 30"W x 74.5"H	Dual					
TWC183078	20.5"D x 30"W x 78.75"H	Triple					
TWC243078	26.5"D x 30"W x 78.75"H	Triple					

the **NEXT** generation of warming cabinets

† Not available for SWC151822
 †† Only available for Triple Chamber Units
 ††† Only available for Dual and Triple Chamber Units
 †††† Only available for Single Chamber Units

* T = Chambers in Thirds

** E = Equal Chambers



Example 1: Model # SWC182464, Glass Door, Left Hinge, Mobile Base, Sloped Top = Order Model # SWC182464-G-LH-MB-ST

Example 2: Model # DWC243074E, Stainless Steel Door, Right Hinge, 4" Base, Celsius, Lower Chamber Door Lock = Order Model # DWC243074E-4B-C-LDL

Note: Additional sizes available for specific requirements. Please contact us for further assistance.



Warranty

All products manufactured by MAC Medical, Inc. will carry a limited lifetime guarantee against product craftsmanship, one year labor and two year parts guarantee. The factory will service all units without cost to the buyer for one year from shipment. After the one year period, replacement of a defective part (labor) will be at buyer's expense. We will exchange all defective parts at no cost to the buyer for a period of two years from shipment. All defective parts must be returned within 30 days to ensure proper credit. An RMA from MAC Medical, Inc. must be obtained prior to items return.



Customer needs are our first priority.

MAC Medical, Inc.
 Printed in USA
 Publication No. FFU-02 Rev D
 August 2016

Information regarding products are subject to change without prior notice.



Corporate Office
 325 West Main Street
 Belleville, IL 62220

Manufacturing Plant
 820 South Mulberry Street
 Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • www.macmedical.com