

D-Series Data Logging Blanket & Fluid Warming Cabinets




Model # DWC183074E

(E = Equal Chambers)


DIMENSIONS

SHOWN WITH
OPTIONAL 4" BASE




DOOR REMOVED
FOR CLARITY


FRONT VIEW



FRONT VIEW



SIDE VIEW



SPECIFICATIONS

Overall Size

20.5"D x 30"W x 74.5"H
(520.7mm x 762mm x 1892mm)
13.04 Cubic Feet Storage
4 Adjustable Shelves

Upper Chamber I.D.

17"D x 26"W x 25"H
(431.8mm x 660.4mm x 635mm)
6.39 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 2° F)

Lower Chamber I.D.

17"D x 26"W x 26"H
(431.8mm x 660.4mm x 660.4mm)
6.65 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 2° F)

Electrical

120V, 11.4 Amp, 50/60 Hz, .82 kWh (Avg.),
2798 BTU/hr (Avg.)

Approximate Capacity

64 Blankets & 72 Solution Bottles (1 Liter Each)

Door Types

Stainless Steel or Glass