

Stryker Prime™ Series

Mobile.
Safe.
Efficient.



Extraordinarily mobile and user-friendly,

the Stryker Prime Series stretchers are designed to greatly reduce the physical strain on clinicians while enhancing patient care and comfort. Advanced mobility options such as the Zoom® Motorized Drive System or Big Wheel® make transporting heavy patients easy, while electric functionality allows patients to adjust their own position without calling a caregiver for assistance. Fully equipped and highly configurable, the Stryker Prime Series allows you to create the ideal solution to meet your specific needs.



Zoom Motorized Drive System on the Stryker Prime stretcher virtually eliminates manual pushing.



Lift Assist™ Backrest uses a drop-seat design which uses the patient's weight to help raise the head of the bed while cradling the patient to help reduce boosting.



Electric functionality includes the ability to give patients control of their own comfort, while Chaperone® Stretcher Exit System helps reduce flight and fall risk.

Advanced mobility, with optional Big Wheel or Zoom Motorized Drive System, helps make transport safe and efficient.

NEW electric-powered base and litter option makes positioning the stretcher and patient as simple as pressing a button.

Industry's first stretcher exit system alerts the nurse to a potential patient fall.

Optional integrated scale system provides accurate, repeatable patient weights.

Pioneer® mattress offers superior pressure redistribution and comfort as a solution for preventative skin care.

Power-washable* design with 700 lb. capacity.

Stryker Prime Series



1105

Standard Features

- Thermoformed ABS base with extra capacity storage tray
- 700 lb. weight capacity
- Less than 21-inch low height
- Integrated utility tray
- Litter-mounted oxygen bottle holder
- Integrated transfer board
- Power-washable design (1105 and 1115 only)
- Four-wheel, center-locking brake system with dual end activators
- Dual, side-mounted foot control hydraulics with uni-lower pedal
- Glideaway™ Siderails
- 4 integrated IV receptacles
- +/- Trendelenburg
- 8-inch Omni Surface™ Casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual end brake/steer controls
- 2-year warranty

1105: All Standard Features Plus

- Fifth Wheel steering

1115: All Standard Features Plus

- Big Wheel steering

1125: All Standard Features Plus

- Zoom Motorized Drive System

Optional Features

- 4-inch (10 cm) Enhanced Comfort Mattress
- 4-inch (10 cm) or 5-inch (13 cm) Ultra Comfort Mattress
- Pioneer Pressure Redistribution Mattress
- Four-sided brake and steer controls (1105 and 1115 only)
- Pop-up push handles (1105 and 1115 only)
- Knee Assist™ No-crank Hydraulic Knee Gatch
- Exclusive Lift Assist Backrest design with drop seat
- Integrated scale system
- Three-sided hydraulic controls (1105 and 1115 only)



1115

Specifications

Model Number

Overall Length

(1105 and 1115)

(1125)

Overall Width

Siderails Up

Siderails Down (1105 and 1115)

Siderails Down (1125)

Patient Surface

Weight Capacity

Height Range

High

Low (1105 and 1115)

Low (1125)

Litter Positioning

Backrest

Knee Gatch

Trend./Reverse Trend. (1105 and 1115)

Trend./Reverse Trend. (1125)

Siderails

Caster Diameter

EMC



1125

26-inch Litter

1105 (Prime with Fifth Wheel) CE

1115 (Prime with Big Wheel) CE

1125 (Prime with Zoom Drive) CE

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

30-inch Litter

1105 (Prime with Fifth Wheel) CE

1115 (Prime with Big Wheel) CE

1125 (Prime with Zoom Drive) CE

85" (216 cm)

86" (218 cm)

34" (86.4 cm)

30.25" (76.8 cm)

30.25" (76.8 cm)

26" x 75.25" (66 x 191 cm)

700 lb. (318 kg)

34" (86.4 cm)

20.75" (53 cm)

23" (58 cm)

0°–90°

0°–40°

±17°

±16°

14" x 57.75" (36 x 147 cm)

8" (20 cm)

IEC 60601-1-2

85" (216 cm)

86" (218 cm)

38" (96.5 cm)

30.5" (77.5 cm)

31" (78.7 cm)

30" x 75.25" (76 x 191 cm)

700 lb. (318 kg)

34" (86.4 cm)

20.75" (53 cm)

23" (58 cm)

0°–90°

0°–40°

±17°

±16°

14" x 57.75" (36 x 147 cm)

8" (20 cm)

Warranty

Two years parts, labor and travel

Electric Functionality

- Electric option (only 30 inches wide)
- Electric patient controls
- Foot end nursing controls
- Powered lift (1105 and 1115 only)
- Chaperone Stretcher Exit System
- 120V or 240V available

Stryker reserves the right to change specifications without notice.

FSC LOGO