

Standard Features

- Fully welded stainless steel construction
- Height adjustment range - 36" - 58"
- 2" Swivel Casters



MYO-2000



MYO-3000



MYO-4000



MYO-5000

Model	Tray Size	Height Adjustment Method	Weight Capacity
MYO-2000	21 1/4" x 16 1/4"	Foot Friction Lock	50 lbs (dual post)
MYO-2001	19 1/8" x 12 5/8"	Foot Friction Lock	50 lbs (dual post)
MYO-2002	25" x 20"	Foot Friction Lock	50 lbs (dual post)
MYO-3000	21 1/4" x 16 1/4"	Hand Friction Lock	50 lbs (dual post)
MYO-3001	19 1/8" x 12 5/8"	Hand Friction Lock	50 lbs (dual post)
MYO-3002	25" x 20"	Hand Friction Lock	50 lbs (dual post)
MYO-4000	21 1/4" x 16 1/4"	Manual Knob	40 lbs (single post)
MYO-4001	19 1/8" x 12 5/8"	Manual Knob	40 lbs (single post)
MYO-5000	21 1/4" x 16 1/4"	Foot Friction Lock	40 lbs (single post)
MYO-5001	19 1/8" x 12 5/8"	Foot Friction Lock	40 lbs (single post)

Optional Feature/ All Models

Mobile Base
(Model # - MB)

