



**Ergonomically
Friendly**

Processing Sinks

Pipe Base Units



STANDARD FEATURES

- Pipe base
- Single, dual, or triple basins
- Fully welded construction with reinforced tops
- Right to left workflow or left to right workflow
- 37" standard height for stationary sinks
- 4" backsplash
- Under shelf
- Adjustable feet
- Drain boards
- Etched gallon markings
- Basin size: 24"W x 17"L x 10"D
- Full CAD rendering of specified sink
- Multiple layouts for individual requirements
- All sinks can be manufactured to individual specs

OPTIONAL FEATURES

- **Ergonomically friendly** electronic adjustable height with four position programming
 - Electronic adjustable height range: 35"-42"
- Wrap around apron
- Mixing valve faucet
- Pre-rinse faucet
- DI water
- Air gun (with 5 tips)
- Water gun (with 8 tips)
- Eye wash
- Shut off valve
- 6" backsplash
- 8" backsplash
- Side splash
- Tubing manifold
- Hand lever drains
- Foot operated drains (not available for adjustable height units)
- Copper drain system (single, dual, and triple)

PS2612035-9

(electronic adjustable height shown with left to right workflow and optional pre-rinse faucet, DI water, air gun, and water gun)



Note: ETL compliance for adjustable height models only.

Standard Pipe Base Models (Stationary)

Base Model #	Overall Dimensions
PS267237-1	26"D x 72"W x 37"H, Single Basin
PS269637-1	26"D x 96"W x 37"H, Dual Basin
PS2612037-1	26"D x 120"W x 37"H, Triple Basin

Note: Additional sizes available for specific requirements. Please contact us for further assistance.

Note: Model # will change based on features added.

Standard Pipe Base Models (Electronic Adjustable Height)

Base Model #	Overall Dimensions
PS267235-9	26"D x 72"W x 35"H, Single Basin
PS269635-9	26"D x 96"W x 35"H, Dual Basin
PS2612035-9	26"D x 120"W x 35"H, Triple Basin

Note: Additional sizes available for specific requirements. Please contact us for further assistance.

Note: Model # will change based on features added.