

The
CureTM
SERIES

20 years of experience in developing Curing Lights that have been engineered and improved over the years, Spring Dental offers the best in LED Curing Light technology! These durable technological marvels are virtually indestructible, have no glass rods, LED is contained in the head, won't chip or break if dropped and are made in the USA!

- **VIRTUALLY INDESTRUCTIBLE • THE LED IS IN THE HEAD**
- **WON'T CHIP OR BREAK IF DROPPED • NO GLASS ROD • EASY TO DISINFECT**

The
CureTM
SERIES TC-3

New!

THE LED CURING LIGHT WITH THE WIDER LENS!

- ONE OF THE LIGHTEST CURING LIGHTS - ONLY 3.1 OZ.
- SLIM & ERGONOMICALLY DESIGNED
- LARGER SIZE 11MM LENS
- DUAL CONTROL SWITCHES
- HIGH/LOW 700-1400 mW/cm² intensities
- RECHARGEABLE LONG-LASTING LITHIUM ION BATTERY (can be easily changed in the office)
- MADE IN USA



TC-3 **\$714⁰⁰**
 MSPR

The
CureTM
SERIES TC-24a

New!

THE LED CURING LIGHT WITH THE WIDER LENS!

- ONE OF THE LIGHTEST CURING LIGHTS ONLY 1.7 OZ.
- SLIM AND ERGONOMICALLY DESIGNED
- LARGER SIZE 11MM LENS
- DUAL CONTROL SWITCHES
- EASY TO DISINFECT
- HIGH/LOW 700-1400 mW/cm² intensities
- AC POWER
- MADE IN USA



TC-24a **\$699⁰⁰**
 MSPR

The
CureTM
SERIES

BEST POWER FOR THE PRICE!

- ONE OF THE SMALLEST AND MOST POWERFUL CORDED LIGHTS
- FEATURES THE LATEST IN HIGH-POWER LED TECHNOLOGY
- ERGONOMIC SHAPE PROVIDES A COMFORTABLE GRIP FROM ANY ANGLE
- 1200 MW/CM2 INTENSITY
- THE LED HAS 4 DIODES AND ALLOWS FOR 6-SECOND CURE
- AC POWER
- MADE IN USA



TC-01 **\$598⁵⁰**
 MSPR

The
Cure²
SERIES

THE SMALLEST, MOST POWERFUL & DURABLE LED CURING LIGHTS!

- ONE OF THE SMALLEST AND MOST POWERFUL CORDLESS LIGHTS
- FEATURES THE LATEST IN HIGH-POWER LED TECHNOLOGY
- 1200 MW/CM2 INTENSITY
- THE LED HAS 4 DIODES AND ALLOWS FOR 6-SECOND CURE
- ERGONOMIC SHAPE PROVIDES A COMFORTABLE GRIP FROM ANY ANGLE
- MADE IN USA



TC-CL II **\$642⁰⁰**
 MSPR



BUY ANY SPRING CURING LIGHT
AND GET 3 - 25/PACKS
OF SPRING REGULAR FG DIAMONDS FREE!!!
MIX & MATCH A \$179 VALUE*

***Redemption Information:**

Email a copy of your invoice dated 7/1/ - 9/30/19 to customercare@springhealthproducts.com. Send your name address and choice of diamonds. Must be received by 10/15/19.