

Turbo FG Diamonds/Regular FG Diamonds

												
Shape No.	368T	863T	880T	850T	855T	850T	855T	856T	847T	847T	837KRT	850KRT
Size (1/10 mm)	023	014	010	014	016	018	018	018	016	018	012	018
Length (mm)	5.0	8.0	7.0	10.0	7.0	10.0	7.0	9.0	9.0	8.0	8.0	9.0
Order No.	285.5T	265.8T	KS0T	771.10T	772.7T	772.10T	773.7T	773.9T	702.9T	703.8T	837KRT012	850KRT018



To order, specify:

Order Number
Grit: (see Grit Legend below)
Package Count: (10pk, 25pk, 50pk, 100pk)
Ex: 285.5 in fine, 25 pack = **285.5F-25pk**

	Inverted Cone			Double Cone		Flat End Cylinder					Pear			Round End Cylinder				
																		
Shape No.	805	805	805	813	806	835	835	835	S835	835	836	837	830	830	838	880	880	881
Size (1/10 mm)	010	014	016	014	012	008	009	010	010	012	012	012	010	014	012	010	012	014
Length (mm)	1.0	2.0	2.0	2.0	2.5	4.0	3.0	4.0	4.0	4.0	5.0	7.0	2.2	3.0	4.0	7.0	7.0	8.0
Order No.	312	313	320.2	1220	1250	514.4	541	835-010	S835-010	835-012	515.5	515.7	331	365.4	GW2	KS0	KS1	KS2
Grit	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■

	Round End Taper											Flat End Taper								
																				
Shape No.	850	855	856	850	855	S855	856	850	855	856	856	856		847	847	848	847	848	847	848
Size (1/10 mm)	012	012	014	014	016	016	016	018	018	018	021	023		012	014	016	016	018	018	023
Length (mm)	11/0.5	7/0.7	8/0.9	10/0.8	7/1.0	7/1.0	8/1.1	10/1.1	7/1.4	9/1.4	8/1.6	7/1.8		9/0.6	9/0.7	11/0.8	9/1.1	10/1.1	8/1.1	10/1.4
Order No.	770.11	771.7	771.8	771.10	772.7	S772.7	772.8	772.10	773.7	773.9	774.8	775.7		700.9	701.9	701.11	702.9	702.10	703.8	703.10
Grit	   	 	  	  	  	  	  	  	  	  	  	  		  	  	  	  	  	  	

																				
Shape No.	879	862	863	863	862	863	863	878K	879K	878	879	878	878K	878K	878K	879K	801	801	801	801
Size (1/10 mm)	018	012	012	014	016	016	018	012	012	012	012	014	014	016	018	018	010	014	018	023
Length (mm)	9.0	8.0	10.0	8.0	8.0	10.0	9.0	8.0	10.0	8.0	10.0	8.0	9.0	8.0	8.0	11.0	1.0	1.4	1.6	2.3
Order No.	256.9	260.8	260.10	265.8	267.8	267.10	270.9	1712.8	1712.10	250.8	250.10	251.8	254SB	257SB	258SB	257SC	115	120	125	135
Grit	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■	■ ■ ■

	Occlusal Reduction / Finishing						Wheel		Needle					Mosquito		Modified Shoulder		Bulk Reducer	Depth Cutter 0.5 mm
																			
Shape No.	379	368	368	368	811	811	909	909	889	S858	858	858	859	852	392	837KR	850KR	6055	834
Size (1/10 mm)	023	016	018	023	033	037	035	043	008	008	008	014	014	023	016	012	018	018	021
Length (mm)	4.5	3.0	5.0	5.0	4.2	7.0	1.1	1.3	3.0	3.0	3.0	8.0	10.0	10.0	5.0	8.0	9.0	8.0	6.0
Order No.	275.5	292.3	280.5	285.5	230	234	860	863	201.3	S1300	1310.3	1314.8	1314.10	209.10	392 016	837KR012	850KR018	ST8	834 021
Grit	  	  	  	  	 	 	 	 	 	 	 	  	  	 	 	 	 	 	

GRIT LEGEND	
Extra Fine 30µ=EF	■
Fine 100µ=F	■
Medium 130µ=M	■
Coarse 145µ=C	■
Super Coarse 180µ=SC	■



Package Count
Available in 10, 25, 50, and 100



Heat Dissipates

