

QuickVue® At-Home COVID-19 Antigen Test

Easy to Use - At-Home Testing

QuickVue® At-Home OTC COVID-19 tests do not require a machine reader and take 10 minutes to get results. The test uses a gentle self-collected anterior nasal (nares) swab sample to determine a positive or negative COVID-19 result. The swab is swirled in a tube of reagent solution, then removed, before a test strip is inserted.

FDA Emergency Use Authorized (EUA)

The QuickVue® At-Home OTC COVID-19 Test has been authorized by the FDA under an Emergency Use Authorization (EUA) for the detection of proteins from SARS-CoV-2.

Available Now

We offer reliable production and delivery schedules. Our team is committed to providing the highest quality service at the best possible prices.



Features


- Identify your current COVID-19 infection status
- Reliable results in 10 minutes
- Anterior nasal specimen collection
- Fast and easy to use at home
- Qualitative detection of SARS-CoV-2 antigens
- Storage conditions: room temperature in a place out of direct sunlight and out of reach of children
- All-in-one package, including collection swabs

Kit Contents

- 25 test devices
- 25 sterile anterior nasal swabs
- 25 extraction reagent tubes
- Package insert with simple instructions

Performance Measures

- 83.5% Sensitivity
- 99.2% Specificity



For more information or to place an order, call 702-340-0562 email invest@rkl2holding.biz

QuickVue® At-Home COVID-19 Antigen Test

Product Specifications and Dimensions

Kit: 25 tests / <1lb

Case: 10 kits / 250 tests / 7 lbs. / 17" x 14" x 16"

Attachments and Documentation (click below)

[FDA EUA Letter](#)

[FDA Instructions for Use](#)

[Package Insert](#)

[Health Care Provider Fact Sheet](#)

[Patient Fact Sheet](#)

For more information or to place an order, call 702-340-0562 email invest@rkl2holding.biz