




# BIONIK P13



Hearing aid type: behind-the-ear hearing aid Battery type: 13th zinc-air battery  
Charging port: / Charging time: / Use time: 300H Single ear mass: 3.8g  
Overall mass: /



## Scope of application



	IEC 60118-7 / ANSI S3.22 - 2003
Peak OSPL90(dB SPL)	≤138dB
OSPL90-HFA(dB SPL)	120±5dB
Peak Gain(dB )	≥65dB
HFA/FOG Gain(dB )	50±4dB
Equivalent input noise(dB)	≤32DB
Frequency response range(Hz)	F1≤750HZ F2≥5000HZ
Harmonic distortion	≤5%
Battery drain	≤2mA

## Main features

- Programmable**  
Fits a variety of hearing loss curves
- Digital Chip**  
Listen to the good and the comfort
- Zinc-air battery**  
Stable discharge, long battery life
- Comfortable to wear**  
Ergonomic Design
- Volume memory**  
Say goodbye to frequent debugging
- HD sound quality**  
Powerful audio processing
- Independent button**  
Separate buttons for volume and program
- Independent switch**  
Put the battery in and start
- Multiple programs**  
Easily cope with various environments

**BIONIK EVOLUTION LLC**  
**30 N. Gould St. Ste. R**  
**Sheridan, WY USA 82801**