

# 320



Comfort Mobility  
Your Family Home Health Care Store  
519-988-1234

# Stannah

## For Outdoor Stairs

## The Stairlift People

The Stannah 320 is a practical solution where you have a straight staircase and need to access your driveway, deck or yard. Quickly fastened to a few of the stair treads, the 320 is easier to install than a ramp and takes up much less room. The 320 was designed to withstand both hot and cold climates, and has a good track record in locations as diverse as New England and Hawaii.

Operating the Stannah 320 couldn't be easier. After sitting in the secure seat (which comes with a lap belt), move the joystick on the arm-rest to glide effortlessly up or down the stairs. At the top of the stairs you can turn the chair to face the landing. To access the great outdoors, it's just the lift you need!



With its narrow rail and slim design, the 320 stairlift provides a space-saving way to access your driveway, deck or yard.



To use the Stannah 320, simply unfold it, sit down and gently move the joystick on the end of the arm-rest in the direction you want to travel.



Landing controls (with key lock) let you call or send your stairlift - handy for carrying groceries.



The 320 is weatherproof, and also comes with a protective cover to keep out leaves, ice and snow.

## Model 320 Key Features

### Weather-resistant construction

The model 320 is made of light gray waterproof, UV-resistant materials.

### Seatbelt

A lap-belt is fitted as standard.

### Directional Control

Able to be mounted on either armrest, the joystick controls the stairlift's movement.

### Manual swivel seat

Lever on either side of the stairlift unlock the swivel seat, allowing the user to turn towards the top landing. The back of the chair then acts as a barrier across the stairs, for safety.

### Folding seat

The seat may be folded away when not in use.



### Folding armrests

Padded, for comfort

### Power-on light

### On/Off switch

### Key lock

If required, the keyswitch can be used to completely immobilize the stairlift.

### Batteries

Two 12-Volt batteries keep the stairlift working in the event of a power outage. Recharging is continuous and automatic, to overcome cold weather.

### Safety edges

Sensitive safety edges stop the stairlift immediately if it meets an obstruction.

### Landing controls

Hard-wired controls (with key locks) on both landings allow the stairlift to be called or sent.

### Rated load

One person, 300lbs.

## REQUIREMENTS

- Minimum staircase width: 27"
- Maximum user weight: 300lbs
- Operating temperature range: 14°F to 105°F (-10°C to 40°C)
- Power outlet: 110V indoors, with a 1/2" sleeve leading to outdoors. Total distance to stairlift should not exceed 100'
- Stair tread construction: Concrete, solid wood, steel or solid stone.

Family Owned & Operated since 2014

 Comfort Mobility  
Home Health Care  
**519-988-1234**  
2707 Temple Drive, Windsor (Off Central Ave North of EC Row)  
Sales, Service & Rentals  
[www.comfortmobility.ca](http://www.comfortmobility.ca)



**CARP**  
APPROVED

**Stannah 150 1867-2017**

**Stannah Canada Inc.**

**Telephone:** (888) 647-8219

**Email:** [Sales@Stannah.ca](mailto:Sales@Stannah.ca)

**Website:** [www.Stannah.ca](http://www.Stannah.ca)

**Office & Showroom:**

25 Valleywood Dr, Unit 11, Markham, ON L3R 5L9