

Locally Owned
& Operated



PRODUCT CATALOG



519-988-1234
comfortmobility.ca



Comfort Mobility Home Health Care

2707 Temple Drive, Windsor
Off Central Avenue just North of EC Row
Phone 519-988-1234 Fax 519-988-1244

Open Monday – Friday 9:00AM – 4:00PM
(Holiday hours may differ)



Rental Coordinator

Lucy Tamer

rentals@comfortmobility.ca



**Compression, Incontinence &
Bracing Specialist**

Bev Barry

info@comfortmobility.ca



General Manager

Jessica Fase Senay

jessica@comfortmobility.ca



**Mobility Consultant Specializing in
Home Assessment**

Patti Osborne

Cell 226-345-5747

sales@comfortmobility.ca



**Mobility Consultant
Specializing in Retirement & LTC**

Jena Simpson

Cell 519-999-7456

jena@comfortmobility.ca



Locally Owned & Operated

Experience exceptional personalized service and reliability with a small family owned and operated home health care store that will treat you like family!



Comfort Mobility Home Health Care
was Founded by the Fasse Family



STAY SAFE & INDEPENDENT AT HOME WITH



Comfort Mobility Home Health Care

Rejuvenate your home to
make it safe for you and
your loved ones as you age.

What we do:

- Mobility Equipment including Scooters, Wheelchairs, Travel Mobility and more.
- Porch Lifts, Ramps, Stair Lifts
- Bathroom Safety
- Bedroom Safety
- Electric HomeCare Beds
- Compression Stockings
- Incontinence Care Supplies
- And More!




2707 Temple Drive, Windsor
519-988-1234
www.comfortmobility.ca
sales@comfortmobility.ca

Patti Osborne
Mobility Consultant
Home Health Care Professional
Specializing in Home Assessment

Toileting

A



Toilet Safety Frame (also known as Versa Frame)

Available with legs under the seat (figure B), to the floor (figure A) or stand alone (left)

B



C



Toilet Seat Riser (figure C)

Installs between your toilet and existing toilet seat to add additional height. A great choice if you want to keep using your original toilet seat.

D



Raised Toilet Seat (figure D)

Clamps on to your existing toilet replacing the original seat to provide additional height. Available with or without armrests. Standard or Elongated options.

E



Universal Raised Toilet Seat (E)

If you are not sure of your needs, the universal raised toilet seat is height adjustable and works on both standard and elongated toilets!

Helpful Tips

Raised Toilet Seats are available in 2", 4" or 6" heights. Make sure you know what height you require and if your toilet is regular or elongated!

Toileting

A



B



C



D



E



Commodes

Stationary or Wheeled

Regular Arms (A) or
Drop Arm (B) for transferring

Standard or
Heavy Duty/Bariatric (C)

Shower Commodes

Multi Use Commodes for over the toilet, at the
bedside and in the shower.

Options:

Standard (D) or Tilting (E)

Add Soft Padded Seat

Standard or Heavy Duty (Bariatric)

Contact your Comfort Mobility Sales
Representative for assistance choosing the
right shower commode for you!

Also Available

Absorbent Pads or powder for bucket to
reduce odour

Splash Guard for over toilet use

Commode liners



Bathing

A



Extra Long Shower Head with Shut Off Valve (figure A)

If you are often seated in the shower or in your tub using a shower head, this is the tool for you!



B



Clamp On Tub Rail (figure B)

Add a handle to the side of your tub for balance and mild support entering and exiting your bath tub. If you require more support than just balance, it is recommended that a Super Pile with horizontal bar is used.

C



Non-Skid Bath Mat (figure C) or Treads (figure D)

Add grip to your shower or bath tub with a non-skid Bath Mat or tread stickers.

D



Grab Bars

Bath Step (figure D)

While these steps can be used anywhere, they are ideal for a bathroom setting. Stackable to create the ideal height for you.

E



Bath Lift (figure E)

Want to take a bath but have trouble getting in and out? A bath lift is a seat that sits inside your tub and raises up and down to raise and lower you in the tub. The transfer flaps make the transition out of the tub a breeze.

Bathroom Safety



Row 1: Bath Chair, Heavy Duty Bath Chair, Bath Chair with Arms **Row 2:** Bath Bench, Folding Bath Bench, Blown Plastic Bath Chair **Row 3:** Swivel Bath Stool, Transfer Bench, Padded Transfer Bench **Row 4:** 2 Transfer Benches with Curtain Guards, Bath Board


Aids for Daily Living

There's something for everyone at Comfort Mobility. From weighted utensils to shoehorns, whatever you are struggling with in your everyday life there is probably a gadget out there to help you with your independence.

Photos below: weighted utensils, large pill organizers, extra long shoe horn, overall bib, specialty utensils, non-slip socks, handybar, pulse oximeter, extra long bath sponge.



Wellness, Therapy & Exercise



Photos: Blood Pressure Monitors, Electrotherapy TENS Machines and replacement electrodes, Pedal exerciser, ROM Board, Light Therapy, Thermometer

Powered Lift Chairs



A great addition to your home for comfort and reliability. The lifting mechanism makes getting up much easier. With a backup battery, never worry about getting stuck in your chair in a power outage.



Call Comfort Mobility to schedule a free consultation

We will measure you for the most appropriate size and discuss options, fabrics and add on features available

Save the Tax with a valid doctor's prescription!



Visit our Online Store for all up to date models and pricing!

Measuring for a Lift Chair

All Lift Chairs are as varied in size and shape as the people who enjoy them. Taking the following detailed measurements prior to a purchase will ensure a proper fit and deliver maximum satisfaction.

The Importance of Sizing

Comfort

Finding the perfect fit for your new lift chair will provide the most comfortable option available. These are not a one-size fits all.

Safety


Your safety is our top priority! If a lift chair is too big or too small, it can become a safety hazard when transferring in and out of the chair. We want to make sure you find the perfect fit for safe transfers every time you enter and exit your chair.

Longevity

Did you know that if a lift chair is too big or too small it can cause mechanical issues over time? We want to help you find the perfect fit to encourage comfort and not headaches due to mechanical breakdown.

All Lift Chairs sold are warrantied through Comfort Mobility

Please Note: No warranty covers labour.



A – Seat to Top of Back _____ in

Measure the vertical distance from the base of a seat to the top of the user's head while seated.

B – Seat Depth _____ in

Measure the horizontal distance from the back of the seat to the back of the user's knee while seated.

C – Floor to Top of Seat _____ in

Measure the vertical distance from the floor to the top of the back of the user's knee while seated.


D – Width Between Arms _____ in

Measure the horizontal distance between the hips while seated (add 2 inches for wiggle room).

User Height: _____ ft _____ inches

User Weight: _____ lbs

Homecare beds



Luxury Harmony Hi Low Bed

For those who may need a little assistance transitioning from sitting to standing to get out of bed, the Harmony Hi-Low can help you stand with the touch of a button. This is also a great option for those who require patient care at home or in an assisted living environment. The Harmony Hi-Low provides the features of a homecare bed without the institutional look. The removable rear section of the bed frame even allows the user to attach their residential headboard making it a beautiful addition to your room.

Features

Lift capacity of 880 lbs.

NEW! Split Frame Design for easier delivery and installation

NEW! Wireless 21 button hand control that is backlit and features a flashlight for nighttime use

NEW! Tilt Feature

Dual Zone Massage

The bed is 10.5" from the top of the deck to the floor in the low position and raises 11" to a high position of 21.5"

Side rails and Casters are available as an option

Adjustable three piece legs and headboard bracket included

Available in Twin XL, Full, Full XL, Queen and Split King

Find out more on our website including specs, warranty information and pricing

Homecare beds



Tröst Bed by Permobil

The Tröst Bed from Permobil is an electric, adjustable bed ideal for long-term care patients and their care providers. It offers a low height of 6.75" and a high height of 30" to meet all your caregiving and client needs. Tröst Bed has a 12-button hand pendant and optional footboard attendant controls, so customized comfort is an easy adjustment.



WeCare by Joerns

The Joerns WeCare® bed offers all the functions and stability of an institutional frame with an easy to deliver and assemble package. The convenient design combines easy, one-person delivery with ultra-low bed height for fall prevention and a full-range of options to deliver convenient, effective care in any environment.

Standard length 80" Width 35" – Optional 39" or 42"
with CareWide Accessory
Low height 8" (Base model 8¾")
High height 28" (Base model 26")
Underbed clearance 18"
Safe working load 450 lb
Product weight 165 lb
Delivery weight 77 lbs or less
Meets safety standards

Medacure by Human Care

The Medacure Full Electric Homecare bed is lightweight and easy to install. The rich black textured finish and grey accents provide a sleek, contemporary look for a standard homecare bed.

3 Function electric bed, hand pendant raises the head, foot and height of the bed. The mattress deck has reinforced steel slats for safety and comfort with built in mattress keepers. Washable bed ends are easy to clean and maintain. Includes emergency hand crank in case of power failure. Half rails included.

Find out more on our website including specs, warranty information and pricing

Sheet Sets

Most HomeCare Beds
measure 36"W x 80"L

200 Thread Count Cotton-Polyester
Sheet Sets: Cream, Denim Blue, Dark
Grey

Flannel 100% Cotton Sheet Sets:
Grey, Blue Brown

Standard White Polyester-Cotton
Fitted Sheet Set includes 2 fitted
sheets



Bedding in a Box Set contains 1
pillow case, 1 fitted sheet, 1 cover
sheet, and 1 blanket
Sheets Cotton-Polyester
Blanket 100% Cotton



Incontinence Supplies

With your commitment to us, Comfort Mobility will keep your supplies in stock and at the ready for you.

We can arrange a regular delivery schedule or you can simply call when you need more and we will deliver them right to your door!

We also carry...

Wipes

Barrier Cream

Cleansing Cream

Pads and Guards

Reusable and Disposable Bed
Pads

Ostomy and Catheterization
Supplies

The TENA logo, featuring the word 'TENA' in blue capital letters with a registered trademark symbol, and a stylized blue and green swoosh graphic.

Find the right fit with our FREE Sample Program

With the correct sizing and absorbency levels, you should never leak or feel uncomfortable in your incontinence wear. With our free sample program, we will provide you with trials of different options available to you until you find the perfect match.



Bracing & Supports

- Braces
 - Splints
 - Wraps
 - Abdominal Binders
 - Hernia Belt
 - Supports with Stays
 - Stabilizers
 - Immobilizers
 - Cast Boots
 - Cast Protectors
 - Belts
 - Hot/Cold Therapy
- Knee
Shoulder
Ankle
Hip
Spine
Wrist and Thumb
Elbow
Back



Compression Stockings

Key facts about Compression Therapy

Compression therapy plays a central role in the treatment of venous edema and lymphatic disorders. It serves two main purposes: to reduce swelling caused by fluid filtration in the tissues and to reduce venous reflux and venous hypertension. The term "compression" refers to a mechanical pressure exerted from the outside onto skin and tissue.

**Certified Compression Stocking
Fitters on staff.**

Call us today for your free fitting!

**How can compression treatment
help me with my venous disease?**

Compression therapy for legs means wearing special medical compression stockings that come in different lengths and styles, to suit a variety of venous diseases and individual needs.

So, why wear compression stockings? Because they help to normalize the blood flow, to close venous valves better again and ease typical problems such as heavy legs.

For this reason, compression stockings may also help to reduce heavy, swollen legs during pregnancy or long days of sitting or standing at work. Compression therapy is also a good solution for people who regularly take long flights or journeys: compression stockings for travel are specially designed for this purpose.

Medical compression for lymphedema

Compression solutions, such as those offered by JOBST, are also an important component in the management of lymphedema.

Lymphedema is a disorder of the lymphatic system, whereby lymph fluid accumulates in the tissue space and is unable to drain properly due to physical damage to or abnormality of the lymphatic system.

This causes swelling and discomfort in the affected limb.

Compression treatment exerts pressure from the outside onto the tissue, increases the interstitial pressure, and supports the fluid transport and recovery in the lymphatic system. Thereby fluid accumulation is counteracted, which helps to maintain the volume and shape of the affected body part of lymphedema patients..

Ask your doctor for more information about wearing compression stockings. They will be happy to provide detailed advice about the application possibilities of compression solutions and further benefits of compression therapy. Visit or call Comfort Mobility to fill your prescription!

Venous disorders are common, just like high blood pressure, heart disease or diabetes, and should be taken seriously.

Compression Levels OTC (Over the Counter)

8-15 mmHg

15-29 mmHg

20-30 mmHg

40 + mmHg





HEEL PROTECTORS

SP Boot Pillows

Designed to relieve pressure on the heel and ankle. The pillow is made of soft, top-grade fiber that conforms to the shape of the foot, providing a comfortable fit. The interior fabric is also antimicrobial.

Available in Blue or Black/Tan



Choice of polyester pile or vented foam with contact closure. One size fits all. Ideal for Protecting the elbow or heel to aid injured and/or bedridden patients.



Mesh Heel/Elbow Protector Knitted stretch stockinette covered comfort pad to aid in use on heel, elbow or knee. One size fits all.



Made with pressure-relieving polyfoam. With a hook and loop closure, the protector is secure and easy to put on and take off. Its over-sized design is to keep the protector in place. Latex Free



Comfort Mobility
Home Health Care
519-988-1234
comfortmobility.ca



HEEL PROTECTORS

Posey Premium Heel Guard The Posey Premium Heel Guard is an

adjustable foam elevation

block that alleviates pressure on the heel. The soft outer lining allows for firm attachment of ankle and foot drop straps and reduces friction against bed sheets. Straps can crisscross over the foot for maximum foot drop support or attach down the side to minimize contact with the foot. Excellent for personal and professional use and care. Interior is soft for long-term and effective use. Measures 16"L x 11"H in size. Latex free



Skil-Care Heel Guard

Provides pressure relief under heel, ankle and lower calf. Made from plush synthetic sheepskin. Large vent holes promote air circulation, to reduce perspiration and moisture buildup. Secures with VELCRO brand closures. Machine washable Latex free



Skil-Care Heel Float II

To prevent and treat pressure sores, the open heel eliminates heel pressure by suspending the heel in the air. The lining is made from a soft microfiber to minimize friction. For individual fit for each user, the dual straps can be adjusted to any position. The height is extended to keep blankets off toes while still reducing the risk of plantar flexion (foot drop).
DVT compatible Machine washable

Small: circ.- 10"-15", foot - 3" wide

Medium: circ. - 11"-16", foot - 4" wide

Large: circ. - 12"-22", foot - 5" wide



Comfort Mobility
Home Health Care

519-988-1234

2707 Temple Dr Windsor
(Off Central Ave just North of EC Row
opposite North Service Road)

comfortmobility.ca






HEEL PROTECTORS

MaxxcareHeel Pro Evolution

Helps protect against pressure injuries and shear forces. The mechanically strong air tubes allow the user to walk short distances while wearing the Heel Pro Evolution to promote increased mobility. It supports wound healing and is 100% washable.



HeelMedix Advanced Heel Protector

Latex Free, Machine Washable
sold as each



Built-in ultra-fresh antimicrobial protection keeps boot fresh and odor-free*



Screen printing of "toes" and "heel" for proper boot application

NEW
Spandex straps

Color-coded straps pull down on sides help protect against foot drop

FOOT DROP PROTECTION
ATTACH TO HOLE

Screen printed instructions on the strap

Patented side hole design minimizes change of DVT tubing touching patient's skin

PETITE

5-11 inches

calf circumference
(12.7-27.9 cm)

STANDARD

11-19 inches

calf circumference
(27.9-48.2 cm)

X-LARGE

19-25 inches

calf circumference
(48.2-63.5 cm)



Comfort Mobility
Home Health Care

519-988-1234

2707 Temple Dr Windsor
(Off Central Ave just North of EC Row
opposite North Service Road)

comfortmobility.ca





Wheelchairs

Get help with funding from the Assistive Devices Program (ADP)

Comfort Mobility is a trusted ADP Authorized Vendor with over 25 years of experience building custom equipment for clients.

Types of Wheelchairs

Transport Wheelchairs
Manual Wheelchairs
Performance Wheelchairs
Tilt and/or Recline Wheelchairs
Power Chairs

We work together with you, and ADP Authorizers to build the most suitable chair for your needs. Direct billing to Green Shield, WSIB, Veteran's and more!



Veterans Affairs
Canada



Mobility Equipment Purchase Guide

Getting Started

Contact Comfort Mobility!

Call Comfort Mobility to get into contact with your Mobility Consultant or let your therapist know which company you choose to deal with for your equipment needs. Your therapist will call one of our Mobility Consultants to schedule an Assessment to fit you for your mobility equipment based on your current and future needs.

Equipment Selection

Work together with your therapist and vendor to find the right equipment for you. With a Registered Nurse with over 30 years of experience in Home Health Care on staff, Comfort Mobility can provide medical and equipment expertise required to ensure that your clinical requirements are matched with quality and reliable products. We are here to maximize your health care potential with top-of-the-line equipment.

Funding

There are funding agents available for most mobility equipment. Comfort Mobility is an ADP Authorized Vendor. We deal directly with the Assistive Devices Program, Green Shield, Blue Cross Veteran's Affairs, Rotary, WSIB, and more.

If you are unsure of what funding is available to you, speak with your Mobility Consultant about your options. Let us help you make this process stress free.

Order, Set Up and Delivery

We pride ourselves on continuity of care from start to finish. Your Mobility Consultant will be with you every step of the way.

Once your equipment has been selected and ordered, your Mobility Consultant will personally deliver and set up your equipment and make adjustments where required. We will answer any questions you have regarding everyday care and maintenance of your new equipment and provide you with all the necessary paperwork required for the life of your product.

Custom Equipment takes time to build, and we understand that! That's why Comfort Mobility provides complimentary loaner equipment when available until your order arrives. If we are awaiting funding approval and the client portion has not been paid, Comfort Mobility has a wide range of affordable rental equipment available to you at a monthly rate.



Ask us for a copy of our
Simplifying ADP Brochure
to learn all you need to know
about the
Assistive Devices Program



Your Mobility Specialists

Patti Osborne

226-345-5747

sales@comfortmobility.ca

Jena Simpson

519-999-7456

jena@comfortmobility.ca



MINISTRY OF HEALTH AND LONG-TERM CARE

Ontario

ADP - Assistive Devices Program

Who Qualifies?

To qualify, you must:

Be an Ontario Resident

Have a valid Ontario
Health Card

Have a disability
requiring a mobility aid
for 6 months or longer

ADP does not consider
your income.

Who Does Not Qualify?

You do not qualify, if you:

Already qualify for or are
receiving financial support
from the

Workplace Safety and
Insurance Board for the
same mobility aid

Are a Group "A" veteran and
already qualify for or are
receiving financial support
from Veterans Affairs
Canada for the same
mobility aid

Types of Mobility Aids Covered

manual wheelchairs, power wheelchairs and power scooters

power add-ons (devices added to a manual wheelchair if
you don't need a power wheelchair)

positioning devices (e.g. cushions, back and head supports,
power tilt and recline)

forearm-crutches

wheeled walkers for adults

pediatric walkers, standers and strollers

Replacements

You can apply to replace your mobility aid if:

Your medical condition and/or functional ability has changed
and your current mobility aid no longer meets your needs

Your body size has changed and your current mobility aid no longer fits

The mobility aid is worn out, no longer covered by warranty
and cannot be repaired at a reasonable cost

ADP does not cover costs to replace a lost mobility aid or to repair aids
damaged through misuse or neglect.

Equipment sold at Comfort Mobility Inc. has a limited warranty against
manufacturer defects.



How Much Is Covered?

The Assistive Devices Program (ADP) will pay 75% of the ADP price for mobility aids. You must pay the remaining 25% of the cost.

When ADP sets the ADP price for each mobility aid, they take into account: Manufacturer's and/or distributor's costs, Complexity of device. A business can't charge you more than the ADP price for the approved device.

Additional accessories, upgrades or items are not covered by the ADP (e.g. add-on accessories, batteries for power wheelchairs, etc.)

You Have A Choice!

You can choose the company you want for your mobility equipment and services whether you reside at home or in a long-term care facility.

Other mobility companies may offer similar equipment but service is where we excel!

Company/Vendor must be a Registered ADP Vendor

So How Do I get Started?

A formal Assessment by an Occupational Therapist or Physiotherapist registered as an ADP Authorizer is required.
(Fee for Assessment to be paid directly to Therapist)

This will confirm that you need a mobility aid and determine which type best suits your needs.

ADP Application must be completed by you and your Therapist. We Comfort Mobility Inc. will then fill your prescription based on the Therapist's application.



COMFORT MOBILITY

MINISTRY OF HEALTH AND LONG-TERM CARE



Ontario

ADP - Assistive Devices Program

Information Sourced from
the Ontario ADP website

More Funding Sources



With an easy online application, terms up to five years and flexible payment options (either bi-weekly or monthly), we take the stress out of big bills by giving you ample time to pay for your purchase.

Contact Comfort Mobility to set up your application and re-approval process!
or do it yourself here:

<https://applicant.myfrontline.app/LendCare/45684KHQXHSV>



Veterans Affairs
Canada

Provides coverage for assistive devices and accessories that help people do everyday tasks and live independently. Repairs and maintenance are also covered.

Examples of devices or accessories covered include:

- Wheelchairs or walking aids (such as canes, walkers, ice grippers/rubber tips);
- bathroom aids (such as grab bars, raised toilet seats, bath boards, bathtub rails);
- self-help aids that assist with getting dressed, preparing food and living safely in your home.



Ontario Works can help with some of the costs that the Assistive Devices Program does not cover. Ontario Works can also help with the cost of batteries and repairs to assistive devices, such as powered wheelchairs.

Ontario Works provides some benefits for health and other costs, depending on your circumstances.

If you need to apply for Ontario Works, you can [find out if you are eligible and how to apply](https://www.ontario.ca/page/social-assistance) online.
<https://www.ontario.ca/page/social-assistance>



March of Dimes Canada's Home and Vehicle Modification Program provides funding for basic home and/or vehicle modifications. By reducing or eliminating physical barriers which create life safety risks, modifications, adaptations, and devices enable children and adults with mobility restrictions to continue living in their homes, avoid job loss, and participate in their communities.



ODSP offers:

- money to help you and your eligible family members with living expenses, including food, and rent
- health benefits, including prescription drugs and vision care
- employment support to help you find and keep a job or advance your career

If you are eligible for ODSP, the amount of money you get will depend on your specific situation.

More Funding Sources



ALS Canada helps people diagnosed with ALS to cope with the daily challenges of decreasing mobility and communication ability and to help with retaining independence. Access to equipment and assistive devices is essential for your safety, comfort, and functioning as a person living with ALS.

To find out more about the program and see if you are eligible, please review the **Program Overview** on their website.

<https://als.ca/wp-content/uploads/2024/03/Equipment-program-overview-2024.pdf>



The Quality of Life Program provides practical assistance to people, helping them access equipment & assistive technology that promotes independence in their daily lives.

You can apply online at mscanada.ca



Assistance with navigation support as well as financial help through the Equipment for Independence Program.

<https://sciontario.org/services/equipment/>



Helping Kids with
Physical Disabilities
Succeed

Easter Seals Ontario is dedicated to helping children, youth, and young adults with physical disabilities carve their own path into the world.

Visit the Easter Seals website for more information

<https://services.easterseals.org/>



Provides wage-loss benefits, medical coverage and support to help people get back to work after a work-related injury or illness.

Funded by premiums paid by Ontario businesses. WSIB provides no-fault collective liability insurance and access to industry-specific health and safety information.



Muscular Dystrophy Canada's Equipment Program has assisted thousands of Canadians to access the equipment needed to enhance their quality of life. By working with other funders and the government, they share resources to best leverage donor dollars to cover the costs of assistive devices, equipment and accessible technology.

Shop Online

Scooters



A mobility scooter is a great way to enjoy the great outdoors. There are many different styles and sizes to choose from. All scooters are backed by manufacturer warranty.

Stop in our showroom to for a test drive, shop online or call us! Models are ever changing and expanding! Ask us for a current scooter catalogue





Pride Mobility Collection

PRIDE
MOBILITY
LIVE YOUR BEST®



GO GO
ELITE TRAVELLER[®] 2

iTurn Technology
300 lbs Weight Cap
Max Speed 7km/hr
16.7 km per charge
18ah Batteries



GO GO
LX WITH CTS SUSPENSION

CTS Suspension
300 lbs Weight Cap
Max Speed 7km/hr
16.9km per charge
18ah Batteries



GO GO
ENDURANCE Li

325 lbs Weight Cap
Max Speed 6.7 km/hr
10 km per charge
8ah Lithium Battery



REVO 2.0

400 lbs Weight Cap
Max Speed 8 km/hr
28.8 km per charge
U1 Batteries



BAJA BANDIT

400 lbs Weight Cap
Max Speed 12.9 km/hr
20ah: 14.5 km/charge
40ah: 46.7km/charge
20 or 40ah Lithium



VICTORY 10.2

400 lbs Weight Cap
Max Speed 8.4 km/hr
40ah: 21.7 km/charge
U1: 17.7 km/charge
40 ah or U1 Batteries



ZERO TURN 10

400 lbs Weight Cap
Max Speed 12 km/hr
38 km per charge
40 ah Batteries



MAXIMA

500 lbs Weight Cap
Max Speed 9.3 km/hr
3W: 31.8km/charge
4W: 31.3 km/charge
58 ah Batteries



RAPTOR

400 lbs Weight Cap
Max Speed 22 km/hr
See range per charge
specs online
75 ah Batteries



PX4

500 lbs Weight Cap
Max Speed 14.8 km/hr
38 km per charge
75 ah Batteries



WRANGLER

450 lbs Weight Cap
Max Speed 18 km/hr
75ah 32.5km/charge
100ah 42.5 km/charge
75 or 100 ah Batteries



PURSUIT 2

400 lbs Weight Cap
Max Speed 15 km/hr
28 km per charge
50 or 100ah Lithium
Batteries

Note: Range per charge specs are only a guideline and vary on user's weight
For more information and specifications on these scooters, contact us or visit our website!

Golden Technologies Collection

GOLDEN
technologies



375 lbs Weight Cap
Max Speed 8km/hr
28.9km per charge
U1 Batteries



350 lbs Weight Cap
Max Speed 8km/hr
28.9km per charge
U1 Batteries



300 lbs Weight Cap
Max Speed 6.4km/hr
16 km charge
12AH Batteries or
option to upgrade



325 lbs Weight Cap
Max Speed 6.4km/hr
14.4km per charge
12AH Batteries or
option to upgrade



300 lbs Weight Cap
Max Speed 6.4km/hr
11.2km per charge
12AH Batteries



300 lbs Weight Cap
Max Speed 5.8km/hr
28.9km per charge
(1) 15AH Battery



300 lbs Weight Cap
Max Speed 8km/hr
28.9km per charge
22AH Batteries



400 lbs Weight Cap
Max Speed 8.5km/hr
30.4km per charge
U1 Batteries or option
to upgrade



500 lbs Weight Cap
Max Speed 11.2km/hr
28.9km per charge
NF22 Batteries



500 lbs Weight Cap
Max Speed 12.8km/hr
37km per charge
NF22 Batteries or
option to upgrade



BUZZ
AROUND
Carry On

Forbes
HEALTH

NAMES THE BUZZAROUND
CARRYON® AMONG
THE BEST LIGHTWEIGHT
MOBILITY SCOOTERS!

Note: Range per charge specs are only a guideline and vary on user's weight
For more information and specifications on these scooters, contact us or visit our website!



Walkers

**Standard Folding
Walker with or
without Wheels**

Available in Standard
and Heavy Duty
(Bariatric)

Height Adjustable
Made of lightweight
aluminum

5" wheels optional



Rollators

**Complete with 4
Wheels, Hand Brakes
and a Seat**

Many styles and brands
available

Funding is available
through the Assistive
Devices Program (ADP)





Type 2
(Indoor)

Unique material
ultimate comfort seat
Ergonomic handgrips
with optimal wrist
positioning
Available in three
colours – Black,
Blue, and Red



Type 3
(Indoor/Outdoor)

The nexuS 3 has innovated the way we think about mobility. With the first-ever cross-folding frame, it set the standard for future rollator design. The nexuS can stand independently in the folded position even with the zippered bag on the rollator. The frame design allows the user ample space for an ergonomic walking gait, or if needed, a slightly folded frame to pass through narrow doorways in smaller areas. Because of this unique design, stability is not compromised, ensuring the nexuS is the most reliable and safe rollator on the market.

STELLA Rollator



The light weight, easy to use STELLA rollator includes features to help people with reduced mobility in their daily activities. Built on an aluminum frame, the STELLA weighs only 14.75 lb* (6,7 kg) in its standard configuration. The oversized soft material backrest offers enhanced comfort and stability when seated.

*Weight does not include optional equipment.

standard features


- Reinforced seat
- Four 8" (205 mm) casters with smooth urethane tires
- Locking system when folded
- Stronger brake cables (25% thicker)
- Only 14.75 lb (6,7 kg)
- Oversized height adjustable backrest strap
- Height-adjustable handles with quick release system

Nexus Rollator Accessories

A



B



A Soft Bag

B Hard Basket

C Clip on Tray

D Slow Down Brakes

C



E Cane Holder

F Curb Climber

G Waterproof Seat Cover

H Thick Seat Pad

D



E



F



G



H



Stella Rollator Accessories

Stella Rollators come with a cup holder and cane holder

Other accessories available are shown below.

The most popular reason our clients choose the Stella is for it's ability to be converted into a Transport Chair!



Transport bag



Slow-down brakes



Ovation 805

Light weight side frame folding aluminum walker complete with 8" casters, seat sling, folding cloth basket, extra large backrest & curb assist feature.

Comes with cane holder and basket

Specifications:

- Seat height: 19" | 21" | 23"
- Seat width: 19"
- Width between handle bars: 19"
 - Overall width: 25"
 - Overall length: 28"
 - Folded width: 12"
- Walker weight: 16 lbs
- Weight capacity: 300 lbs
- Color choice: Black, Blue, Red or Silver



Evolution
TECHNOLOGIES*

With the new padded seat cushion and wide back strap pad, users can sit in comfort and walk in style. The Trillium comes with an all-new braking system, allowing for effortless handle height adjustment and easy maintenance.



Type 2
(Indoor)



Type 3 (Indoor/Outdoor)

Features


- Removable Seat Cushion 1.5" thick, quilted & padded
- Contoured Backrest Extra-wide for added comfort
- Concealed-cable brakes Easy to adjust and maintain
- Collapsible Basket Removable & foldable
- Durable Wheels 8" PU tires
- Offset Rear Forks for wider gait

Trillium Rollator Accessories



A

A Cup Holder



B

B Slow Down Brakes

C *Cane Holder



C



D

D *Sunglasses Holder

E Accessory Bracket – required for all 'holders'*

F Tray



E

G Curb Climber

H *Phone Holder



F



G

H



Type 2 Options

Note: All specs are in inches & lbs

| Manufacturer | ADP CODE | Model | Cable Free? | Seat W (WBH) | Seat Ht (in) | Overall W (in) | Wt Cap (lbs) | Wt (lbs) | HH Min (in) | HH Max (in) |
|--------------|-----------|---------------------|-------------|--------------|--------------|----------------|--------------|----------|-------------|-------------|
| Parsons | MW2DG0010 | 4200 DX LOW | No | 17.5 | 21.75 | 23.5 | 250 | 12 | 32.5 | 37 |
| Parsons | MW2DG0010 | 4200 DX SUPER LOW | No | 16.5 | 18.75 | 22.5 | 250 | 12 | 30 | 34 |
| Parsons | MW2DG0005 | 4200 LX ONE SIZE | No | 17.5 | 21.75 | 23.5 | 250 | 11 | 32 | 36 |
| Evolution | MW2EV0025 | CHALLENGER GRANDE | No | 17 | 21 | 24 | 250 | 13 | 33 | 38 |
| Evolution | MW2EV0010 | CHALLENGER MINI | No | 16 | 18 | 22 | 250 | 13 | 30 | 35 |
| Human Care | AA0001810 | NEXUS ONE SUPER LOW | Yes | 17.3 | 17 | 23 | 353 | 17.9 | 28.4 | 31.5 |
| Human Care | AA0001810 | NEXUS ONE LOW | Yes | 17.3 | 20 | 23 | 353 | 14.4 | 29 | 35 |
| Human Care | AA0001810 | NEXUS ONE STANDARD | Yes | 17.3 | 22 | 23 | 353 | 17 | 31 | 38 |
| Evolution | AA0001306 | TRILLIUM 2 MINI | Yes | 18 | 18 | 25 | 300 | 13 | 29 | 32 |
| Evolution | AA0001306 | TRILLIUM 2 REGULAR | Yes | 18 | 21 | 25 | 300 | 13 | 31.5 | 34.5 |
| Evolution | AA0001306 | TRILLIUM 2 TALL | Yes | 18 | 24 | 25 | 300 | 15 | 34.5 | 37.5 |

Type 3 Options

All specs are in inches or lbs

*WBH = Width Between Handles

| Manufacturer | ADP CODE | Model | Cable Free? | Seat W (WBH) | Seat Ht (in) | Overall W (in) | Wt Cap (lbs) | Wt (lbs) | HH Min (in) | HH Max (in) |
|--------------|-----------------|------------------------|--------------|--------------|--------------|----------------|--------------|----------|-------------|-------------|
| Parsons | MW3DG0035 + AA8 | 4200 SWX ONE SIZE | No | 24 | 22 | 27 | 440 | 27 | 30.5 | 36.5 |
| Parsons | MW3DG0030 | 4200 WX LOW | No | 17.5 | 22.5 | 24.5 | 250 | 14.5 | 34.5 | 38.5 |
| Parsons | MW3DG0030 | 4200 WX SUPER LOW | No | 17.5 | 19 | 24.5 | 250 | 14.5 | 30.5 | 35 |
| Parsons | MW3DG0030 + AA8 | 4200 WX WIDE LOW | No | 20.5 | 22.5 | 27 | 350 | 15.5 | 34.5 | 38.5 |
| Parsons | MW3DG0030 + AA8 | 4200 WX WIDE SUPER LOW | No | 20.5 | 19 | 27 | 350 | 15.5 | 34.5 | 38.5 |
| Evolution | MW3EV0035 + AA8 | CHALLENGER LOW WIDE | No | 21 | 19 | 26 | 400 | 16 | 33 | 37 |
| Evolution | MW3EV0030 + AA8 | CHALLENGER WIDE | No | 21 | 21 | 26 | 400 | 16 | 33 | 37 |
| Triumph | AA0001084 | ESCAPE LOW | Yes | 18.5 | 21 | 25 | 275 | 15.6 | 29.5 | 37.5 |
| Triumph | AA0001084 | ESCAPE STANDARD | Yes | 18.5 | 24 | 25 | 275 | 15.8 | 32 | 40 |
| Triumph | AA0001084 | ESCAPE SUPER LOW | Yes | 18.5 | 19 | 25 | 275 | 15.4 | 28 | 36 |
| Drive | MW3DGM0065 | EXCURSION HEMI X18 | No | 18 | 18 | 25 | 300 | 17 | 29.5 | 35 |
| Drive | MW3DGM0065 | EXCURSION REGULAR X20 | No | 18 | 20 | 25 | 300 | 17 | 30.5 | 37 |
| Drive | MW3DGM0065 | EXCURSION TALL X23 | No | 18 | 23 | 25 | 300 | 17 | 33.5 | 39.5 |
| Drive | AA0000124 + AA8 | EXCURSION XWD | No | 19.5 | 20 | 27 | 400 | 18 | 30 | 36.5 |
| Drive | AA0000694 + AA8 | FUSION 2-IN-1 HEMI X18 | No | 19.5 | 18 | 27 | 300 | 17 | 29.5 | 35 |
| Drive | AA0000694 + AA8 | FUSION 2-IN-1 REG X20 | No | 19.5 | 20 | 27 | 300 | 17 | 30 | 36.5 |
| Drive | AA0000694 + AA8 | FUSION 2-IN-1 TALL X23 | No | 19.5 | 23 | 27 | 300 | 19 | 33.5 | 39.5 |
| Sunrise | AA0001082 | GEMINO 30 | No | 18.5 | 24 | 23.6 | 330 | 15.5 | 30.7 | 39.4 |
| Sunrise | AA0001082 | GEMINO 30 M | No | 18.5 | 21 | 23.6 | 285 | 15.5 | 27.16 | 34.6 |
| Sunrise | AA0001082 | GEMINO 30 S | No | 16.5 | 18 | 21.6 | 275 | 14.5 | 25.6 | 30.3 |
| Sunrise | AA0001236 + AA8 | GEMINO 60 | No | 20 | 24 | 25 | 330 | 18.7 | 29.9 | 37.4 |
| Sunrise | AA0001236 + AA8 | GEMINO 60 M | No | 20 | 21 | 25 | 330 | 18 | 28.7 | 32.6 |
| Human Care | MW3DG0045 | NEXUS 3 LOW | Yes | 18 | 21 | 24 | 441 | 15 | 30.5 | 36.5 |
| Human Care | MW3DG0045 | NEXUS 3 STANDARD | Yes | 18 | 24 | 24 | 441 | 16 | 32.5 | 39.5 |
| Human Care | MW3DG0045 | NEXUS 3 SUPER LOW | Yes | 18 | 18.5 | 24 | 441 | 15 | 30.5 | 33.5 |
| Drive | AA0001238 + AA8 | NITRO HEAVY DUTY | Yes | 23 | 23.5 | 28.75 | 450 | 20.5 | 36 | 41 |
| Drive | AA0001193 | NITRO HEMI | Yes | 18 | 18 | 27.75 | 300 | 20.5 | 28 | 31 |
| Drive | AA0001036 | NITRO REGULAR | Yes | 18 | 20.5 | 23 | 300 | 17.5 | 33.5 | 38.25 |
| Drive | AA0001194 | NITRO TALL | Yes | 18 | 23.6 | 23 | 300 | 17.5 | 36 | 41 |
| Nighthawk | MW3NM1005 | OVATION 805 HIGH-RIDER | Yes | 18 | 23 | 25 | 300 | 16 | 35 3/4 | 39 3/4 |
| Nighthawk | MW3NM1005 | OVATION 805 LOW-RIDER | Yes | 18 | 19 | 25 | 300 | 16 | 33.5 | 37.5 |
| Nighthawk | MW3NM1005 | OVATION 805 MID-RIDER | Yes | 18 | 21 | 25 | 300 | 16 | 33.5 | 37.5 |
| Nighthawk | MW3NM1010 | OVATION 806 HIGH-RIDER | Yes | 18 | 23 | 26 | 300 | 16 | 35 1/4 | 39 1/4 |
| Nighthawk | MW3NM1010 | OVATION 806 MID-RIDER | Yes | 18 | 21 | 26 | 300 | 16 | 33.5 | 37.5 |
| Amyllor | AA0001580 | STELLA 18 | Yes | 17.5 | 18 | 24 | 300 | 14.5 | 29.5 | 33.5 |
| Amyllor | AA0001580 | STELLA 20 | Yes | 17.5 | 20 | 24 | 300 | 14.5 | 30.5 | 34.5 |
| Amyllor | AA0001580 + AA8 | STELLA 20 HD | Yes | 19.5 | 20 | 26 | 400 | 14.5 | 30.5 | 34.25 |
| Amyllor | AA0001580 | STELLA 23 | Yes | 17.5 | 23 | 24 | 300 | 14.5 | 33.5 | 37.5 |
| Amyllor | AA0001580 + AA8 | STELLA 23 HD | Yes | 19.5" | 23" | 26" | 400 | 14.5 | 33.5 | 38.5 |
| Evolution | AA0001307 | TRILLIUM 3 MINI | Yes | 18 | 18 | 26 | 300 | 14 | 29 | 32 |
| Evolution | AA0001307 | TRILLIUM 3 REGULAR | Yes | 18 | 21 | 26 | 300 | 14 | 31.5 | 34.5 |
| Evolution | AA0001307 | TRILLIUM 3 TALL | Yes | 18 | 24 | 26 | 350 | 14 | 34.5 | 37.5 |
| Evolution | MW3EV0040 | XPRESSO LITE CF MINI | Either Avail | 18 | 18 | 22.5 | 300 | 15 | 30.5 | 35.5 |
| Evolution | MW3EV0040 | XPRESSO LITE CF REG | Either Avail | 18 | 21 | 22.5 | 300 | 15 | 30.5 | 35.5 |
| Evolution | MW3EV0040 | XPRESSO LITE CF TALL | Either Avail | 18 | 24 | 22.5 | 350 | 15 | 33.5 | 40 |
| Evolution | MW3EV0040 | XPRESSO MINI | Yes | 18 | 18 | 22.5 | 300 | 18 | 30.5 | 33.5 |
| Evolution | MW3EV0040 | XPRESSO REGULAR | Yes | 18 | 21 | 22.5 | 300 | 18 | 30.5 | 34.5 |
| Evolution | MW3EV0040 | XPRESSO TALL | Yes | 18 | 24 | 22.5 | 350 | 18 | 33.5 | 40 |
| Evolution | AA0001283 + AA8 | XPRESSO WIDE | Yes | 20 | 22 | 24.8 | 450 | 31 | 32.5 | 37.5 |
| Evolution | AA0001283 + AA8 | XPRESSO WIDE LOW | Yes | 20 | 19 | 24.8 | 450 | 30 | 30.5 | 35.5 |
| Evolution | AA0001282 | XPRESSO ZERO MINI | Yes | 17.5 | 18 | 26 | 300 | 16 | 28.5 | 31.5 |
| Evolution | AA0001282 | XPRESSO ZERO REG | Yes | 17.5 | 21 | 26 | 300 | 16 | 31 | 35 |
| Evolution | AA0001282 | XPRESSO ZERO TALL | Yes | 17.5 | 24 | 27 | 350 | 17 | 33 | 37 |
| Human Care | AA0001895 | Crossover | Yes | 18 | 21 | 26 | 300 | 17 | 33.8 | 38.5 |

Aluminum Folding Walkers



Aluminum Folding Walkers are lightweight folding walkers that are available with no wheels or two front wheels.

Additional accessories such as rear glides are available for ease of use on carpets and other terrains. Aluminum folding walkers do not have brakes or a seat.



Walker Accessories



A



B

A Over handle
Tray

B Front flip down
Tray

C Wire basket

D Zippered Bags

E, F Glide Skis

G Tennis Ball Glide

i Swivel Locking

J 5" Wheels

K 3" Swivel Locking

H Bariatric Wheels



C



D



E



F




G




H



i



J



K

Transport Chairs

Transport chairs are meant for exactly that, transport!

These lightweight folding chairs are ideal for a quick stroll in the park, attending doctors appointments, etc. Ideal for putting in your vehicle for travel.



Basic Transport Chair

- Lightweight aluminum construction (only 20 lbs)
- Flared handles and handy back-release hinge allow the back to fold down flat
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Push-to-lock rear wheel locks
- Composite 8" wheels with solid rubber tires are lightweight and maintenance free
- 300 lbs Weight Capacity



Transport Chair with Handle Brakes

- Weighs only 19 lbs
- Available in blue or red frames
- Folds flat for easy transport
- Handy lever on chair back makes it easy and convenient to fold back down
- Comes with seat belt
- Companion activated wheel locks
- 12" rear "flat-free" wheels
- 300 lbs Weight Capacity



CGT Deluxe Transport Chair

- Dynamic brakes
- 12" back wheels
- Quick-release, fold-down backrest
- Flip-up, removable, and dual-function armrests (full length and desk length)
- Wheel locking mechanism
- Swing-away and length adjustable footrests
- Available Options:**
 - Width 17", 19" or 22"
 - Burgundy or blue



Luxury Travel Mobility



**Experience newfound freedom
with a folding power chair**

Buy with confidence at Comfort Mobility
Home Health Care

Supplier of only trusted products backed by
reliable service and warranty, you can count on
us for quality equipment for competitive prices!

Call, Stop in or visit our
website for more
information!



Jazzy Carbon

300 lbs Weight Capacity
43 lbs, 39 without battery

Limited Suspension

Lithium Ion Battery

24" Turn Radius

Up to 6km/hr Speeds

Seat Size 17x17

8.5" Rear Wheels

7" Caster Wheels

24V 12Ah Battery Size

Warranty: Limited Lifetime on Frame,
1 year limited electronics and batteries,
1 year drive motors

Jazzy Ultra Lite

300 lbs Weight Capacity
38-41 lbs, 33 without battery

Suspension Caster Forks

Lithium Ion Battery

24" Turn Radius

Up to 6km/hr Speeds

Seat Size 18x18

10.71" Rear Wheels

8" Caster Wheels

Battery Size Options: 10Ah and 20Ah

Warranty: Limited Lifetime on Frame, 1
year drivetrain and electronics, 13 months
on batteries



AS
LIGHT AS
**52.6
lbs.!***



 **GOLDEN**

 **GOLDEN**
Stride 

Standard Features:

- Aluminum Frame
- 300 lb. Weight Capacity
- Removable Compact Joystick with USB Charger
- 50 Amp Controller
- Comfort Spring Suspension
- Removable Padded Seat
- Convenient Under-Seat Storage & Back Seat Pocket
- Offboard Charging Capability
- Airline-approved Lithium-ion Battery
- Easy One-Step Folding
- Stands on End for Storage

*Excludes battery, joystick, cushion and bag.



 **GOLDEN**
CRICKET

Standard Features:

- 100% Carbon Fiber Frame
- 300 lb. Weight Capacity
- Removable Compact Joystick with USB Charger
- 35 Amp Controller
- Convenient Under Seat Storage Bag
- Removable Padded Seat
- Offboard Charging Capability
- Easy One-Step Folding
- Stands on End for Space-saving Storage
- Airline-approved Lithium-ion Battery

* Excludes battery, joystick, cushion and bag.



**ONLY
33 lbs.***
to carry
but holds up to
300 lbs.!





Aerolux Carbon Travel Buggy

29 lbs, 32 lbs with battery
6 Km/Hr Max Speed
18km per charge
Travel Safe lithium Ion Battery



Dash Travel Buggy

43 lbs, 46 lbs with battery
7 Km/Hr Max Speed
15km per charge
Travel Safe lithium Ion Battery



City 2 Plus Travel Buggy

Reclining Seat, Adjustable
Footplate and Removable Rear
Wheels for ease of travel

40 lbs, 65 lbs with battery and
Rear Wheels
7 Km/Hr Max Speed
25km per charge
Travel Safe lithium Ion Battery



Did you know that cruise ships are now enforcing new rules about mobility scooters? As of March 7th, 2024, they've updated their rules around scooter storage.

Now, scooters must be stored inside the stateroom and cannot be stored in corridors or public areas. This might have you worried but don't be. We've got you covered! The BuzzAround Carry On Power Chair is "Cruise Approved", as it was designed for compact storage either under a bed or in a closet.

- 300lb weight capacity
- Center Wheel Drive for Optimal Control & Maneuverability
- 21" Turning Radius
- Independent 4-Wheel Suspension
- Springless Articulating Chassis
- USB Port
- Front & Rear Lights
- Width Adjustable Armrests
- Textured Soft Touch Armrests
- Angle Adjustable Foot Plate
- Easy Access, Airline-Friendly, 12.5AH Lithium-Ion Battery with Travel Case
- PowerLink Connectivity for Easy Assembly/Disassembly
- Removeable Storage Basket Under Seat
- Easily Interchangeable Colour Shrouds (Red installed/Blue in box)





The Go Chair[®] MED is #1 for portability! With convenient feather-touch disassembly, the Go Chair MED allows you to experience lightweight travel and independence on the go. Compact, portable and perfect for travel, the Go Chair[®] MED is a feature-rich power chair that disassembles easily in just a few steps! A large footplate gives you the comfort you need.

Go Chair[®] MED Features:

- 300 lb. weight capacity
- Speeds up to 3.7 mph
- Compact, 25.5" Turning radius
- Intuitive controller and easy to operate
- Range of up to 9.41 miles per charge
- Compact, foldable black vinyl seat with 5-point swivel seat with latch release available in 18" x 17" or 20" x 18" and 60" lap belt
- Compact turning radius
- Large footplate for added comfort
- One hand feather-touch disassembly for quick transport and storage
- The heaviest piece is only 36 lbs.



The Robooter X40 Power Wheelchair


features a sleek, modern design with clean lines and a streamlined silhouette. Available in an elegant white finish, it combines a fine matte surface with a lightweight yet durable magnesium-aluminum alloy frame. The ergonomic design includes a 2.75" breathable memory foam cushion and a 15-degree backrest tilt for superior comfort. This stylish and functional design ensures that the wheelchair is both aesthetically pleasing and highly practical for everyday use.



Features:

Hands-free automatic folding and unfolding
Advanced dual suspension system for a smoother ride

Remote control via a user-friendly mobile app
Removable battery for convenient charging
Ergonomically designed backrest for maximum comfort



Omnidirectional wheel for seamless navigation on diverse surfaces



The Robooter E40 is a sleek, foldable indoor electric wheelchair that seamlessly combines performance, portability, and style. Designed for easy daily use with comfortable leg space and a small turning radius, it is powered by a robust 20AH lithium-ion battery, covering an impressive 14.3 miles on a single charge. With smooth maneuverability on various floor surfaces, its lightweight design makes it ideal for transportation in cars, airplanes, and more.

Features:

Easy to Fold
Lightweight Design (55lbs)
Small Turning Radius
Adjustable Armrest
Removable Battery for convenient charging



Patient Lifts

The Hoyer Advance® is compact, yet able to perform an outstanding range of transfers. Whether from the floor or onto a healthcare bed, the Advance copes with ease. The Advance's unique tapered leg design allows the lift to get even closer to the widest of obstacles. This is especially useful when encountering large chairs, bulky commodes, and wheelchairs. The tapered design also ensures that the resident feels safe and secure during the transfer.

Span America F450T standing transfer aid is designed to make transport quick and simple while encouraging weight bearing as well as maintenance of trunk and upper body strength. F450T stand-up patient lift is the ideal solution for patients with the ability to participate in their own transfers. F450T standing transfer aid has swing-away padded seat for secure, comfortable transport and the underbed clearance allows use with most low beds. The dual, contoured shin pads offers added stability and non-slip footplate provides stable base for transfer.

FLOOR LIFTS STAND-AID LIFTS SLINGS

Also Available:
Porch Lifts
Stair Lifts
Lift Chairs



Usage requirements: The F450T is designed to provide transfer assistance to patients with the necessary strength to Bear weight in the standing position, Pull themselves up into standing position, Lower themselves safely into the seated position For patients who meets above criteria, use of the F450T encourages activity and helps maintain strength.

Joerns®
HEALTHCARE

SPAN America
Innovative Solutions.

Contact your Mobility Sales Consultant for assistance choosing the correct lifting device and/or sling. At Comfort, we only sell brands you can trust!

Stair Lifts

How a Stair Lift Works

Stair lifts employ a battery-powered rack and pinion drive system, constantly charged by a nearby outlet at the top or bottom of the stairs. This design allows the stair lift to carry its rated load up and downstairs. Safety features, like the over-speed governor and emergency braking system, ensure that the stair lift will never travel faster than it is supposed to. Typically, stair lifts travel at about 18 ft per minute.

There are three basic stair lift types: straight, curved and outdoor. Each model is suited to a different type of staircase and comes with a unique set of attributes and benefits.

Contact John, our stair lift expert for more information and pricing!





WHO CAN BENEFIT FROM A STAIRLIFT?

While the list is truly endless, here is an idea of our common clientele

- Seniors
- Those with limited mobility
- Those living with chronic conditions or recovering from injury affecting mobility

[LEARN MORE](#)

 **CONTACT US:**
519-988-1234



Founded in 1867, Stannah Stairlifts is a fifth-generation family business, where customer service, quality and safety come first. We have been challenging and setting standards in stairlift design for more than 40 years. With over 750,000 satisfied customers, Stannah is a name you can trust!

Jonathan Fase
Trusted Service
Technician and
Professional Installer



When to Consider a Stairlift

If stairs are becoming a problem, consider getting in touch with us. You shouldn't have to limit yourself to the first floor or plan your life around a few difficult trips up and down the stairs. Our dedicated team of experienced professionals can help you regain the full use of your home, so you can carry on enjoying life in the place you know and love.

Take the first step today. We're standing by and ready to help. Let us share our story with you and show you why, after so many years, Stannah continues to stand out among the world's most trusted stairlift brands.



Standard Features You'll Love

Stannah's exceptional warranty and reliability will bring you peace of mind. Our thoughtfully designed key features make your Stannah stairlift safe and easy to use. Features below apply to all indoor models unless otherwise specified.

Easy-to-Use Controls

The light-touch control allows you to operate the stairlift with any part of your hand.

A Chair that Swivels

Turn to a safe position at the top or bottom of your stairs, with manual swivel levers, or optional one or two-way powered swivel.

Armrest Safety

For your safety, the downside armrest must be lowered before riding.

Safety Sensors

Stop the stairlift if it meets an obstruction.

350lb Weight Capacity

Keyed Lock (not shown)

Prevents unauthorized use.

Powered Footrest

For Starla and Siena 600, eliminates the need to bend over to unfold or fold the footrest.



Retractable Seatbelt

Easy to reach retractable seatbelt, requires only one hand to use (Starla & Siena 600).

Siena 600 Chair Width Options

Available in two width options depending upon your needs.

On/Off Switch

Conveniently placed, prevents the chair from moving unexpectedly.

Adjustable Seat Height

Set during installation for a comfortable and tailored fit.

Battery Operated

Automatically recharging, long-life batteries (not shown) keep your stairlift working in the event of a power outage.

Two Remotes

Call and send your stairlift or operate the powered swivel.

Power-On Light (not shown)

Options



Powered Swivel

The Stannah powered swivel turns the chair to face away from the stairs, for ease-of-use and safety. The powered swivel can be set to swivel at both the top and bottom of the stairs. Stannah offers the best installation solution when there is a door at the top of the stairs.



Retractable Rail

For straight stairs where the rail could be a tripping hazard at the bottom of the stairs, the Stannah retractable rail raises and lowers automatically as the chair moves up and down the stairs.

Seatbelt Options

Our lap-diagonal and five-point seatbelts provide extra restraint over the shoulders and between the legs.



Rudder Control

The rudder control option is helpful for those who may have difficulty with dexterity, for example due to severe arthritis or tremors.



Stay safe & independent at home with a Stannah Stairlift

Working with Stannah, Comfort Mobility can offer one of the world's most popular stairlift brands.

If you need a safe and reliable stairlift to help make your home safer, get in touch with Comfort Mobility today, and stay in the home you love!



Comfort Mobility
Home Health Care

\$1,000 CREDIT FROM STANNAH
FOR VETERAN'S WHEN YOU
PURCHASE A CURVED UNIT

Stannah

LIFETIME
WARRANTY ON
MOTOR &
GEARBOX
27 MONTHS ON
ALL PARTS



comfortmobility.ca



Call Jonathan for a quote
519-988-1234



519-988-1234

Rehabilitation Equipment - Sales & Rentals

KNEE REPLACEMENT?

Knee Immobilizers
Cold Therapy Machines
Crutches and Walkers
Compression Stockings
Knee Braces
Bathroom Safety Equipment
'ROM or Slider' Boards
& More!

We direct bill to Green Shield!

2707 Temple Drive, Windsor
info@comfortmobility.ca

comfortmobility.ca



VISIT OUR NEW & IMPROVED WEBSITE

Learn about new product, find out what funding is available to you
shop our online store
Stay up to date with the newest equipment and latest deals


Free Local Delivery
when you spend \$85



comfortmobility.ca



519-988-1234



Follow Us on Social Media

Let's Get Connected for Our Latest News & Updates
On our Facebook page, we do giveaways and announce
our latest sale flyers



on X (formerly Twitter) @CMWindsor



on Instagram @comfortmobility519



on Facebook @ComfortMobilityWindsor



on Youtube @comfortmobilityinc.3711



COMFORT MOBILITY

EST 2014

Your Mobility Journey Begins Here

Visit us online at
www.comfortmobility.ca

