

# PrevaMatt® Series Therapeutic Support Surfaces



Revolutionizing Care Together

## Superior Technology = Superior Performance

The Joerns RecoverCare PrevaMatt® series has undergone a complete transformation to create a portfolio of foam mattresses designed for comfort and safety all while providing pressure ulcer prevention. The expanded portfolio gives healthcare Providers a range of options to improve outcomes, reduce risks to both the caregiver and patient, and meet financial objectives.

- **IMPROVE CLINICAL OUTCOMES**  
Excellent pressure redistribution aids in wound prevention
- **REDUCE RISK OF INJURY**  
Standard safety features reduce fall risk and promote mobility
- **REDUCE FINANCIAL BURDEN**  
Broad and diverse portfolio designed for quality and durability



GREENER & MORE SUSTAINABLE PORTFOLIO

### Variable Pressure Foaming (VPF™) process

- Provides greatest consistency and highest quality mattress
- Virtually zero-emissions manufacturing process that exceeds Clean Air Act requirements

### Sustainability

- Plant-based raw materials
- Elimination of unnecessary chemicals results in purer material and can lead to improved durability





**Open cell memory foam technology**

- Allows greater air flow resulting in more breathability and a cooler, drier patient experience
- Materials less affected by fatigue – better height and support retention with faster recovery time

**MANIPULATION OF FOAM PROPERTIES**

**Surface Modification Technology (SMT™)**

- Alterations in pattern, size of pattern, spacing between patterns and depth of cut to protect high risk areas
- Zoned pressure redistribution with improved air circulation and maximized comfort



**TEMPUR® Option**



Partnership with Tempur-Pedic® to provide genuine TEMPUR® material on the PrevaMatt Alleviate

- Reduces wound incidence rates
- Treats up to Stage III uncomplicated pressure ulcers



**SAFETY AND SECURITY**

**Firm Perimeter**




























- Provides edge of bed stability
- Aids ingress and egress

**Optional Raised Perimeter**

- Aids in fall prevention
- Reduces risk of entrapment

When used as part of a Joerns bed system, our PrevaMatt mattresses are guaranteed to meet the guidelines set forth in the Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment developed by the Food and Drug Administration (FDA) and the Hospital Bed Safety Workgroup (HBSW).

**PrevaMatt® Series of Prevention Surfaces**

	Basic	Ease	Defend	Console	Alleviate	Enmerse
VPF™						
SMT®						
Natural Ingredients						
Cover Material (Nylon, 2-way, 4-way stretch)						
Open Cell Visco						
TEMPUR® Option						
2nd Layer SMT®						
Clinical Indications	Prevention	Prevention	Prevention to Stage I	Prevention to Stage I	Prevention to Stage II <sup>†</sup>	Prevention to Stage II

<sup>†</sup>TEMPUR® option is suitable for treatment of up to Stage III uncomplicated pressure ulcers

\*VPF is a trademark and SMT is a registered trademark of FXI®

\*\*TEMPUR is a registered trademark of Tempur-Pedic®