




Global solutions Matachana  
for pathology and morgue  
**Anatomical Pathology**

# Autopsy rooms, morgue & laboratory equipment

MATACHANA carries out complete equipment projects for autopsy rooms in hospitals, funeral parlours and forensic medicine buildings.

We apply our wide experience to every project, implementing comprehensive solutions for autopsy rooms and morgues, areas which must be completely aseptic and hygienic at all times and kept strictly controlled and organized.

Our portfolio includes a full line of products for equipping different areas such as dissecting, necropsy and mortuary, as well as lab furniture and transport carts.



- |   |   |    |  |
|---|---|----|--|
| 1 | Lifting and transport systems             | 6  | Lighting and video systems   |
| 2 | Walk-in refrigerator systems              | 7  | Autopsy tables   |
| 3 | Mortuary refrigeration and freezing units | 8  | Grossing stations /Cutting tables                                    |
| 4 | Corpse washing area                       | 9  | Stainless steel furniture  |
| 5 | Forensic anthropology                     | 10 | Equipment for washing, disinfection and sterilization of instruments |



## Heat sterilizer system for bio-contaminated effluents

In addition, as manufacturers of steam sterilizers, we provide devices and solutions for the treatment of bio-contaminated waste effluents produced in all those areas. Unlike the chemical disinfection equipment, chemical residues are not generated by the steam sterilization. Therefore, it does not arise environmental problems. It is simple, clean (only steam is used) quick and effective, as it eliminates all pathogens capable of generating human infections.



matachana


## Mortuary refrigerators and deep freezing units

MATACHANA provides the best solutions for corpse storage and preservation needs with our Refrigeration Units and Walk-in Refrigerator Systems. We also offer the possibility to execute customized projects with top quality materials.





Complete range in preservation/  
freezing chambers for standard  
and obese bodies, with all  
possible configurations.  
Capacities from 2 to 20 bodies,  
and possibility of a tailored-  
made study.



## Equipment for autopsy rooms

MATACHANA designs, projects, manufactures and equips any facility for the most advanced and fully equipped mortuary, laboratory or autopsy room.


Our long experience handling and working stainless steel for over 50 years, allows us to provide a wide range of dissection and cutting tables, laboratory furniture and other elements entirely manufactured in stainless steel.

## Autopsy tables

MATACHANA dissection tables can be configured under the following variations: fixed & adjustable height, rotating with fixed and adjustable height.

Each device can be customized with the following list of options: aspiration, waster, mixer for hot/cold water, extensible showers, electric sockets, water exhauster, irrigation systems; and accessories, including instrument trays, organ table, body holders, neckrest holders, scales, etc.





Autopsy rooms  
**Anatomical Pathology**





## MATACHANA GROUP

# Innovation is the way forward

Since our foundation more than 50 years ago, our mission has been to provide the best service, bringing our knowledge and field experience to our customers to facilitate their daily work, allowing them to be efficient in the production whereas keeping rigorously the quality.

Our history has taken us from the beginnings in 1962 to the transformation in 2000 in MATACHANA GROUP, with its subsidiaries in Germany, France, USA, Argentina, Malaysia and China, as well as the collaboration of direct distributors in over 110 countries.

In MATACHANA GROUP we are aware that Training and Service contribute decisively to achieve customer satisfaction. For this reason, we have always invested in the development of these two areas of activity, which enables a direct contact with customers and help us to develop together a continuous improvement process.

### Technical Engineering Support

An assistance provided by engineers, highly skilled expert technicians and support staff, all committed to ensure the proper equipment operation and condition.

### Training Academy



Committed to the responsibility to obtain the optimum operation from our equipment, in MATACHANA GROUP we invest in the training of future users.

In addition, we offer training courses regularly in the 5 continents, providing this service to all our customers.

### Environmentally friendly

The sterilizers are designed and manufactured using the latest technologies on the market to achieve the best results in terms of energy savings and reduction in water consumption.

Our Production Center complies with the ISO 14001 Environmental Management System and ISO 50001 Energy Management.

### Quality

The MATACHANA devices are developed, manufactured and tested within a strict quality control according to the International Standard ISO 9001.

www.matachana.com

CA\_AP\_01\_00\_en 0617

All specifications herein are provided for information purposes only and may be modified without notice.

 **matachana** | **+50** YEARS | Experience that improves lives

**MATACHANA GROUP**  
Headquarters · Almogàvers, 174  
E-08018 Barcelona

**Spain**  
Tel. (+34) 93 486 87 00  
hospitales@matachana.com

**International**  
Tel. (+34) 93 300 85 04  
int@matachana.com

**Spain · Germany · France · USA · Argentina · Malaysia · China**