



Hemraj
INDIA

LAB 2 ADVANCED LABORATORY SPRAY DRYER

Capacity : 3 - 5 kg / hr water evaporation

**SS 316L construction; High purity design; Ideal for low volumes
Various pharmaceutical grade finishes; Sanitary type ports are standard
Solvent and Aqueous feeds; HEPA filters for gas streams
Pressure and two fluid nozzle systems; Co - Current spray; High yields
Compact design ideal for laboratories and kilo labs.**

Recent Customers :

- 3M Corp., USA (2 Units)
- Misr Café Co., Egypt
- Taste Master Flavors, India
 - Lupin Ltd., India
 - Praj-Matrix, India
 - Givaudan, India
 - ITG, Israel
- Vedant Dyestuffs Intermediates Pvt. Ltd., India



Fully automated PLC based touch screen operator panel with recipe function and USB data logger. Safe, flexible and easy to handle.



Hemraj
INDIA

SS SMALL SCALE SPRAY DRYER

Capacity : 5 - 8 kg / hr water evaporation

SS 316L construction; Robust and high purity design per GMP standards; Various pharmaceutical grade finishes; Sanitary type ports are standard; Solvent and Aqueous Feeds with flameproof design; Uses of inert process gas or air as the drying medium; HEPA filters for gas streams; Ideal for establishing process data for commercial production; Easy to dismantle for cleaning and fast product switch.



Recent Customers :

- **Dr. Reddy's Laboratories Ltd., Vizag (3 units)**
- **Ranbaxy Laboratories Ltd., India**
- **Saint Gobain Industries, France**
- **Labotrade SL, Spain**
- **IFF Ltd., India**

Additional Features :

- **Simple and skid mounted arrangement.**
- **Available with rotary atomizer or pressure and two fluid nozzle systems.**
- **Available with auto / manual mode and PLC based system with touch screen operation panel, online trends and USB data logging of all parameters.**
- **CIP systems available.**