



**Hemraj**  
INDIA

## **SS SMALL SCALE SPRAY DRYER**

**Capacity : 5 - 10 kg / hr water evaporation**

**SS 316L construction; Robust and high purity design per GMP standards; Various pharmaceutical grade finishes; Sanitary type ports are standard; Solvent and Aqueous Feeds with flameproof design; Uses of inert process gas or air as the drying medium; HEPA filters for gas streams; Ideal for establishing process data for commercial production; Easy to dismantle for cleaning and fast product switch.**



### **Recent Customers :**

- Dr. Reddy's Laboratories Ltd., India
- Ranbaxy Laboratories Ltd., India
- Saint Gobain Industries, France
- Labotrade SL, Spain
- IFF Ltd., India

### **Additional Features :**

- Simple and skid mounted arrangement.
- Available with rotary atomizer or pressure and two fluid nozzle systems.
- Available with auto / manual mode and PLC based system with touch screen operation panel, online trends and USB data logging of all parameters.
- CIP systems available.