


Multi-Band Ligator

Indications

This device is used in conjunction with endoscopes for ligation of esophageal varices or ligation of anorectal hemorrhoids.

Features

- ※ The ligation ring is installed at the back, the carrier is made of super transparent material, and the ligation area is bright and clear
- ※ 0.2mm double wire release; wider field of view, more conducive to capturing tissue, releasing ligature ring
- ※ Multi-ring burst, one operation with multiple points, faster and more convenient
- ※ Medical rubber material has good biocompatibility and physical properties
- ※ Humanized handle design, more comfortable and convenient to operate
- ※ Both tactile and audible loop release cues
- ※ WARNING: Gives an eye-catching reminder when there is only one ligature ring left
- ※ The square cross-section collar design makes the ligation stronger and safer



Specifications

Code	Loop No.	Code	Loop No.
BL18-4A	4 bands	BL18-4M	4 bands
BL18-5A	5 bands	BL18-5M	5 bands
BL18-6A	6 bands	BL18-6M	6 bands
BL18-7A	7 bands	BL18-7M	7 bands
BL18-8A	8 bands	BL18-8M	8 bands
BL18-9A	9 bands	BL18-9M	9 bands