




**THE
RIGHT
CHOICE
FOR
RONGEURS**

CHOOSE VARIETY & PRECISION

From Straight to Curved, Micro to Spurling, Stainless to Coated, having the proper Rongeurs is critical to positive surgical outcomes. That is why we offer the largest selection of Rongeurs on the market, all of them expertly manufactured in Germany using the highest quality materials and state-of-the-art craftsmanship. All BOSS Rongeurs are backed by our Lifetime Warranty* against defects in materials and workmanship, and you can count on the support of our sales and customer service team for quotations, cross-references, order tracking, and post-sales service.

Our commitment to providing the widest variety of precision Rongeurs, combined with our best-in-class service, makes BOSS the clear choice for your Rongeur needs.

* Please refer to page 18 for warranty information.



ANATOMY OF A BOSS KERRISON RONGEUR

Micro Style

The most popular style of Kerrison Rongeur for both lumbar and cervical procedures. Micro Kerrison Rongeurs feature a 9mm bite opening and a handle reach 30% shorter than a Spurling.

APPROXIMATE DIMENSIONS:

- A** - Jaw Opening - 9 millimeters
- B** - Handle Width - 100 millimeters
- C** - Handle Height - 120 millimeters




Spurling Style

The original handle style for all Kerrison Rongeurs, Spurling Kerrison Rongeurs feature a 15mm bite opening and a long handle reach.

APPROXIMATE DIMENSIONS:

- A** - Jaw Opening - 15 millimeters
- B** - Handle Width - 140 millimeters
- C** - Handle Height - 145 millimeters



Tarpon™ Hybrid Style

Tarpon™ from BOSS offers the 15mm bite opening of a Spurling with a smaller handle size similar to the popular Micro.

APPROXIMATE DIMENSIONS:

- A** - Jaw Opening - 15 millimeters
- B** - Handle Width - 100 millimeters
- C** - Handle Height - 130 millimeters

PERFORMANCE COATINGS

All of our Kerrison Rongeurs are crafted in high-quality, hardened, German stainless steel. However, to maintain a sharp edge and smooth action over frequent usage, BOSS Kerrison Rongeurs are available with two performance coatings, ensuring your Rongeur is at peak performance.



PANTHER™

Using an Aluminum Titanium Nitride (AlTiN) coating bonded to the steel at high heat and pressure, Panther Rongeurs feature a low-glare black finish, smooth cutting action, and increased cutting edge hardness.



SUPER-COAT

Super-Coat utilizes the strength and lubricity of Titanium Nitride bonded to the instrument surface. This special coating gives the Rongeur its distinctive gold color.


APEX KERRISON RONGEURS

Clearly the clean choice.

What was once one of the most difficult to clean instruments is now one of the easiest. With the simple push of a button, the top rail of our APEX Kerrison Rongeur lifts up for easy inspection and cleaning. For quick reassembly, the top rail locks back into place with a simple squeeze of the handle.

APEX Kerrison Rongeurs are available in stainless steel or with a Panther™ performance coating.


APEX Kerrison Rongeurs allow you to save both time and expense associated with cleaning by using a simple, innovative design. Make the clearly clean choice and switch to BOSS APEX Kerrison Rongeurs.



MICRO



Shorter handle reach
Smaller bite opening




Pin ejector



Shown with
Panther™ coating



T-bar ejector



MICRO KERRISONS small handle

40°▲ 9mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
7" 175mm	Non-Ejector	70-0160	70-0230	70-0232	70-0234	70-0236	70-0238
	Pin Ejector	—	70-1700	70-1701	70-1702	70-1703	70-1704
8" 200mm	Non-Ejector	70-0162	70-0330	70-0332	70-0334	70-0336	70-0338
	Pin Ejector	—	70-0690	70-0691	70-0692	70-0693	70-0694
9" 225mm	Non-Ejector	70-0164	70-0352	70-0354	70-0356	70-0358	70-0360
	Pin Ejector	—	70-0730	70-0731	70-0732	70-0733	70-0734
10" 250mm	Non-Ejector	70-0166	70-0364	70-0366	70-0368	70-0370	70-0372
	Pin Ejector	—	70-0750	70-0751	70-0752	70-0753	70-0754
11" 275mm	Non-Ejector	70-0168	70-0376	70-0378	70-0380	70-0382	70-0384
	Pin Ejector	—	70-0770	70-0771	70-0772	70-0773	70-0774

APEX MICRO KERRISONS small handle

40°▲ 9mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
7" 175mm	T-bar Ejector	70-2780E	70-2782E	70-2784E	70-2785E	70-2786E	70-2787E
8" 200mm	T-bar Ejector	70-2800E	70-2802E	70-2804E	70-2805E	70-2806E	70-2807E
9" 225mm	T-bar Ejector	70-2820E	70-2822E	70-2824E	70-2825E	70-2826E	70-2827E
11" 275mm	T-bar Ejector	70-2860E	70-2862E	70-2864E	70-2865E	70-2866E	70-2867E

* Additional sizes may be available as made-to-order. Extended lead times required.


For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

 **P = Panther™** |  **SC = Super-Coat**

TARPON™



Shorter handle reach
Larger bite opening




Pin ejector



Shown with
Panther™ coating



T-bar ejector



TARPON™ KERRISONS

midsize handle

40°▲ 15mm

Length	Style	1mm	2mm	3mm	4mm	5mm	6mm
7" 175mm	Non-Ejector	70-1831	70-1832	70-1833	70-1834	70-1835	70-1836
	Pin Ejector	–	70-1842	70-1843	70-1844	70-1845	70-1846
8" 200mm	Non-Ejector	70-1851	70-1852	70-1853	70-1854	70-1855	70-1856
	Pin Ejector	–	70-1862	70-1863	70-1864	70-1865	70-1866
9" 225mm	Non-Ejector	70-1951	70-1952	70-1953	70-1954	70-1955	70-1956
	Pin Ejector	–	70-1882	70-1883	70-1884	70-1885	70-1886
10" 250mm	Non-Ejector	70-1957	70-1958	70-1959	70-1960	70-1961	70-1950
	Pin Ejector	–	70-1970	70-1972	70-1974	70-1976	70-1978
11" 275mm	Non-Ejector	70-2151	70-2152	70-2153	70-2154	70-2155	70-2156
	Pin Ejector	–	70-1981	70-1982	70-1983	70-1984	70-1985

APEX TARPON™ KERRISONS

midsize handle

40°▲ 15mm

Length	Style	1mm	2mm	3mm	4mm	5mm	6mm
8" 200mm	T-bar Ejector	70-2910E	70-2912E	70-2914E	70-2915E	70-2916E	70-2917E
9" 225mm	T-bar Ejector	70-2920E	70-2922E	70-2924E	70-2925E	70-2926E	70-2927E
11" 275mm	T-bar Ejector	70-2910-11E	70-2912-11E	70-2914-11E	70-2915-11E	70-2916-11E	70-2917-11E

* Additional sizes may be available as made-to-order. Extended lead times required.

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

◆ P = Panther™ | ◆ SC = Super-Coat

SPURLING



Spurling


Longer handle reach
Larger bite opening



Shown with
Panther™ coating



T-bar ejector



SPURLING KERRISONS large handle

40°▲ 15mm

Length	Style	2mm 	3mm 	4mm 	5mm 	6mm 
7" 175mm	Non-Ejector	70-0530	70-0532	70-0534	70-0536	70-0538
8" 200mm	Non-Ejector	70-0630	70-0632	70-0634	70-0636	70-0638
9" 225mm	Non-Ejector	70-0640	70-0642	70-0644	70-0646	70-0648
10" 250mm	Non-Ejector	70-0650	70-0652	70-0654	70-0656	70-0658
11" 275mm	Non-Ejector	70-0660	70-0662	70-0664	70-0666	70-0668

APEX SPURLING KERRISONS large handle

40°▲ 15mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
8" 200mm	T-bar Ejector	70-2940E	70-2942E	70-2944E	70-2945E	70-2946E	70-2947E

* Additional sizes may be available as made-to-order. Extended lead times required.

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

 **P = Panther™** |  **SC = Super-Coat**



MICRO THIN FOOTPLATE

small handle

Length	Style	1.5mm	2mm	3mm	4mm	5mm
7" 175mm	Non-Ejector	70-0070	70-0071	70-0072	70-0073	70-0074
	Pin Ejector	—	70-0071E	70-0072E	70-0073E	70-0074E
8" 200mm	Non-Ejector	70-0080	70-0082	70-0084	70-0086	70-0088
	Pin Ejector	—	70-0082E	70-0084E	70-0086E	70-0088E

Made-to-order. Additional lead time may be required



Shown with
Panther™ coating

APEX MICRO THIN FOOTPLATE

small handle

Length	Style	1mm	2mm	3mm	4mm	5mm
7" 175mm	T-bar Ejector	70-2780TFE	70-2782TFE	70-2784TFE	70-2785TFE	70-2786TFE
8" 200mm	T-bar Ejector	70-2800TFE	70-2802TFE	70-2804TFE	70-2805TFE	70-2806TFE

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

 **P = Panther™** |  **SC = Super-Coat**



FORAMINOTOMY

streamlined hingeless design, delicate shaft curved up

Length	Style	1mm	2mm	3mm	4mm	5mm	6mm
8" 200mm	Standard	–	70-0950	70-0951	70-0952	70-0953	70-0954



FORAMINOTOMY

hinged design for extra strength, shaft available curved up, left and right

Length	Jaw Angle	Style	2mm	3mm	4mm	5mm	6mm
8" 200mm	40°	Curved Up	70-0682	70-0683	70-0684	70-0685	70-0686
	90°	Curved Left	–	70-0670	70-0671	70-0672	–
	90°	Curved Right	–	70-0674	70-0675	70-0676	–
10" 250mm	40°	Curved Up	70-0682L	70-0683L	70-0684L	–	–

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

◆ P = Panther™ | ◆ SC = Super-Coat



ROTATING KERRISONS small handle (sold separately as 70-2390)

Length	1mm	1.5mm	2mm	3mm	4mm	5mm	6mm
40° 15mm							
7" 175mm	70-2450	70-2451	70-2452	70-2454	70-2455	70-2456	70-2457
8" 200mm	70-2460	70-2461	70-2462	70-2464	70-2465	70-2466	70-2467
90° 15mm							
7" 175mm	70-2400	70-2401	70-2402	70-2404	70-2405	70-2406	70-2407
8" 200mm	70-2410	70-2411	70-2412	70-2414	70-2415	70-2416	70-2417

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

 **P = Panther™** |  **SC = Super-Coat**



BAYONET MICRO small handle

Length	Jaw Size and Angle	Style	1mm	2mm	3mm	4mm	5mm	6mm
6" 150mm	40°▲ 10mm	Non-Ejector	70-2251	70-2252	70-2253	70-2254	70-2255	70-2256
		Pin Ejector	–	70-2280	70-2281	70-2282	70-2283	70-2284
	90°▲ 10mm	Non-Ejector	70-2241	70-2242	70-2243	70-2244	70-2245	70-2246
		Pin Ejector	–	70-2286	70-2287	70-2288	70-2289	70-2290



Shown with
Panther™ coating

BAYONET SPURLING large handle

Length	Jaw Size and Angle	Style	1mm	2mm	3mm	4mm	5mm	6mm
6" 150mm	40°▲ 15mm	Non-Ejector	70-2320	70-2322	70-2323	70-2324	70-2325	70-2326
		Pin Ejector	–	70-2340	70-2342	70-2344	70-2346	70-2348
	90°▲ 15mm	Non-Ejector	70-2310	70-2312	70-2313	70-2314	70-2315	70-2316
		Pin Ejector	–	70-2350	70-2352	70-2354	70-2356	70-2358

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC).

◆ P = Panther™ | ◆ SC = Super-Coat

IVD RONGEURS



Shown with
Panther™ coating



Length	Style						
		Cushing 1.5mm	Cushing 2mm	Love- Gruenwald 3mm	Spurling 4mm	Spurling 5mm	Cloward 6mm
5" 125mm	Straight	–	70-0980	70-1030	70-1080	–	70-1250
	Up	–	70-0982	70-1032	70-1081	–	70-1251
	Down	–	70-0984	70-1034	70-1082	–	70-1252
6" 150mm	Straight	70-0967	70-0990	70-1040	70-1090	–	70-1430
	Up	70-0968	70-0992	70-1042	70-1091	–	70-1432
	Down	70-0969	70-0994	70-1044	70-1092	–	70-1434
7" 175mm	Straight	70-0970 ♦	70-1008 ♦	70-1053 ♦	70-1100 ♦	–	–
	Up	70-0971 ♦	70-1009 ♦	70-1054 ♦	70-1101 ♦	–	–
	Down	70-0972 ♦	70-1010 ♦	70-1055 ♦	70-1102 ♦	–	–
	Straight w/Teeth	–	70-1120 ♦	70-1123 ♦	70-1126	–	–
	Up w/Teeth	–	70-1121 ♦	70-1124 ♦	70-1127	–	–
	Down w/Teeth	–	70-1122 ♦	70-1125 ♦	70-1128	–	–
9" 225mm	Straight	70-0976	70-1022	70-1070 ♦	70-1222	70-1227	–
	Up	70-0977	70-1024	70-1072 ♦	70-1224	70-1228	–
	Down	70-0978	70-1026	70-1074 ♦	70-1226	70-1229	–
	Straight w/Teeth	70-1021	–	–	–	–	–
10" 250mm	Straight	–	70-1018	70-1066	70-1218	70-1236	–
	Up	–	70-1019	70-1067	70-1219	70-1237	–
	Down	–	70-1020	70-1068	70-1220	70-1238	–
	Straight w/Teeth	–	–	70-1069	70-1221	70-1239	–
11" 275mm	Straight	–	70-1027	70-1076	70-1230	70-1231	–
	Up	–	70-1028	70-1077	70-1232	70-1233	–
	Down	–	70-1029	70-1078	70-1234	70-1235	–
12" 300mm	Straight	–	70-1150	70-1155	70-1160	70-1165	–
	Up	–	70-1151	70-1156	70-1161	70-1166	–
	Down	–	70-1152	70-1157	70-1162	70-1167	–
	Straight w/Teeth	–	–	–	–	70-1239-12	–



Available with Panther™ performance coating. Please add "P" to the end of the item number (e.g. 70-1120P / 70-1350P).



DRAGON IVD RONGEURS

Length	Style	 2mm x 12mm	 3mm x 12mm	 4mm x 14mm	 5mm x 14mm	 6mm x 16mm
5½" 140mm	Straight	70-1281	70-1282	70-1283	–	–
	Up	–	70-1290	70-1291	–	–
	Down	–	70-1296	70-1297	–	–
6¼" 160mm	Straight	70-1285 ♦	70-1286 ♦	70-1287 ♦	70-1288 ♦	70-1289
	Up	–	70-1293	70-1294	70-1295	–
	Down	–	70-1298	70-1299	–	–
7¼" 185mm	Straight	70-1256 ♦	70-1257 ♦	70-1258 ♦	70-1259 ♦	70-1260
	Up	70-1262 ♦	70-1263 ♦	70-1264 ♦	70-1265	–
	Down	70-1272	70-1273	70-1274	–	–



Available with Panther™ performance coating. Please add "P" to the end of the item number (e.g. 70-1120P / 70-1350P).



FERRIS-SMITH IVD RONGEURS

						
Length	Style	2mm	3mm	4mm	5mm	6mm
7" 175mm	Straight	70-1350 ♦	70-1351 ♦	70-1352 ♦	–	70-1354
	Up	70-1360 ♦	70-1361 ♦	70-1362 ♦	–	70-1363
	Down	70-1364	70-1365	70-1366	–	–
10" 250mm	Straight	70-1540	70-1541	70-1542	70-1543	70-1544
	Up	70-1550	70-1551	70-1552	70-1553	70-1554
	Down	70-1555	70-1556	70-1557	70-1558	70-1559
12" 300mm	Straight	–	70-1562	70-1563	70-1564	70-1565
	Up	–	70-1567	70-1568	70-1569	70-1570
	Down	–	70-1572	70-1573	70-1574	70-1575

NOTE: Any Ferris-Smith IVD Rongeur can be ordered With Teeth by adding "WT" to the item number (e.g. 70-1575WT)



Available with Panther™ performance coating. Please add "P" to the end of the item number (e.g. 70-1120P / 70-1350P).



BAYONET IVD RONGEUR

Length	Style	2mm 	4mm 
6 1/2" 165mm	Straight	70-2750P	70-2754P
	Up	70-2760P	70-2764P



Only available with Panther™ performance coating.

BOSS INSTRUMENT ACCESSORIES

Visit bossinstruments.com for the full range of BOSS Instrument accessories.

Rongeur Cases

KERRISON RONGEUR CASES

All cases hold 8 Rongeurs

00-1220 - 10 x 12 x 6" (for up to 9" Rongeurs)

00-1230 - 10 x 16 x 6" (for up to 13" Rongeurs)

FORAMINOTOMY RONGEUR CASE


Extra height accommodates the curve of the foraminotomy Rongeur

00-1240 - 10 x 12 x 7.5" (extra height for curve)

IVD RONGEUR CASE

Holds up to 16 Rongeurs

00-1250 - 9 x 14 x 5"



00-1220

LIFETIME WARRANTY

MATERIALS AND DESCRIPTIONS

Unless otherwise stated, all instruments in this catalog are stainless steel. Instrument photographs accurately represent the instrument; however, they may not be actual size. Our instruments are hand-crafted, so measurements noted may vary slightly. We reserve the right to alter our patterns and make modifications to any of our products at any time and without notice.

WARRANTY

All BOSS™ instruments are guaranteed to be free of defects in material and workmanship. Any instrument failing as a result of such defect shall be replaced free of charge at the sole discretion of BOSS INSTRUMENTS, LTD. No other warranty or affirmation of fact expressed or implied is made or authorized. Immediate cancellation of this guarantee shall result from the alteration, modification or repair of any product by persons not authorized by BOSS INSTRUMENTS, LTD., or by the physical or chemical abuse of any product.

The buyer and/or user assumes full responsibility for determining the suitability of the application of all products herein listed and for all risks and liability associated with their use.



BOSS

INSTRUMENTS LTD

104 Sommerfield Drive
Gordonsville, VA 22942 USA

Toll Free: 800-210-2677
Phone: 540-832-5000
Fax: 540-832-5515
Email: inquiry@bossinstruments.com
bossinstruments.com

©2024 BOSS INSTRUMENTS, LTD. REV G L-34
BOSS AND BOSS INSTRUMENTS, LTD. ARE REGISTERED TRADEMARKS