




nxt™ non-stick bipolar forceps



Contents



06-29

nxt™ bipolar forceps
with nxt™ non-stick technology

nxt™ Bipolarpinzetten
mit nxt™ Non-Stick Technologie



30-33

nxt™ irrigation bipolar forceps
with nxt™ non-stick technology

nxt™ Spülpinzetten
mit nxt™ Non-Stick Technologie



34-39

nxt™ suction bipolar forceps

nxt™ Saugpinzetten



40-47

Lehmburg Keyhole™ bipolar forceps
with nxt™ non-stick technology


Lehmburg Keyhole™ Bipolarpinzetten
mit nxt™ Non-Stick Technologie



48-55

Round handle bipolar forceps
with nxt™ non-stick technology

Rundgriff Bipolarpinzetten
mit nxt™ Non-Stick Technologie



56-57

bipolar cables

Bipolarkabel

Superior Non-Stick Technology. Made in Germany.




adeor nxt™ non-stick bipolar forceps with instrument tips made of solid silver (not coating) are characterised by very good heat transfer. This prevents overheating of the tips during coagulation. Independent tests confirm their superior non-stick properties. The non-stick properties remain intact throughout the entire lifespan of the instruments.


adeor nxt™ Non-Stick Bipolarpinzetten mit Instrumentenspitzen aus solidem Silber (keine Beschichtung) zeichnen sich durch einen sehr guten Wärmetransfer aus. Dadurch wird eine Überhitzung der Pinzettenspitzen während der Koagulation verhindert. Unabhängige Tests bestätigen die überlegenen Antihafteigenschaften. Die Non-Stick Eigenschaften bleiben über die gesamte Lebensdauer der Pinzetten erhalten.

nxt™ Non-Stick Bipolar Forceps

vs.

Standard Bipolar Forceps






Optimal balance and weight. Fatigue free operation.

Optimal balance and weight, good tactile feedback,
fatigue free operation.


*Optimale Balance und Gewicht, gutes taktiler Feedback,
ermüdungsfreies Arbeiten.*



nano™ non-stick bipolar forceps



nano™ tips with 4 mm tip length



classic tines vs. nano™ tines
Ø 3.5 mm vs. Ø 2.0 mm




TruPar™ pin guide



Wide ergonomic handle



US pin or EU flat connectors



nano™ non-stick



Ø 0.2 x 4 mm nano™ tip straight

nano™ series / working / total length 55 / 170 mm

BIP001EU eu flat connector

BIP001US us pin connector



Ø 0.4 x 4 mm nano™ tip straight

nano™ series / working / total length 55 / 170 mm

BIP002EU eu flat connector

BIP002US us pin connector




Ø 0.7 x 4 mm nano™ tip straight

nano™ series / working / total length 55 / 170 mm

BIP003EU eu flat connector

BIP003US us pin connector



nano™ non-stick



Ø 0.2 x 4 mm nano™ tip straight
nano™ series / working / total length 85 / 200 mm
BIP004EU eu flat connector
BIP004US us pin connector




Ø 0.4 x 4 mm nano™ tip straight
nano™ series / working / total length 85 / 200 mm
BIP005EU eu flat connector
BIP005US us pin connector



Ø 0.7 x 4 mm nano™ tip straight
nano™ series / working / total length 85 / 200 mm
BIP006EU eu flat connector
BIP006US us pin connector



Ø 1.0 x 4 mm nano™ tip straight
nano™ series / working / total length 85 / 200 mm
BIP007EU eu flat connector
BIP007US us pin connector



nanoTM non-stick



Ø 0.2 x 4 mm nanoTM tip straight

nanoTM series / working / total length 115 / 230 mm

BIP008EU eu flat connector

BIP008US us pin connector



Ø 0.4 x 4 mm nanoTM tip straight

nanoTM series / working / total length 115 / 230 mm

BIP009EU eu flat connector

BIP009US us pin connector



Ø 0.7 x 4 mm nanoTM tip straight

nanoTM series / working / total length 115 / 230 mm

BIP010EU eu flat connector

BIP010US us pin connector




Ø 1.0 x 4 mm nanoTM tip straight

nanoTM series / working / total length 115 / 230 mm

BIP011EU eu flat connector

BIP011US us pin connector



nano™ non-stick



Ø 0.2 x 4 mm nano™ tip straight

nano™ series / working / total length 135 / 250 mm

BIP012EU eu flat connector

BIP012US us pin connector



Ø 0.4 x 4 mm nano™ tip straight

nano™ series / working / total length 135 / 250 mm

BIP013EU eu flat connector

BIP013US us pin connector



Ø 0.7 x 4 mm nano™ tip straight

nano™ series / working / total length 135 / 250 mm

BIP014EU eu flat connector

BIP014US us pin connector




Ø 1.0 x 4 mm nano™ tip straight


nano™ series / working / total length 135 / 250 mm

BIP015EU eu flat connector


BIP015US us pin connector




nano™ non-stick




Ø 0.4 x 4 mm nano™ tip 80° upwards angled
nano™ series / working / total length 85 / 200 mm
BIP016EU eu flat connector
BIP016US us pin connector




Ø 0.4 x 4 mm nano™ tip 80° downwards angled
nano™ series / working / total length 85 / 200 mm
BIP017EU eu flat connector
BIP017US us pin connector




Ø 0.4 x 4 mm nano™ tip 45° upwards angled
nano™ series / working / total length 85 / 200 mm
BIP018EU eu flat connector
BIP018US us pin connector



Ø 0.4 x 4 mm nano™ tip 45° downwards angled
nano™ series / working / total length 85 / 200 mm
BIP019EU eu flat connector
BIP019US us pin connector




nano™ non-stick



Ø 0.4 x 4 mm nano™ tip 80° upwards angled
nano™ series / working / total length 115 / 230 mm
BIP020EU eu flat connector
BIP020US us pin connector



Ø 0.4 x 4 mm nano™ tip 80° downwards angled
nano™ series / working / total length 115 / 230 mm
BIP021EU eu flat connector
BIP021US us pin connector



Ø 0.4 x 4 mm nano™ tip 45° upwards angled
nano™ series / working / total length 115 / 230 mm
BIP022EU eu flat connector
BIP022US us pin connector




Ø 0.4 x 4 mm nano™ tip 45° downwards angled
nano™ series / working / total length 115 / 230 mm
BIP023EU eu flat connector
BIP023US us pin connector





classic non-stick bipolar forceps



Tips with 6 or 8 mm tip length



TruPar™ pin guide



Wide ergonomic handle



US pin or EU flat connectors

classic non-stick



Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 30 / 110 mm

BIP024EU eu flat connector

BIP024US us pin connector



Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 30 / 110 mm

BIP025EU eu flat connector

BIP025US us pin connector



Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 30 / 110 mm

BIP026EU eu flat connector

BIP026US us pin connector



Ø 0.7 x 8 mm standard tip 45° angled

standard series / working / total length 30 / 110 mm

BIP027EU eu flat connector

BIP027US us pin connector



Ø 1.0 x 8 mm standard tip straight


standard series / working / total length 30 / 110 mm

BIP028EU eu flat connector


BIP028US us pin connector



1:1



classic non-stick




Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 45 / 185 mm

BIP029EU eu flat connector

BIP029US us pin connector




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 45 / 185 mm

BIP030EU eu flat connector

BIP030US us pin connector




Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 45 / 185 mm

BIP031EU eu flat connector

BIP031US us pin connector




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 45 / 185 mm

BIP032EU eu flat connector

BIP032US us pin connector




Ø 1.0 x 8 mm standard tip 45° angled


standard series / working / total length 45 / 185 mm

BIP033EU eu flat connector

BIP033US us pin connector



classic non-stick




Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 60 / 200 mm

BIP034EU eu flat connector

BIP034US us pin connector




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 60 / 200 mm

BIP035EU eu flat connector

BIP035US us pin connector




Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 60 / 200 mm

BIP036EU eu flat connector

BIP036US us pin connector




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 60 / 200 mm

BIP037EU eu flat connector

BIP037US us pin connector




Ø 1.0 x 8 mm standard tip 45° angled

standard series / working / total length 60 / 200 mm

BIP038EU eu flat connector

BIP038US us pin connector




Ø 2.0 x 8 mm standard tip 45° angled


standard series / working / total length 60 / 200 mm

BIP039EU eu flat connector


BIP039US us pin connector



classic non-stick



Ø 1.5 x 8 mm standard tip straight
standard series / working / total length 90 / 230 mm
BIP040EU eu flat connector
BIP040US us pin connector



Ø 2.0 x 8 mm standard tip straight
standard series / working / total length 90 / 230 mm
BIP041EU eu flat connector
BIP041US us pin connector

TruPar™ pin guide

wide handle



Ø 1.0 x 8 mm standard tip 45° angled
standard series / working / total length 120 / 260 mm
BIP042EU eu flat connector
BIP042US us pin connector



Ø 2.0 x 8 mm standard tip 45° angled
standard series / working / total length 120 / 260 mm
BIP043EU eu flat connector
BIP043US us pin connector

TruPar™ pin guide

wide handle

1:1

classic non-stick



Ø 2.5 x 8 mm standard tip straight

standard series / working / total length 140 / 280 mm

BIP044EU eu flat connector

BIP044US us pin connector

TruPar™ pin guide

wide handle

1:1

classic non-stick




Ø 2.0 x 8 mm standard tip straight


standard series / working / total length 50 / 180 mm

BIP045EU eu flat connector

BIP045US us pin connector



classic non-stick




Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 40 / 155 mm

BIP046EU eu flat connector

BIP046US us pin connector




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 40 / 155 mm

BIP047EU eu flat connector

BIP047US us pin connector




Ø 0.7 x 6 mm micro tip straight


standard series / working / total length 40 / 155 mm

BIP048EU eu flat connector

BIP048US us pin connector



classic non-stick




Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 55 / 170 mm

BIP049EU eu flat connector

BIP049US us pin connector




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 55 / 170 mm

BIP050EU eu flat connector

BIP050US us pin connector




Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 55 / 170 mm

BIP051EU eu flat connector

BIP051US us pin connector



Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 55 / 170 mm

BIP052EU eu flat connector

BIP052US us pin connector


TruPar™ pin guide

wide handle



1:1

classic non-stick




Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 85 / 200 mm

BIP053EU eu flat connector

BIP053US us pin connector




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 85 / 200 mm

BIP054EU eu flat connector

BIP054US us pin connector




Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 85 / 200 mm

BIP055EU eu flat connector

BIP055US us pin connector




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 85 / 200 mm

BIP056EU eu flat connector

BIP056US us pin connector




Ø 1.2 x 8 mm standard tip straight

standard series / working / total length 85 / 200 mm

BIP057EU eu flat connector

BIP057US us pin connector



Ø 2.0 x 8 mm standard tip straight

standard series / working / total length 85 / 200 mm

BIP058EU eu flat connector

BIP058US us pin connector


TruPar™ pin guide

wide handle



1:1

classic non-stick



Ø 1.2 x 8 mm standard tip 45° downwards angled
standard series / working / total length 85 / 200 mm
BIP059EU eu flat connector
BIP059US us pin connector



Ø 1.2 x 8 mm standard tip 45° upwards angled
standard series / working / total length 85 / 200 mm
BIP060EU eu flat connector
BIP060US us pin connector


classic non-stick

Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 115 / 230 mm

BIP061EU eu flat connector

BIP061US us pin connector




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 115 / 230 mm

BIP062EU eu flat connector

BIP062US us pin connector




Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 115 / 230 mm

BIP063EU eu flat connector

BIP063US us pin connector




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 115 / 230 mm

BIP064EU eu flat connector

BIP064US us pin connector




Ø 1.2 x 8 mm standard tip straight

standard series / working / total length 115 / 230 mm

BIP065EU eu flat connector

BIP065US us pin connector



TruPar™ pin guide

wide handle


1:1

TruPar™ pin guide


wide handle

1:1


classic non-stick




Ø 0.7 x 8 mm standard tip 45° angled upwards
standard series / working / total length 115 / 230 mm
BIP066EU eu flat connector
BIP066US us pin connector



Ø 1.2 x 8 mm standard tip 45° angled upwards
standard series / working / total length 115 / 230 mm
BIP067EU eu flat connector
BIP067US us pin connector



Ø 0.7 x 8 mm standard tip 45° angled downwards
standard series / working / total length 115 / 230 mm
BIP068EU eu flat connector
BIP068US us pin connector



Ø 1.2 x 8 mm standard tip 45° angled downwards
standard series / working / total length 115 / 230 mm
BIP069EU eu flat connector
BIP069US us pin connector

TruPar™ pin guide

wide handle



1:1

classic non-stick



Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 135 / 250 mm

BIP070EU eu flat connector

BIP070US us pin connector



Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 135 / 250 mm

BIP071EU eu flat connector

BIP071US us pin connector



Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 135 / 250 mm

BIP072EU eu flat connector

BIP072US us pin connector



Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 135 / 250 mm

BIP073EU eu flat connector

BIP073US us pin connector



Ø 1.2 x 8 mm standard tip straight

standard series / working / total length 135 / 250 mm

BIP074EU eu flat connector


BIP074US us pin connector

TruPar™ pin guide

wide handle

1:1

classic non-stick




Ø 0.7 x 8 mm standard tip 45° angled upwards

standard series / working / total length 135 / 250 mm

BIP075EU eu flat connector

BIP075US us pin connector




Ø 1.2 x 8 mm standard tip 45° angled upwards

standard series / working / total length 135 / 250 mm

BIP076EU eu flat connector

BIP076US us pin connector




Ø 0.7 x 8 mm standard tip 45° angled downwards

standard series / working / total length 135 / 250 mm

BIP077EU eu flat connector

BIP077US us pin connector



Ø 1.2 x 8 mm standard tip 45° angled downwards


standard series / working / total length 135 / 250 mm

BIP078EU eu flat connector

BIP078US us pin connector



non-stick irrigation bipolar forceps




Integrated irrigation channel




TruPar™ Pin guide




Luer Lock irrigation connector




non-stick irrigation




Ø 1.5 x 8 mm standard tip straight
standard series / working / total length 100 / 200 mm
BIP001EU-IR eu flat connector
BIP001US-IR us pin connector




non-stick irrigation



Ø 0.5 x 6 mm micro tip straight
standard series / working / total length 90 / 200 mm
BIP002EU-IR eu flat connector
BIP002US-IR us pin connector



Ø 1.0 x 8 mm standard tip straight
standard series / working / total length 90 / 200 mm
BIP003EU-IR eu flat connector
BIP003US-IR us pin connector



non-stick irrigation




Ø 0.5 x 6 mm micro tip straight
standard series / working / total length 120 / 230 mm
BIP004EU-IR eu flat connector
BIP004US-IR us pin connector



Ø 1.0 x 8 mm standard tip straight
standard series / working / total length 120 / 230 mm
BIP005EU-IR eu flat connector
BIP005US-IR us pin connector



suction bipolar forceps




Integrated suction channel




Wide ergonomic grip




Standardized suction connector



suction forceps




Ø 1.4 mm standard tip straight
standard series / working / total length 90 / 200 mm
BIP001EU-SUC eu flat connector
BIP001US-SUC us pin connector



Ø 2.0 mm standard tip straight
standard series / working / total length 90 / 200 mm
BIP002EU-SUC eu flat connector
BIP002US-SUC us pin connector



suction forceps




Ø 1.4 mm standard tip straight

standard series / working / total length 120 / 230 mm

BIP003EU-SUC eu flat connector

BIP003US-SUC us pin connector




Ø 2.0 mm standard tip straight


standard series / working / total length 120 / 230 mm

BIP004EU-SUC eu flat connector


BIP004US-SUC us pin connector



suction forceps



Ø 2.0 standard tip straight
standard series / working / total length 75 / 190 mm
BIP005EU-SUC eu flat connector
BIP005US-SUC us pin connector



meuser suction forceps




Ø 3.0 standard tip straight


standard series / working / total length 75 / 210 mm

BIP006EU-SUC eu flat connector

BIP006US-SUC us pin connector



meuser suction forceps



Ø 3.2 standard tip straight
standard series / working / total length 75 / 210 mm
BIP007EU-SUC eu flat connector
BIP007US-SUC us pin connector

Lehmberg Keyhole™

non-stick bipolar forceps

nxt™ non stick - solid silver tips

Lehmberg Keyhole™ non-stick bipolar forceps

The Lehmberg Keyhole™ non-stick bipolar forceps were developed in close collaboration with Prof. Dr. med. Jens Lehmberg. Engineered for precision in narrow surgical spaces, Lehmberg Keyhole™ bipolar forceps combine the tactile feedback of traditional forceps with innovative functionality. Overcompression of the forceps reopens the tips which allows for working in very narrow spaces.

Lehmberg Keyhole™ non-stick bipolar forceps with instrument tips made of solid silver (not coating) are characterised by very good heat transfer. This prevents overheating of the tips during coagulation. The non-stick properties remain intact throughout the entire lifespan of the instruments.

Die Lehmberg Keyhole™ Non-Stick Bipolarpinzetten wurden in enger Zusammenarbeit mit Prof. Dr. med. Jens Lehmberg entwickelt. Das Entwicklungsziel bestand darin, eine Keyhole Bipolarpinzette zu entwickeln, welche das Arbeiten unter sehr beengten Platzverhältnissen ermöglicht, ohne auf das taktile Feedback einer klassischen Bipolarpinzette verzichten zu müssen. Erreicht wurde dieses Ziel durch ein spezielles Schenkeldesign, welches die Pinzettenspitzen wieder öffnet, sobald die Pinzette überdrückt wird.

Alle Lehmberg Keyhole™ Non-Stick Bipolarpinzetten sind mit Instrumentenspitzen aus solidem Silber (keine Beschichtung) ausgestattet. Dadurch wird eine Überhitzung der Pinzettenspitzen während der Koagulation verhindert. Die Non-Stick Eigenschaften bleiben über die gesamte Lebensdauer der Pinzetten erhalten.



Prof. Dr. med. Jens Lehmberg, MHBA

Chair, Department of Neurosurgery

Munich Clinic Bogenhausen & Munich Clinic Schwabing

Munich Municipal Hospital Group (MMHG), Germany



Position 1

Open Forceps.
Offene Pinzette.



Position 2

Closed Forceps. In this state, the forceps can be inserted into very narrow spaces (keyhole surgery).
Geschlossene Pinzette. In diesem Zustand können die Pinzetten in sehr enge Zugänge gebracht werden.



Position 3

Overcompression of the forceps reopens the tips which allows for working in very narrow spaces.
Überdrücken öffnet die Pinzettenspitzen. Dadurch wird das Arbeiten unter sehr beengten Platzverhältnissen ermöglicht.

Lehmberg Keyhole™

non-stick bipolar forceps



Lehmberg Keyhole™ non-stick bipolar forceps



nano™ tips (3 mm length)



nano™ Ø 2.5 mm tines
slightly curved upwards for optimal
tip visibility under the microscope



Multi position wide handle



US pin or EU flat connectors

Lehmberg Keyhole™
non-stick bipolar forceps



Lehmberg Keyhole™

non-stick bipolar forceps



Lehmberg Keyhole™



Ø 0.2 x 3 mm nano™ tip straight

nano™ series / working / total length 85 / 200 mm

LF01EU eu flat connector

LF01US us pin connector



Ø 0.4 x 3 mm nano™ tip straight

nano™ series / working / total length 85 / 200 mm

LF02EU eu flat connector

LF02US us pin connector



Ø 0.2 x 3 mm nano™ tip 40° up

nano™ series / working / total length 85 / 200 mm

LF03EU eu flat connector

LF03US us pin connector



Ø 0.4 x 3 mm nano™ tip 40° up

nano™ series / working / total length 85 / 200 mm

LF04EU eu flat connector

LF04US us pin connector

Lehmborg Keyhole™

non-stick bipolar forceps



Lehmborg Keyhole™



Ø 0.2 x 3 mm nano™ tip straight

nano™ series / working / total length 110 / 225 mm

LF05EU eu flat connector

LF05US us pin connector



Ø 0.4 x 3 mm nano™ tip straight

nano™ series / working / total length 110 / 225 mm

LF06EU eu flat connector

LF06US us pin connector



Ø 0.2 x 3 mm nano™ tip 40° up

nano™ series / working / total length 110 / 225 mm

LF07EU eu flat connector

LF07US us pin connector



Ø 0.4 x 3 mm nano™ tip 40° up

nano™ series / working / total length 110 / 225 mm

LF08EU eu flat connector

LF08US us pin connector

Lehmborg Keyhole™

non-stick bipolar forceps



Lehmborg Keyhole™



Ø 0.2 x 3 mm nano™ tip straight
nano™ series / working / total length 130 / 245 mm
LF09EU eu flat connector
LF09US us pin connector



Ø 0.4 x 3 mm nano™ tip straight
nano™ series / working / total length 130 / 245 mm
LF10EU eu flat connector
LF10US us pin connector



Ø 0.2 x 3 mm nano™ tip 40° up
nano™ series / working / total length 130 / 245 mm
LF11EU eu flat connector
LF11US us pin connector



Ø 0.4 x 3 mm nano™ tip 40° up
nano™ series / working / total length 130 / 245 mm
LF12EU eu flat connector
LF12US us pin connector



— nano™ tips 3 mm length

— nano™ tines 2.5 mm diameter

slightly curved upwards for optimal
tip visibility under the microscope

— multi position
wide handle

1:1

Lehmborg Keyhole™

non-stick bipolar forceps

Lehmborg Keyhole™



Ø 0.2 x 3 mm nano™ tip straight

nano™ series / working / total length 150 / 265 mm

LF13EU eu flat connector

LF13US us pin connector



Ø 0.4 x 3 mm nano™ tip straight

nano™ series / working / total length 150 / 265 mm

LF14EU eu flat connector

LF14US us pin connector



Ø 0.2 x 3 mm nano™ tip 40° up

nano™ series / working / total length 150 / 265 mm

LF15EU eu flat connector

LF15US us pin connector



Ø 0.4 x 3 mm nano™ tip 40° up


nano™ series / working / total length 150 / 265 mm

LF16EU eu flat connector


LF16US us pin connector




round handle non-stick bipolar forceps




nano™ tips (4 mm length)



Round-handle design



TruPar™ Pin guide



US pin or EU flat connectors



Ø 0.5 x 4 mm straight

working / total length 50 / 150 mm

BIP-S141EU eu flat connector

BIP-S141US us pin connector



Ø 0.8 x 4 mm straight

working / total length 50 / 150 mm

BIP-S142EU eu flat connector

BIP-S142US us pin connector



Ø 0.5 x 4 mm angled

working / total length 50 / 150 mm

BIP-S146EU eu flat connector

BIP-S146US us pin connector



Ø 0.8 x 4 mm angled

working / total length 50 / 150 mm

BIP-S147EU eu flat connector

BIP-S147US us pin connector



Ø 0.5 x 4 mm angled

working / total length 80 / 180 mm

BIP-S246EU eu flat connector

BIP-S246US us pin connector



Ø 1.2 x 4 mm angled

working / total length 80 / 180 mm

BIP-S248EU eu flat connector

BIP-S248US us pin connector



Ø 0.5 x 4 mm straight

working / total length 100 / 200 mm

BIP-S341EU eu flat connector

BIP-S341US us pin connector



Ø 0.8 x 4 mm straight

working / total length 100 / 200 mm

BIP-S342EU eu flat connector

BIP-S342US us pin connector



Ø 0.5 x 4 mm angled

working / total length 100 / 200 mm

BIP-S346EU eu flat connector

BIP-S346US us pin connector



Ø 0.8 x 4 mm angled

working / total length 100 / 200 mm

BIP-S347EU eu flat connector

BIP-S347US us pin connector



Ø 0.5 x 4 mm straight

working / total length 90 / 200 mm

BIP-S441EU eu flat connector

BIP-S441US us pin connector



Ø 0.8 x 4 mm straight

working / total length 90 / 200 mm

BIP-S442EU eu flat connector

BIP-S442US us pin connector



Ø 1.2 x 4 mm straight


working / total length 90 / 200 mm

BIP-S443EU eu flat connector

BIP-S443US us pin connector




Ø 0.5 x 4 mm down angled
working / total length 100 / 195 mm
BIP-S444EU eu flat connector
BIP-S444US us pin connector




Ø 0.5 x 4 mm up angled
working / total length 90 / 200 mm
BIP-S446EU eu flat connector
BIP-S446US us pin connector



Ø 0.8 x 4 mm up angled
working / total length 90 / 200 mm
BIP-S447EU eu flat connector
BIP-S447US us pin connector



Ø 1.2 x 4 mm up angled
working / total length 90 / 200 mm
BIP-S448EU eu flat connector
BIP-S448US us pin connector



Ø 0.5 x 4 mm straight

working / total length 140 / 240 mm

BIP-S541EU eu flat connector

BIP-S541US us pin connector



Ø 0.8 x 4 mm straight

working / total length 140 / 240 mm

BIP-S542EU eu flat connector

BIP-S542US us pin connector



Ø 0.5 x 4 mm angled

working / total length 140 / 240 mm

BIP-S546EU eu flat connector

BIP-S546US us pin connector



Ø 0.8 x 4 mm angled

working / total length 140 / 240 mm

BIP-S547EU eu flat connector

BIP-S547US us pin connector



Ø 0.2 x 4 mm straight

working / total length 75 / 220 mm

BIP-S641EU eu flat connector

BIP-S641US us pin connector



Ø 0.5 x 4 mm straight

working / total length 75 / 220 mm

BIP-S642EU eu flat connector

BIP-S642US us pin connector



Ø 0.8 x 4 mm straight

working / total length 75 / 220 mm

BIP-S643EU eu flat connector

BIP-S643US us pin connector



Ø 1.2 x 4 mm straight

working / total length 75 / 220 mm

BIP-S644EU eu flat connector

BIP-S644US us pin connector



Ø 2.0 x 4 mm straight

working / total length 75 / 220 mm

BIP-S645EU eu flat connector

BIP-S645US us pin connector



Ø 0.5 x 4 mm up angled

working / total length 80 / 220 mm

BIP-S646EU eu flat connector

BIP-S646US us pin connector



Ø 0.8 x 4 mm up angled

working / total length 80 / 220 mm

BIP-S647EU eu flat connector

BIP-S647US us pin connector



Ø 1.2 x 4 mm up angled

working / total length 80 / 220 mm

BIP-S648EU eu flat connector

BIP-S648US us pin connector



Ø 0.5 x 4 mm straight

working / total length 110 / 250 mm

BIP-S741EU eu flat connector

BIP-S741US us pin connector



Ø 0.8 x 4 mm straight

working / total length 110 / 250 mm

BIP-S742EU eu flat connector

BIP-S742US us pin connector



Ø 1.2 x 4 mm straight

working / total length 110 / 250 mm

BIP-S743EU eu flat connector

BIP-S743US us pin connector



Ø 0.5 x 4 mm up angled

working / total length 110 / 250 mm

BIP-S746EU eu flat connector

BIP-S746US us pin connector



Ø 0.8 x 4 mm up angled

working / total length 110 / 250 mm

BIP-S747EU eu flat connector

BIP-S747US us pin connector



Ø 1.2 x 4 mm up angled


working / total length 110 / 250 mm

BIP-S748EU eu flat connector

BIP-S748US us pin connector

BIPOLAR CABLES

RE-USABLE NON-STERILE



bipolar cables*

*not an adeor product.

*kein adeor Produkt.

Generator side



Forceps side

Forces side	EU flat	US pin
		
Valleylab, Conmed, Codman, Olympus - Banana plug*	ST-92021 3 m	ST-92114 3 m
Valleylab, Conmed, Codman, Olympus - Banana plug*	ST-92022 5 m	ST-92093 5 m
Aesculap, Martin, Berchtold, Wolf	ST-92017 3 m	ST-92113 3 m
Aesculap, Martin, Berchtold, Wolf	ST-92018 5 m	ST-92092 5 m
Erbe, Downs, Siemens	ST-92015 3 m	ST-92110 3 m
Erbe, Downs, Siemens	ST-92016 5 m	ST-92111 5 m
Valleylab, EMD, Dolley - 2-pin distance 29 mm	ST-92019 3 m	ST-92115 3 m
Valleylab, EMD, Dolley - 2-pin distance 29 mm	ST-92020 5 m	ST-92094 5 m

② multi use
autoclavable 134°C

*not available within the EU and various other markets.

*nicht in der EU und einigen anderen Märkten verfügbar.

adeor® is a registered trademark of adeor medical AG.

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing device are reserved and require our written approval.

Copyright 2024 adeor medical AG

01-1224-D/E

Printed in Germany

December 2024

Spetzler and Malis are trademarks of Stryker Inc.

Sutter is a trademark of Sutter Medizintechnik GmbH

adeor is not affiliated with any of the above.



adeor Germany (HQ)

adeor medical AG

Martinshof 5

D-83626 Valley

Germany

Phone +49 (0) 8024 47717 0

Fax +49 (0) 8024 47717 99

Email office@adeor.com

Web www.adeor.com

adeor USA

adeor medical USA INC.

1100 Arrowhead Court | Suite A

Crown Point | IN 46307

USA

Phone +1 (219) 361 1020

Fax +1 (219) 361 1030

Email usa@adeor.com

Web www.adeor.com

Printed on recycled paper.

Gedruckt auf Recyclingpapier.

