



Sugita AVM Microclip II

Manufactured by:



Mizuho Corporation

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN
<https://www.mizuho.co.jp>

Distributed by:

Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Caution


- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

Sugita AVM Microclip II

A complete solution for Cerebral Arteriovenous Malformation and other small vessel occlusions
Permanent Microclips / Appliers / Tray

Sugita AVM Microclip II

- Permanent clip supplied in dedicated microclip case
- The clip case is designed for easy grip and extraction of microclips, enabling more efficient surgical procedures
- Compatible with MRI scanning at up to 3-tesla
- Made from cobalt-chrome alloy
- Ideal for AVM and other small vessels
- A closing force within a range of 50-70gf




Product No.	Description
0 7 - 9 3 2 - 1 2	Sugita AVM Microclip II 2 mm
0 7 - 9 3 2 - 1 3	Sugita AVM Microclip II 3 mm
0 7 - 9 3 2 - 1 4	Sugita AVM Microclip II 4 mm
0 7 - 9 3 2 - 1 5	Sugita AVM Microclip II 5 mm



Sugita AVM Tray

- Designed for storage, transport and sterilization of microclips and standard appliers*
- Fully autoclavable


*See IFU for details



Product No.	Description	Component part
0 7 - 9 3 2 - 0 9	SUGITA AVM Microclip II CASE Complete Set	①,②,③
0 7 - 9 3 2 - 1 0	SUGITA AVM Microclip II CASE A	①,③
0 7 - 9 3 2 - 1 1	SUGITA AVM Microclip II CASE B	②

Standard Applier

- Standard locking applier optimized for AVM procedures
- Made from stainless steel with cobalt-chrome alloy jaws




Product No.	Description	Overall Length	Working Length
0 7 - 9 3 5 - 3 0	Sugita AVM Microclip II Applier 100 mm	212 mm	100 mm
0 7 - 9 3 5 - 3 1	Sugita AVM Microclip II Applier 50 mm	165 mm	50 mm

Single Shaft Applier

- Ultra low-profile single shaft for improved visualization during clip insertion
- Ideal for vessel occlusions in remote regions such as skull base

*Not compatible with Sugita AVM Tray



Product No.	Description	Overall Length	Working Length
0 7 - 9 3 5 - 5 2	Sugita AVM Microclip II Single Shaft Applier 90 mm, Str, w/o Catches	244 mm	90 mm
0 7 - 9 3 5 - 5 3	Sugita AVM Microclip II Single Shaft Applier 90 mm, 30 deg, w/o Catches		
0 7 - 9 3 5 - 5 4	Sugita AVM Microclip II Single Shaft Applier 110 mm, Str, w/o Catches	264 mm	110 mm
0 7 - 9 3 5 - 5 5	Sugita AVM Microclip II Single Shaft Applier 110 mm, 30 deg, w/o Catches		
0 7 - 9 3 5 - 5 6	Sugita AVM Microclip II Single Shaft Applier 150 mm, Str, w/o Catches	304 mm	150 mm
0 7 - 9 3 5 - 5 7	Sugita AVM Microclip II Single Shaft Applier 150 mm, 30 deg, w/o Catches		