



Dr. Youssef's Cranial Nerve Dissection Set

Manufactured by:

MIZUHO AMERICA



Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com

Distributed by:

Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Caution

- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

DR. YOUSSEF'S CRANIAL NERVE DISSECTION SET

Comprehensive, high quality, specialized set for nerve dissection

Comprehensive Set for Tumor Treatment



Spatula Dissectors
1.8 mm, 2.2 mm, & 2.4 mm
Tip Sizes




Round Dissectors
2.0 mm Tip Size with 30° &
45° Up Angles



Putter Dissectors
60° Up, Left & Right Angles



Teardrop Probes
15°, 45°, & 90° Up Angles
with Ball Tips



Needle Dissectors
Straight, 45°, & 90° Up Angles



Cup Curettes
45° & 90° Up Angles,
100° Reverse Angle,
2.0 x 2.5 mm & 1.5 x 2.0 mm
Cup Sizes



Teardrop Back-cutting Knife
90° Up Angle with Ball Tip



Recommended intended use for:
acoustic neuroma, meningiomas, other schwannomas.

Dr. Youssef's Cranial Nerve Dissection Set

SUPERIOR DESIGN

Comprehensive set for skull base and ENT tumor treatment

Lower profile instruments to accommodate minimally invasive procedures

Extended working length to adapt to a full range of skull base procedures



Uniquely designed tips to mitigate cranial nerve trauma

Vividly colored handles for better instrument identification

Gentle manipulation of tumors

