



Microscissors

Master Edition

Manufactured by:

MIZUHO AMERICA



Mizuho America Inc.
30057 Ahern Avenue Union City, CA 94587
www.mizuho.com

Distributed by:

Mizuho America Inc.
30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Caution

- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

LawtonElite Microscissors

Ultra-thin design with sharp, curved blades. Non-slip counterbalanced ergonomic handles. In the closed position, microscissors can be used as a micro dissector. Available in Reusable or Reposable options.

Titanium, Straight

- 6.0 cm Working Length
- 18.0 cm Total Length
- Reusable



Titanium, Bayonet

- 7.6 cm Working Length
- 18.5 cm Total Length
- Curved Left or Right
- Reusable



Reposable, Straight

- 5.0 cm Working Length
- 18.0 cm Total Length
- Reposable - always sharp for the life of the product



Reposable, Bayonet

- 10.0 cm Working Length
- 21.5 cm Total Length
- Curved Left or Right
- Reposable – always sharp for the life of the product



Microscissor Tips



LawtonElite Curved Tips

MICROSCISSORS FOR Neurovascular Bypass Applications

Specially crafted microscissors for varying depths of cranial bypass procedures.



Lawton Straight



Lawton Curved



Micro Arteriotomy 60° Up



Micro Sylvian Fissure Straight Tips



Micro Sylvian Fissure 30° Up Tips



Lawton Microscissors

- 6.5 cm Working Length
- 18.0 cm Total Length
- Straight or Curved
- Round or Flat Handle



Micro Arteriotomy Microscissors

- 3.0 cm & 6.0 cm Working Lengths
- 14.5 cm & 17.5 cm Total Lengths
- 60° Up Extra Fine Tip
- Round Handles
- Guillotine Blade Action



Micro Sylvian Fissure Microscissors

- 3.0 cm, 6.0 cm, & 7.0 cm Working Lengths
- 14.5 cm, 18.0 cm, 19.0 cm, & 21.0 Total Lengths
- Straight, 30° Up, 30° Angled Blades & 30° Angled Blades with Irrigation
- Round Handle
- Flush Port for Cleaning
- Closed Blades Options

Microscissors

- 9.5 cm Working Length
- 20.5 cm Total Length
- 60° Up
- Finger Pads



Designed in collaboration with
MICHAEL T. LAWTON, MD
President and CEO, Barrow Neurological Institute
Professor and Chair, Neurosurgery
Chief, Neurovascular Surgery

MICROSCISSORS FOR Multiple-use Cranial and Spinal Applications


Designed for every aspect of neurosurgery, Mizuho America’s portfolio of Microscissors is one of the most comprehensive line-ups available in the industry. Uniquely designed slightly curved microscissors allow greater visualization at the operative site. Serrated blades allow for grasping tissue with a clean-cut, and they stay sharper, longer.

Fukushima/Giannotta Traditional Microscissors

- 6.5 cm, 8.5 cm, 9.0 cm, 11.0 cm, 13.0 cm, & 15.0 cm Working Lengths
- 20.0 cm, 20.5 cm, 22.0 cm, 23.0 cm, 24.0 cm, 25.0 cm, 25.5 cm, 26.0 cm, 26.5, & 27.5 Total Lengths
- Straight or Curved Tips
- Round or Flat Handles
- Bayonet Design




Traditional Microscissors Straight Tips



Traditional Microscissors Curved Tips

Fukushima/Giannotta Ultra Fine Microscissors

- 6.5 cm, 8.5 cm, 11.0 cm, 12.5 cm, 13.0 cm, & 15.5 cm Working Lengths
- 15.0 cm, 20.0 cm, 20.5 cm, 22.0 cm, 22.5 cm, 23.0 cm, 24.0 cm, 25.0 cm, 26.0 cm, 26.5 cm & 27.5 cm Total Lengths
- Straight or Curved Tips
- Round or Flat Handle
- Bayonet Design
- Ideal for Dissection of Arachnoid
- Extremely Thin Blades with Gold Springs



Ultra Fine Microscissors Straight Tips



Ultra Fine Microscissors Curved Tips

Metzenbaum Style Microscissors

Longer blades with blunt tips and flat handles. Designed for blunt and sharp dissection.



Meningioma Microscissors

- 9.0 cm & 11.5 cm Working Lengths
- 21.0 cm & 22.5 cm Total Lengths
- Straight or Curved Tips
- Flat Handle
- Bayonet Design
- Short Blades with Sharp Tips & Black Springs
- Ideal for Cutting Thick Vascular & Scar Tissue



Meningioma Microscissors Straight Tips



Meningioma Microscissors Curved Tips

Blunt Microscissors

- 9.0 cm Working Length
- 20.5 cm Total Length
- Straight and Curved Tips
- Flat Handle
- Bayonet Design
- Longer Blades
- Ideal for Cutting Thick Vascular & Scar Tissue



Blunt Microscissors Straight Tips




Blunt Microscissors Curved Tips

Endonasal Microscissors


Designed to minimize visual obstruction and maximize maneuverability.

Pistol - Grip Microscissors

- 14.0 cm Working Length
- 21.0 cm & 21.5 cm Total Lengths
- Straight & Curved, Left & Right
- Single Shaft



Pistol Grip Microscissors Straight Tips




Pistol Grip Microscissors Curved Tips



Pistol Grip Microscissors Curved Up Tips

Endoscopic Microscissors

- 12.0 cm Working Length
- 22.0 cm Total Length
- Straight, Curved Right & Left
- Pistol-Grip Handle
- Single Shaft



Endoscopic Microscissors Straight Tips



Endoscopic Microscissors Curved Tips


XPS SERIES MICROSCISSORS FOR Sharp Dissection

Extremely thin blades with gold-plated spring, well balanced and comfortably weighted flat handle scissors.


Features & Benefits

- 8.5 cm, 9.0 cm & 10.0 cm Working Lengths
- 17.5 cm & 19.5 cm Total Lengths
- Curved or Straight Non-Serrated Blades
- Gold Plated Springs
- Thin, Sharp Blades
- Designed in collaboration with Dr. Aaron Cohen-Gadol, MD
- Straight & Bayonet Flat Handles

Microscissor Tips



Curved Lateral



Curved Left



Curved Right



Curved Up



Straight



360° EVANS ROTATABLE MICROSCISSORS FOR Endoscopic & Micro Neurosurgical Procedures

360° fully rotatable shafts for endonasal, transcranial endoscopic and microscopic surgery.



Detachable Flush Port

Features & Benefits

- 9.5 cm, 11.5 cm, 12.0 cm, 15.0 cm, 15.5 cm, & 16.0 cm Working Lengths
- Options: 45° Up, Straight, Curved, 125° Back Cutting, & Long Blades
- Single-Shaft Handles
- Axial Bayonet Handles
- Detachable Flush-Port for Improved Visibility & Cleaning
- Safely Prevents Rotation when Jaws are Closed
- Well Balanced for Comfort
- Serrated Blades for Sharp Dissection
- Long Blade Options

Microscissor Tips



Straight



45° Up



Curved



Curved, Extended Blades



125° Back Cutting

Ordering Information

LAWTON ELITE MICROSCISSORS

Product #	Description
10203MA	LawtonElite Microscissors Straight, 18cm
10204MA	LawtonElite Microscissors Bayonet, Curved Right Jaw, 7.6cm WL
10214MA	LawtonElite Microscissors Straight, Reposable, 18cm (Box of 3)
10214MA-1	LawtonElite Microscissors Straight, Reposable, 18cm
10215MA	LawtonElite Microscissors Bayonet, Curved Left Jaw, 7.6cm WL
10294MA	LawtonElite Microscissors Bayonet, Reposable, Curved Right Blades, 10cm WL
10294MA-3	LawtonElite Microscissors Bayonet, Reposable, Curved Right Blades, 10cm WL (Box of 3)
10295MA	LawtonElite Microscissors Bayonet, Reposable, Curved Left Blades, 10cm WL
10295MA-3	LawtonElite Microscissors Bayonet, Reposable, Curved Left Blades, 10cm WL (Box of 3)
07-827-30	Microscissors Straight, Straight Jaw, Round Handle, 18cm
07-827-30F	Microscissors Straight, Straight Jaw, Flat Handle, Ultrafin Blades, 18cm
07-827-31	Microscissors Straight, Curved Jaw, Round Handle, 18cm
07-827-31F	Microscissors Straight, Curved Jaw, Flat Handle, 18cm
07-827-33	Microscissors Straight, Angled Jaw, Finger Pad, 9.5cm WL
07-827-34	Tenotomy Scissors, Straight Jaw, Ring, 11.5cm
07-827-35	Tenotomy Scissors, Curved Jaw, Ring, 11.5cm
07-827-48	Microscissors Straight, Straight Jaw, Round Handle, 3cm WL
07-827-49	Microscissors Straight, Straight Jaw, Round Handle, Irrigation, 3cm WL
07-827-50	Microscissors Straight, Angled Jaw, 60 Degrees Up, Round Handle, 6cm WL
07-827-51	Microscissors Straight, Angled Jaw, 30 Degrees Up, Round Handle, 3cm WL
07-827-52	Microscissors Straight, Angled Jaw, 60 Degrees Up, Round Handle, 3cm WL
07-CVM-02	Microscissors Straight, Straight Jaw, Round Handle, Closed Blades, 6cm WL
07-CVM-03	Microscissors Straight, Straight Jaw, Round Handle, 7.5mm WL
07-CVM-32	Microscissors Straight, Straight Jaw, Round Handle, Closed Blades, Irrigation, 6cm WL
07-CVM-33	Microscissors Straight, Straight Jaw, Round Handle, Closed Blades, Irrigation, 3cm WL

TRADITIONAL MICROSCISSORS

Product #	Description
07-827-01	Microscissors Bayonet, Straight Jaw, Round Handle, Fine Blades, 6.5cm WL
07-827-01F	Microscissors Bayonet, Straight Jaw, Flat Handle, Fine Blades, 8.5cm WL
07-827-02	Microscissors Bayonet, Curved Jaw, Round Handle, Fine Blades, 6.5cm WL
07-827-02F	Microscissors Bayonet, Curved Jaw, Flat Handle, Fine Blades, 8.5cm WL
07-827-03	Microscissors Bayonet, Straight Jaw, Round Handle, Fine Blades, 8.5cm WL
07-827-03F	Microscissors Bayonet, Straight Jaw, Flat Handle, Fine Blades, 11cm WL
07-827-04	Microscissors Bayonet, Curved Jaw, Round Handle, Fine Blades, 8.5cm WL
07-827-04F	Microscissors Bayonet, Curved Jaw, Flat Handle, Fine Blades, 11cm WL
07-827-05	Microscissors Bayonet, Straight Jaw, Round Handle, Fine Blades, 11cm WL
07-827-11	Microscissors Bayonet, Curved Right Jaw, Flat Handle, 9.5cm WL
07-827-12	Microscissors Straight, Curved Jaw, Flat Handle, 17cm
07-827-15	Microscissors Straight, Straight Jaw, Flat Handle, Fine Blades, 15cm
07-827-18F	Microscissors Straight, Straight Jaw, Flat Handle, Fine Blades, 18cm
07-827-19F	Microscissors Straight, Curved Jaw, Flat Handle, Fine Blades, 18cm

ULTRA – FINE MICROSCISSORS

Product #	Description
07-827-20	Microscissors Straight, Straight Jaw, Round Handle, 15cm
07-827-20F	Microscissors Straight, Straight Jaw, Flat Handle, Ultrafin Blades, 15cm
07-827-21	Microscissors Straight, Curved Jaw, Round Handle, 15cm
07-827-21F	Microscissors Straight, Curved Jaw, Flat Handle, 15cm
07-827-22	Microscissors Bayonet, Straight Jaw, Round Handle, Ultrafin Blades, 6.5cm WL
07-827-22F	Microscissors Bayonet, Straight Jaw, Flat Handle, Ultrafin Blades, 8.5cm WL
07-827-23	Microscissors Bayonet, Curved Jaw, Round Handle, Ultrafin Blades, 6.5cm WL
07-827-23F	Microscissors Bayonet, Curved Jaw, Flat Handle, Ultrafin Blades, 8.5cm WL
07-827-24	Microscissors Bayonet, Straight Jaw, Round Handle, Ultrafin Blades, 8.5cm WL
07-827-24F	Microscissors Bayonet, Straight Jaw, Flat Handle, Ultrafin Blades, 11cm WL
07-827-25	Microscissors Bayonet, Curved Jaw, Round Handle, Ultrafin Blades, 8.5cm WL
07-827-25F	Microscissors Bayonet, Curved Jaw, Flat Handle, Ultrafin Blades, 11cm WL
07-827-26	Microscissors Bayonet, Straight Jaw, Round Handle, Ultrafin Blades, 11cm WL
07-827-26F	Microscissors Bayonet, Straight Jaw, Flat Handle, Ultrafin Blades, 13cm WL
07-827-27	Microscissors Bayonet, Curved Jaw, Round Handle, Ultrafin Blades, 11cm WL
07-827-27F	Microscissors Bayonet, Curved Jaw, Flat Handle, Ultrafin Blades, 13cm WL
07-827-28	Microscissors Bayonet, Straight Jaw, Round Handle, Ultrafin Blades, 12.5cm WL
07-827-28F	Microscissors Bayonet, Straight Jaw, Flat Handle, Ultrafin Blades, 15.5cm WL
07-827-29	Microscissors Bayonet, Curved Jaw, Round Handle, Ultrafin Blades, 12.5cm WL

METZENBAUM STYLE MICROSCISSORS

Product #	Description
09-061-10	Mayo Scissors, 17.5cm Straight
09-071-10	Metzenbaum Scissors, 18cm Straight
09-071-11	Metzenbaum Scissors, 18cm Curved
09-071-20F	Microscissors Bayonet, Straight Jaw, Flat Handle, Round Blades, 9cm WL
09-071-21F	Microscissors Bayonet, Curved Jaw, Flat Handle, Round Blades, 9cm WL
09-071-22F	Microscissors Bayonet, Straight Jaw, Flat Handle, Round Blades, 11.5cm WL
09-071-23F	Microscissors Bayonet, Curved Jaw, Flat Handle, Round Blades, 11.5cm WL
09-081-10	Metzenbaum Scissors, 14.5cm Curved
07-827-09	Microscissors Bayonet, Straight Jaw, Flat Handle, 9cm WL
07-827-10	Microscissors Bayonet, Curved Jaw, Flat Handle, 9cm WL

SINGLE SHAFT PISTOL GRIP MICROSCISSORS

Product #	Description
07-134-10	Endoscopic Microscissors, Straight Jaw, Pistol-Grip, 14cm WL
07-134-11	Endoscopic Microscissors, Curved Right Jaw, Pistol-Grip, 14cm WL
07-134-12	Endoscopic Microscissors, Curved Left Jaw, Pistol-Grip, 14cm WL
07-134-13	Endoscopic Microscissors, Curved Up Jaw, Pistol-Grip, 14cm WL
07-181-17MA	Endoscopic Microscissors, Straight Jaw, Pistol-Grip, 12cm WL
07-181-18MA	Endoscopic Microscissors, Curved Right Jaw, Pistol-Grip, 12cm WL
07-181-19MA	Endoscopic Microscissors, Curved Left Jaw, Pistol-Grip, 12cm WL

XPS SERIES MICROSCISSORS

Product #	Description
99-600-58	Microscissors Straight, Curved Lateral Jaw, Flat Handle, 17.5cm
99-600-59	Microscissors Bayonet, Curved Left Jaw, Flat Handle, Slim Blades, 9cm WL
99-600-60	Microscissors Bayonet, Curved Right Jaw, Flat Handle, Slim Blades, 9cm WL
99-600-61	Microscissors Bayonet, Curved Up Jaw, Flat Handle, Slim Blades, 9cm WL
99-600-63	Microscissors Bayonet, Straight Jaw, Flat Handle, 10cm WL
07-827-70MA	Microscissors Straight, Straight Jaw, Flat Handle, Short Blades, 20cm
07-827-71MA	Microscissors Bayonet, Curved Left Jaw, Flat Handle, Short Blades, 8.5cm WL
07-827-72MA	Microscissors Bayonet, Curved Right Jaw, Flat Handle, Short Blades, 8.5cm WL
07-827-73MA	Microscissors Bayonet, Curved Left Jaw, Flat Handle, Short Blades, 9.5cm WL
07-827-74MA	Microscissors Bayonet, Curved Right Jaw, Flat Handle, Short Blades, 9.5cm WL

E3 EVANS ELITE 360° ROTATABLE MICROSCISSORS

Product #	Description
07-827-61	Rotatable Microscissors, Straight Jaw, Single Shaft Handle, 9.5cm WL
07-827-62	Rotatable Microscissors, Straight Jaw, Single Shaft Handle, 11.5cm WL
07-827-63	Rotatable Microscissors, Curved Jaw, Single Shaft Handle, 10cm WL
07-827-64	Rotatable Microscissors, Angled Jaw, 45 Degrees Up, Single Shaft Handle, 12cm WL
07-827-65	Rotatable Microscissors, Angled Jaw, 45 Degrees Lateral, Single Shaft Handle, 9.5cm WL
07-827-66	Rotatable Microscissors, Angled Jaw, 45 Degrees Lateral, Single Shaft Handle, 11.5cm WL
07-829-41	Rotatable Microscissors, Angled Jaw, 45 Degrees Up, Single Shaft Handle, 11.5cm WL
07-829-42	Rotatable Microscissors, Angled Jaw, 45 Degrees Up, Single Shaft Handle, 15cm WL
07-829-43	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, 12cm WL
07-829-44	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, 15cm WL
07-829-45	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, 12cm WL
07-829-46	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, 15.5cm WL
07-829-47	Rotatable Microscissors, Angled Jaw, 125 Degrees Up, Single Shaft Handle, 12cm WL
07-829-48	Rotatable Microscissors, Angled Jaw, 125 Degrees Up, Single Shaft Handle, 16cm WL
07-829-49MA	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, Long Blades, 12cm WL
07-829-50MA	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, Long Blades, 16cm WL



Notes

[illegible]