



LawtonElite Neurovascular Bypass Instruments

Manufactured by:

MIZUHO AMERICA



Mizuho America Inc.
30057 Ahern Avenue Union City, CA 94587
www.mizuho.com

Distributed by:

Mizuho America Inc.
30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Caution

- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

LawtonElite Neurovascular Bypass Instruments

- Available in multiple lengths with flat or round handles
- Accommodates varying depths of the surgical field
- Useful in vascular, tumor, and spine cases
- Low profile allows for control and ease of use
- Strong, Very fine, Extra fine, and Long tips
- Carbide coated jaws
- 30° & 60° up extra fine scissor tips

Reusable & Reposable Microscissors



Microdissection Set



Designed in collaboration with

MICHAEL T. LAWTON, MD

President and CEO, Barrow Neurological Institute
Professor and Chair, Neurosurgery
Chief, Neurovascular Surgery



Micro Arteriotomy Scissors

- 3.0 cm & 6.0 cm shaft lengths
- 60° up extra fine tip
- Round handles
- Guillotine blade action



Micro Sylvian Fissure Scissors

- 3.0 cm shaft length
- Straight & 30° up
- Round handle
- Flush port for cleaning



Micro Biopsy Cup Forceps

- 10.0 cm straight shaft length
- 1.0 mm spoon



Micro Scissors

- 9.5 cm working length
- 20.5 cm total length
- 60° up
- Finger pads



Lawton Scissors

- 6.5 cm working length
- 18.0 cm total length
- Straight or curved
- Round or flat handle



Micro Titanium Needle Holders

- 15.0 cm, 18.0 cm, & 21.0 cm total lengths
- 9.0 mm x 0.4 mm jaws in straight or curved
- No catches
- Round or flat handle
- Diamond dusted jaws prevent needle magnetism/rotation
- Low profile tapered tips for visibility of needle & suture



Micro Dilator Forceps

- 15.0 cm, 18.0 cm, & 21.0 cm total lengths
- Straight & 45° up
- Round or flat handle



Micro Tying Forceps

- 15.0 cm, 18.0 cm, & 21.0 cm total lengths
- 6.0 mm x 0.4 mm tying platform in straight or curved
- Round or flat handle
- Diamond dusted jaws prevent needle magnetism/rotation



Micro Forceps with DeBakey Jaw

- 15.0 cm, 18.0 cm, & 21.0 cm total lengths
- 1.5 mm jaw with DeBakey teeth
- Round handle



Micro Jeweler Forceps

- 11.0 cm, 12.0 cm, & 14.0 cm total lengths
- Strong, very fine, extra fine, & long tips



Lawton Neurovascular Bypass Organization Tray

- Two trays store instruments by size and function
- Laser-etching provides easy instrument identification
- Secure fit protects delicate instrument tips



Ordering Information

Product #	Description
07-866-00MA	Lawton Neurovascular Bypass Organization Tray
07-827-30	Microscissors Straight, Straight Jaw, Round Handle, 18cm
07-827-31	Microscissors Straight, Curved Jaw, Round Handle, 18cm
07-827-33	Microscissors Straight, Angled Jaw, Finger Pad, 9.5cm WL
07-827-34	Tenotomy Scissors, Straight Jaw, Ring Handle, 11.5cm
07-827-35	Tenotomy Scissors, Curved Jaw, Ring Handle, 11.5cm
07-827-48	Microscissors Straight, Straight Jaw, Round Handle, 3cm WL
07-827-49	Microscissors Straight, Straight Jaw, Round Handle, Irrigation, 3cm WL
07-827-50	Microscissors Straight, Angled Jaw, 60 Degrees Up, Round Handle, 6cm WL
07-827-51	Microscissors Straight, Angled Jaw, 30 Degrees Up, Round Handle, 3cm WL
07-827-52	Microscissors Straight, Angled Jaw, 60 Degrees Up, Round Handle, 3cm WL
07-862-01MA	Micro Needle Holder, Straight Jaw, 15cm
07-862-02MA	Micro Needle Holder, Curved Jaw, 15cm
07-862-03	Micro Needle Holder, Straight Jaw, 18cm
07-862-04	Micro Needle Holder, Curved Jaw, 18cm
07-862-05	Micro Needle Holder, Straight Jaw, 20.5cm
07-862-06	Micro Needle Holder, Curved Jaw, 20.5cm
07-863-01	Micro Jeweler Forceps #1, 12cm
07-863-02	Micro Jeweler Forceps #2, 12cm
07-863-03	Micro Jeweler Forceps #3, 12cm -
07-863-04	Micro Jeweler Forceps #4, 11cm
07-863-05	Micro Jeweler Forceps #5, 11cm
07-863-08	Micro Jeweler Forceps #8, 14cm
07-863-10MA	Micro Jeweler Forceps, Long Tip, 18cm
07-866-01	Micro Tying Forceps, Straight Jaw, 15cm
07-866-02	Micro Tying Forceps, Curved Jaw, 15cm
07-866-03	Micro Tying Forceps, Straight Jaw, 18cm
07-866-04	Micro Tying Forceps, Curved Jaw, 18cm
07-866-05	Micro Tying Forceps, Straight Jaw, 21cm
07-866-06	Micro Tying Forceps, Curved Jaw, 21cm
07-868-01	Micro Dilator Forceps, 15cm Straight, 0.2mm Tip
07-868-02	Micro Dilator Forceps, 15cm 45 Degrees, 0.2mm Tip
07-868-03	Micro Dilator Forceps, 18cm Straight, 0.2mm Tip
07-868-04	Micro Dilator Forceps, 18cm 45 Degrees, 0.2mm Tip
07-868-05	Micro Dilator Forceps, 21cm Straight, 0.2mm Tip
07-868-06	Micro Dilator Forceps, 21cm 45 Degrees, 0.2mm Tip
07-869-02	Biopsy Forceps, 1mm Spoon, 10cm WL
07-870-03	DeBaKey Forceps, Straight Jaw, 15cm
07-870-04	DeBaKey Forceps, Straight Jaw, 18cm
07-870-05	DeBaKey Forceps, Straight Jaw, 21cm
07-870-05	Micro Forceps, 21cm Straight, 1.5mm DeBaKey Jaw

FLAT-HANDLED INSTRUMENTS AVAILABLE

**Do not fit in set's tray*

07-827-30F	Microscissors Straight, Straight Jaw, Flat Handle, Ultrafine Blades, 18cm
07-827-31F	Microscissors Straight, Curved Jaw, Flat Handle, 18cm
07-862-01F	Micro Needle Holder, Straight Jaw, Flat Handle, 15cm
07-862-02F	Micro Needle Holder, Curved Jaw, Flat Handle, 15cm
07-862-03F	Micro Needle Holder, Straight Jaw, Flat Handle, 18cm
07-862-04F	Micro Needle Holder, Curved Jaw, Flat Handle, 18cm
07-862-05F	Micro Needle Holder, Straight Jaw, Flat Handle, 21cm
07-862-06F	Micro Needle Holder, Curved Jaw, Flat Handle, 21cm

Notes

This image shows a full page of blank, lined paper. It features approximately 20 evenly spaced horizontal grey lines across the entire width of the page, providing a guide for handwriting or typing. The background is a solid off-white color.