



Kelly Endonasal Instrument Set

All-in-One Endonasal Transsphenoidal

Manufactured by:

MIZUHO AMERICA



Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Mizuho Corporation

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN
<https://www.mizuho.co.jp>

Distributed by:

Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Caution

- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

Designed to support surgical needs for MICROSCOPIC APPROACHES WITH ENDOSCOPIC ASSISTANCE

- Low-profile devices that maximizes visibility at the surgical site
- Comprehensive set saves time with less hassle
- Designed in collaboration with Dr. Daniel F. Kelly, M.D.

Product Line-up

- | | | |
|-------------------------|------------------------------|--------------------------------|
| • Cottle Nasal Speculum | • Decker Cup Forceps | • Spoon Forceps |
| • Papavero Speculum | • Yasargil Forceps | • Feather Blades and Handle |
| • Papavero Crank Handle | • Pistol Grip Micro Scissors | • Kelly Endonasal Probe, 8Mhz |
| • Bayonetted Chisel | • Ring curettes | • Titanium Speculum (optional) |
| • Mallet | • Hemostat | |
| • Kerrison Rongeurs | • Weitlaner Retractor | |
| • Mucosal Dissectors | • Mayo scissors | |
| • Suctions | • Metzenbaum scissors | |



KELLY ENDONASAL INSTRUMENT SET TRAY 1 - OPENING



Cottle Nasal Speculum

- Short & Medium
- 76x7 mm & 90x7 mm Jaws



Papavero Speculum

- Pediatric, Short, Medium & Long
- 65x10 mm, 70x13 mm, 80x13 mm, & 90x13 mm Jaws




Bayonetted Chisel & Mallet

- 24.0 cm Working Length
- 3.0 mm Tip
- Mallet with Interchangeable Head



Kerrison Rongeurs

- 180.0 mm Working Length
- 1.0 mm, 2.0 mm & 3.0 mm Up Biting
- 40° Up Angle



Mucosal Dissectors


- Available in Left & Right Curve
- Spear-like Tip



One-Piece Suctions

- 5Fr, 6Fr, 7Fr, 8Fr, 9Fr & 10Fr
- Long (150.0 mm) & Extra Long (170.0 mm) Working Lengths

07-200-01MA
Kelly Endonasal Set
Organization Tray #1



TRAY 2 - OPERATING & CLOSING



Spatula Dissectors

- 2.4 mm Tip
- Curved Up & Down



Decker Cup Forceps

- Oval & Round Cups
- 2.0 mm, 2x5 mm & 2x6 mm Cup Sizes
- Straight & 45° Up Angle



Yasargil Forceps

- 11.5 cm Working Length
- 24.0 cm Total Length
- Bayonet



Adson Tissue Forceps

- 12.0 cm Total Length



Dressing Forceps

- 16.0 cm Total Length



Meningioma Scissors

- 9.0 cm & 11.5 cm Working Lengths
- 21.0 cm & 22.5 cm Total Lengths
- Curved & Straight Tips



Pistol Grip Micro Scissors

- 14.5 cm Working Length
- 21.0 cm Total Length
- Straight & Curved Left & Right
- Serrated Blades



Spoon Forceps

- Cushing, Love-Gruenwald, Spurling Types
- 16.0 cm Working Length
- 2x10 mm, 3x10 mm & 4x10 mm Cup Sizes
- Straight Angle



Hemostat

- 16.0 cm Total Length



Dissectors

Teardrop Probe

- 2.5 mm Tip
- Straight, 45°, 60° & 90°
- Straight, Up & Down Angles



Ring Curettes

- 12.5 cm Working Length
- 1.5 mm, 3.0 mm, 5.0 mm & 7.0 mm Tip Sizes
- 45° & 90° Up Angle
- 45° Down Angle



Mezenbaum Scissors

- 18.0 cm Total Length
- Straight & Curved




Mayo Scissors

- 17.5 cm Total Length
- Straight



Weitlaner Retractor

- 13.0 cm Total Length



07-200-02MA
Kelly Endonasal Set
Organization Tray #2

Ordering Information

Product #	Description	TRAY #1
07-200-01MA	Kelly Endonasal Set Organization Tray #1	
07-177-10	Cottle Nasal Speculum, 76x7mm Jaw	
07-177-11	Cottle Nasal Speculum, 90x7mm Jaw	
07-907-04	Papavero Speculum, 65x10mm Jaw	
07-907-05	Papavero Speculum, 70x13mm Jaw	
07-907-06	Papavero Speculum, 80x13mm Jaw	
07-907-07	Papavero Speculum, 90x13mm Jaw	
07-907-09	Papavero Crank Handle	
01-413-10	Mallet, 18cm	
02-034-10	Spear-Like Dissector, Left Curve, Bayonet, 12.5cm WL	
02-034-11	Spear-Like Dissector, Right Curve, Bayonet, 12.5cm WL	
02-201-01	Kerrison rongeur, upward, 40 degree, 180mm(L)x1mm(W)	
02-201-02	Kerrison rongeur, upward, 40 degree, 180mm(L)x2mm(W)	
02-201-03	Kerrison rongeur, upward, 40 degree, 180mm(L)x3mm(W)	
07-923-01	Chisel, Bayonet, 3mm Jaw, 24cm TL	
09-081-10	Metzenbaum Scissors, 14cm Curved	
07-960-05L	Tapered Suction Tube, One Piece, 5Fr, 11cm WL	
07-960-05XL	Tapered Suction Tube, One Piece, 5Fr, 14.5cm WL	
07-960-06L	Tapered Suction Tube, One Piece, 6Fr, 11cm WL	
07-960-06XL	Tapered Suction Tube, One Piece, 6Fr, 14.5cm WL	
07-960-07L	Tapered Suction Tube, One Piece, 7Fr, 11cm WL	
07-960-07XL	Tapered Suction Tube, One Piece, 7Fr, 14.5cm WL	
07-960-08L	Tapered Suction Tube, One Piece, 8Fr, 11cm WL	
07-960-08XL	Tapered Suction Tube, One Piece, 8Fr, 14.5cm WL	
07-960-09L	Tapered Suction Tube, One Piece, 9Fr, 11cm WL	
07-960-09XL	Tapered Suction Tube, One Piece, 9Fr, 14.5cm WL	
07-960-10L	Tapered Suction Tube, One Piece, 10Fr, 11cm WL	
07-960-10XL	Tapered Suction Tube, One Piece, 10Fr, 14.5cm WL	

Product #	Description	OR SUITE
108910	8MHz Vascular Doppler Transceiver, Battery Only	
138210	8MHz Kelly Endonasal Probe (Box of 4)	
K-5310	Straight Blade 1mm (Box of 5)	
K-5520	Hook Blade 30° 1mm (Box of 5)	

Product #	Description	ADDITIONAL OPTIONS
07-917-04MA	Papavero Speculum, Titanium, 65x13mm Jaw	
07-917-05MA	Papavero Speculum, Titanium, 70x13mm Jaw	
07-917-06MA	Papavero Speculum, Titanium, 80x13mm Jaw	
07-917-07MA	Papavero Speculum, Titanium, 90x13mm Jaw	
07-927-02	Trapezoidal Speculum, 45x12mm Jaw	
07-927-03	Trapezoidal Speculum, 60x13mm Jaw, Down	
07-927-04	Trapezoidal Speculum, 60x13mm Jaw, Up	
07-927-05MA	Trapezoidal Speculum, 70x13mm Jaw, Down	
07-927-06MA	Trapezoidal Speculum, 70x13mm Jaw, Up	
07-962-05D	Teardrop Suction Tube, 5Fr, Atraumatic, Angled Down, 10.5cm WL	
07-962-05U	Teardrop Suction Tube, 5Fr, Atraumatic, Angled Up, 10.5cm WL	
07-962-07D	Teardrop Suction Tube, 7Fr, Atraumatic, Angled Down, 10.5cm WL	
07-962-07U	Teardrop Suction Tube, 7Fr, Atraumatic, Angled Up, 10.5cm WL	

Product #	Description	TRAY #2
07-200-02MA	Kelly Endonasal Set Organization Tray #2	
02-034-12	Spatula Dissector, 2.4mm, Curved Up, Bayonet, 12.5cm WL	
02-034-13	Spatula Dissector, 2.4mm, Curved Down, Bayonet, 12.5cm WL	
02-082-09	Decker Rongeur, Oval 2x6mm, 15cm WL	
02-082-20	Decker Rongeur, Round 2mm, 45 Degrees Up, 14cm WL	
02-082-21	Decker Rongeur, Round 2mm, 14cm WL	
02-082-22	Decker Rongeur, Oval 2x5mm, 45 Degrees Up, 14cm WL	
02-082-23	Decker Rongeur, Oval 2x5mm, 45 Degrees Up, 14.5cm WL	
02-202-07	Cushing Rongeur, Oval 2x10mm, 16cm WL	
02-202-08	Love-Gruenwald Rongeur, Oval 3x10mm, 16cm WL	
02-202-09	Spurling Rongeur, Oval 4x10mm, 16cm WL	
07-110-10	Yasargil Forceps, Bayonet, 11.5cm WL	
07-134-10	Endoscopic Microscissors, Straight Jaw, Pistol-Grip, 14cm WL	
07-134-11	Endoscopic Microscissors, Curved Right Jaw, Pistol-Grip, 14cm WL	
07-134-12	Endoscopic Microscissors, Curved Left Jaw, Pistol-Grip, 14cm WL	
07-134-13	Endoscopic Microscissors, Curved Up Jaw, Pistol-Grip, 14cm WL	
07-827-09	Microscissors Bayonet, Straight Jaw, Flat Handle, 9cm WL	
07-827-10	Microscissors Bayonet, Curved Jaw, Flat Handle, 9cm WL	
07-820-24	Teardrop Probe, 2.5mm Tip, Straight, Bayonet, 12.5cm WL	
07-820-26	Teardrop Probe, 2.5mm Tip, 45 Degrees Up, Bayonet, 12.5cm WL	
07-820-27	Teardrop Probe, 2.5mm Tip, 90 Degrees Up, Bayonet, 12.5cm WL	
07-820-29	Teardrop Probe, 2.5mm Tip, 60 Degrees Down, Bayonet, 12.5cm WL	
07-822-21	Ring Curette, 1.5mm, 45 Degrees Up, Bayonet, 12.5cm WL	
07-822-23	Ring Curette, 3mm, 45 Degrees Up, Bayonet, 12.5cm WL	
07-822-25	Ring Curette, 5mm, 45 Degrees Up, Bayonet, 12.5cm WL	
07-822-27	Ring Curette, 7mm, 45 Degrees Up, Bayonet, 12.5cm WL	
07-822-31	Ring Curette, 1.5mm, 90 Degrees Up, Bayonet, 12.5cm WL	
07-822-33	Ring Curette, 3mm, 90 Degrees Up, Bayonet, 12.5cm WL	
07-822-35	Ring Curette, 5mm, 90 Degrees Up, Bayonet, 12.5cm WL	
07-822-37	Ring Curette, 7mm, 90 Degrees Up, Bayonet, 12.5cm WL	
07-822-41	Ring Curette, 1.5mm, 45 Degrees Down, Bayonet, 12.5cm WL	
07-822-43	Ring Curette, 3mm, 45 Degrees Down, Bayonet, 12.5cm WL	
07-822-45	Ring Curette, 5mm, 45 Degrees Down, Bayonet, 12.5cm WL	
07-822-47	Ring Curette, 7mm, 45 Degrees Down, Bayonet, 12.5cm WL	
07-822-97	Ring Curette, 7mm, Left to Right, Bayonet, 12.5cm WL	
06-674-00	Hemostat, 15cm	
06-826-00	Mayo-Hegar Needle Holder, 16cm	
07-057-10	Weitlaner Retractor, 13cm	
07-120-10	#3 Scalpel Handle	
07-181-01	Adson Tissue Forceps, 12cm	
07-181-10	Dressing Forceps, 16cm	
09-061-10	Mayo Scissors, 17cm Straight	
09-071-10	Metzenbaum Scissors, 18cm Straight	
09-071-11	Metzenbaum Scissors, 18cm Curved	
09-214-10	Allis Forceps, Straight Jaw, 15cm	
09-071-20F	Microscissors Bayonet, Straight Jaw, Flat Handle, Round Blades, 9cm WL	
09-071-21F	Microscissors Bayonet, Curved Jaw, Flat Handle, Round Blades, 9cm WL	
09-071-22F	Microscissors Bayonet, Straight Jaw, Flat Handle, Round Blades, 11.5cm WL	
09-071-23F	Microscissors Bayonet, Curved Jaw, Flat Handle, Round Blades, 11.5cm WL	
NF-100BA	Titanium Blade Handle Rotatable Bayonet 100mm	

Notes

[illegible]