



Kelly Keyhole Instrument Set

Manufactured by:

MIZUHO AMERICA



Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Mizuho Corporation

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN
<https://www.mizuho.co.jp>

Distributed by:

Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Caution

- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

Designed to support surgical needs for

KEYHOLE APPROACHES

Useful for:

- Supraorbital “Eyebrow”
- Supracerebellar-transtentorial
- Retromastoid
- Mini-pterional
- Endonasal
- Designed in collaboration with Dr. Daniel F. Kelly, M.D.

Product Line-up

- Suctions
- Duckbill Spatulas
- Spatulas
- Tear Drop Nerve Hook
- Rhoton Dissectors
- Arachnoid Knives
- Micro Forceps
- Ring curettes
- Tumor Forceps
- Evans Rotatable Micro Scissors
- Decker Rongeurs

Tray 1



Day-Bailes Suctions and Handle

- Flat Tear-drop Handle
- 5Fr, 6Fr, 7Fr, 8Fr, 9Fr
- Medium and Long Working Length



Kelly Endonasal Suctions

- 5Fr, 7Fr
- 15cm Working Length
- Up & Down Angle



Ring Curettes


- 9cm & 11cm Working Lengths
- 1.5mm, 3mm, & 5mm Diameter
- 45° & 90° Up
- 45° & 90° Down

Tray 2



Duckbill Spatulas

- 9cm & 11cm Working Lengths
- Up & Down Angle




Rhoton #3 Dissectors

- 3mm Tip
- 9cm & 11cm Working Lengths
- 30° Up & Down



Spatulas

- 1.8mm Tip
- 9cm Working Length
- Angled Up & Down



Arachnoid Knives

- 11cm Working Length
- Angle Up & Down



Tear Drop Nerve Hooks

- 9cm & 11cm Working Lengths
- Straight
- 45° Up & Down
- 90° Up & Down



Tumor Forceps

- 9cm Working Length
- 3mm & 5mm Diameter
- Bayonet, flat handle design




Micro Forceps

- 9cm & 11cm Working Lengths
- Bayonet Handle
- 0.4mm (sharp) & 0.9mm (blunt) Tip Size




Evans Rotatable Micro Scissors

- 9cm & 11cm Working Lengths
- Round Handle Design
- Straight
- 45° Up
- Curved, Double Action



Decker Rongeurs

- 11cm Working Length
- 2mm & 2x5mm Cup Sizes
- Round & Oval Cups
- 45° & Straight Tips



Ordering Information

Product #	Description	TRAY #1	Product #	Description	TRAY #2
07-822-00MA	Kelly Keyhole Instruments Organization Tray 1		07-822-01MA	Kelly Keyhole Instruments – Organization Tray 2	
07-119-01SC	Flat Teardrop Suction Handle w/Safety Coupling		02-034-01MA	Rhoton #3, 9cm WL, 3mm 30 Degrees Up	
07-822-48MA	Ring Curette 9cm WL, 5mm 45 Degrees Up		02-034-02MA	Rhoton #3, 9cm WL, 3mm 30 Degrees Down	
07-822-49MA	Ring Curette 9cm WL, 5mm 90 Degrees Up		02-034-03MA	Rhoton #3, 11cm WL, 3mm 30 Degrees Up	
07-822-50MA	Ring Curette 9cm WL, 5mm 45 Degrees Down		02-034-04MA	Rhoton #3, 11cm WL, 3mm 30 Degrees Down	
07-822-51MA	Ring Curette 11cm WL, 5mm 45 Degrees Up		02-034-05MA	Duckbill Spatula, 9cm WL, Up Angle	
07-822-52MA	Ring Curette 11cm WL, 5mm 90 Degrees Up		02-034-06MA	Duckbill Spatula, 9cm WL, Down Angle	
07-822-53MA	Ring Curette 11cm WL, 5mm 45 Degrees Down		02-034-07MA	Duckbill Spatula, 11cm WL, Up Angle	
07-822-54MA	Ring Curette 9cm WL, 3mm 45 Degrees Up		02-034-08MA	Duckbill Spatula, 11cm WL, Down Angle	
07-822-55MA	Ring Curette 9cm WL, 3mm 90 Degrees Up		02-034-14MA	Spatula 1.8mm, 125mm WL up	
07-822-56MA	Ring Curette 9cm WL, 3mm 45 Degrees Down		02-034-15MA	Spatula 1.8mm, 125mm WL dwn	
07-822-57MA	Ring Curette 11cm WL, 3mm 45 Degrees Up		02-082-40	Micro Forceps, 9cm WL Straight, Bayonet Sharp 0.4mm, Flat Handle	
07-822-58MA	Ring Curette 11cm WL, 3mm 90 Degrees Up		02-082-41	Micro Forceps, 9cm WL Straight, Bayonet Blunt 0.9mm, Flat Handle	
07-822-59MA	Ring Curette 11cm WL, 3mm 45 Degrees Down		02-082-42	Micro Forceps, 11cm WL Straight, Bayonet Sharp 0.4mm, Flat Handle	
07-822-60MA	Ring Curette 9cm WL, 1.5mm 45 Degrees Up		02-082-43	Micro Forceps, 11cm WL Straight, Bayonet Blunt 0.9mm, Flat Handle	
07-822-61MA	Ring Curette 9cm WL, 1.5mm 90 Degrees Up		02-082-44	Yasargil Tumor Forceps, 9cm WL, Bayonet 3mm Dia, Flat Handle	
07-822-62MA	Ring Curette 9cm WL, 1.5mm 45 Degrees Down		02-082-45	Yasargil Tumor Forceps, 9cm WL, Bayonet 5mm Dia, Flat Handle	
07-822-63MA	Ring Curette 11cm WL, 1.5mm 45 Degrees Up		02-082-46	Decker Rongeur, 11cm WL, 2mm Round Cup, 45 Degrees Up	
07-822-64MA	Ring Curette 11cm WL, 1.5mm 90 Degrees Up		02-082-47	Decker Rongeur, 11cm WL, 2mm Round Cup, Straight	
07-822-65MA	Ring Curette 11cm WL, 1.5mm 45 Degrees Down		02-082-48	Decker Rongeur, 11cm WL, 2x5mm Oval Cup, Straight	
07-962-05D	Kelly Endonasal Suction Tube, 5Fr, 15cm Down		02-082-49	Kelly Keyhole Decker Rongeur, 11cm WL, 2x5mm Oval Cup, 45 Degrees Up	
07-962-05U	Kelly Endonasal Suction Tube, 5Fr, 15cm Up		07-820-25MA	Teardrop Nerve Hook, 9cm WL Straight	
07-962-07D	Kelly Endonasal Suction Tube, 7Fr, 15cm Down		07-820-28MA	Teardrop Nerve Hook, 9cm WL 45 Degrees Up	
07-962-07U	Kelly Endonasal Suction Tube, 7Fr, 15cm Up		07-820-30MA	Teardrop Nerve Hook, 9cm WL 90 Degrees Up	
07-DAY-05LSC	Day-Bailes Locking Suction 5Fr Long		07-820-31MA	Teardrop Nerve Hook, 9cm WL 45 Degrees Down	
07-DAY-05MSC	Day-Bailes Locking Suction 5Fr Medium		07-820-32MA	Teardrop Nerve Hook, 11cm WL Straight	
07-DAY-06LSC	Day-Bailes Locking Suction 6Fr Long		07-820-33MA	Teardrop Nerve Hook, 11cm WL 45 Degrees Up	
07-DAY-06MSC	Day-Bailes Locking Suction 6 Fr Medium, 125MM W/STYLET		07-820-34MA	Teardrop Nerve Hook, 11cm WL 90 Degrees Up	
07-DAY-07LSC	Day-Bailes Locking Suction 7Fr Long		07-820-35MA	Teardrop Nerve Hook, 11cm WL 45 Degrees Down	
07-DAY-07MSC	Day-Bailes Locking Suction 7Fr Medium		07-821-14MA	Arachnoid Knife CP-2, 11cm WL Up Angle	
07-DAY-08LSC	Day-Bailes Locking Suction 8Fr Long		07-821-15MA	Arachnoid Knife CP-2, 11cm WL Down Angle	
07-DAY-08MSC	Day-Bailes Locking Suction 8Fr Medium		07-821-16MA	Arachnoid Knife CP-4, 11cm WL Up Angle	
07-DAY-09LSC	Day-Bailes Locking Suction 9Fr Long		07-821-17MA	Arachnoid Knife CP-4, 11cm WL Down Angle	
07-DAY-09MSC	Day-Bailes Locking Suction 9Fr Medium		07-827-61	Kelly Keyhole Rotatable Microscissors, 23cm, 9cm WL Straight	
			07-827-62	Kelly Keyhole Rotatable Microscissors, 25cm, 11cm WL Straight	
			07-827-63	Kelly Keyhole Rotatable Microscissors, 23cm, 9cm WL Curved, Double-Action	
			07-827-64	Kelly Keyhole Rotatable Microscissors, 25cm, 11cm WL Curved, Double-Action	
			07-827-65	Kelly Keyhole Rotatable Microscissors, 23cm, 9cm WL 45 Degrees Up	
			07-827-66	Kelly Keyhole Rotatable Microscissors, 25cm, 11cm WL 45 Degrees Up	

Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.