

Designed in collaboration with

James J. Evans, MD

Professor
Chief, Brain Tumor and Stereotactic Radiosurgery Division
Director, Cranial Base and Pituitary Surgery
Director, Cranial Base and Endoscopic Surgery Fellowship



Dr. Evans is the Chief of the Brain Tumor Division at Thomas Jefferson University and Past President of the North American Skull Base Society. He is also Co-Director of the Minimally Invasive Cranial Base Surgery and Endoscopic Neurosurgery Center at Jefferson. Dr. Evans performs complex cranial, skull base, and endoscopic minimally invasive neurosurgery to treat a wide range of tumors and disorders of the brain and skull base, with a special emphasis on pituitary tumors. He founded the cranial base tumor bank at Jefferson and performs molecular genetics research on pituitary and cranial base tumors. In addition, he is the Director of the Cranial Base and Endoscopic Surgery Fellowship, a training program that has produced many new leaders and contributors to the field. Dr. Evans has authored three textbooks, published over 200 manuscripts and chapters, and has been named in Philadelphia “Top Docs” annually for over the past 10 years.



Ordering Information

360° Rotatables

Product #	Description
02-082-00MA	Evans Rotatable Instruments Organization Tray
07-829-01	Rotatable Micro Needle Holder, Straight Jaw, 12cm WL
07-829-02	Rotatable Micro Needle Holder, Straight Jaw, 15.5cm WL
07-829-03	Rotatable Micro Needle Holder, Curved Jaw, 12cm WL
07-829-04	Rotatable Micro Needle Holder, Curved Jaw, 15.5cm WL
07-829-11	Rotatable Biopsy Forceps, 2mm Spoon, 12cm WL
07-829-12	Rotatable Biopsy Forceps, 2mm Spoon, 15cm WL
07-829-13	Rotatable Biopsy Forceps, 45 Degrees Up, 2mm Spoon, 12cm WL
07-829-14	Rotatable Biopsy Forceps, 45 Degrees Up, 2mm Spoon, 15cm WL
07-829-15	Rotatable Biopsy Forceps, 45 Degrees Up, 2x5mm Spoon, 12cm WL
07-829-16	Rotatable Biopsy Forceps, 45 Degrees Up, 2x5mm Spoon, 15cm WL
07-829-17	Rotatable Biopsy Forceps, 2x5mm Spoon, 12cm WL
07-829-18	Rotatable Biopsy Forceps, 2x5mm Spoon, 15cm WL
07-829-19	Rotatable Biopsy Forceps, 4mm Spoon, 12cm WL
07-829-20	Rotatable Biopsy Forceps, 4mm Spoon, 16cm WL
07-829-21	Rotatable Biopsy Forceps, 45 Degrees Up, 4mm Spoon, 12cm WL
07-829-22	Rotatable Biopsy Forceps, 45 Degrees Up, 4mm Spoon, 16cm WL
07-829-31	Rotatable Grasping Forceps, Straight Jaw, 11.5cm WL
07-829-32	Rotatable Grasping Forceps, Straight Jaw, 15.5cm WL
07-829-33	Rotatable Grasping Forceps, Angled Jaw, 11.5cm WL
07-829-34	Rotatable Grasping Forceps, Angled Jaw, 15.5cm WL
07-829-41	Rotatable Microscissors, Angled Jaw, 45 Degrees Up, Single Shaft Handle, 11.5cm WL
07-829-42	Rotatable Microscissors, Angled Jaw, 45 Degrees Up, Single Shaft Handle, 15cm WL
07-829-43	Rotatable Microscissors, Straight Jaw, Single Shaft Handle, 12cm WL
07-829-44	Rotatable Microscissors, Straight Jaw, Single Shaft Handle, 15cm WL
07-829-45	Rotatable Microscissors, Curved Jaw, Single Shaft Handle, 12cm WL
07-829-46	Rotatable Microscissors, Curved Jaw, Single Shaft Handle, 15.5cm WL
07-829-47	Rotatable Microscissors, Angled Jaw, 125 Degrees Up, Single Shaft Handle, 12cm WL
07-829-48	Rotatable Microscissors, Angled Jaw, 125 Degrees Up, Single Shaft Handle, 16cm WL
07-829-49MA	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, Long Blades, 12cm WL
07-829-50MA	Rotatable Microscissors, Angled Jaw, Single Shaft Handle, Long Blades, 16cm WL
P07-829-00	Flushing Port, Detachable, Replacement

Evans Elite Endoscopic Set

10293MA	Evans Endonasal Set Tray
---------	--------------------------

Nasal Forceps

07-181-13MA	Nasal Forceps, Straight Jaw, 12cm WL
07-181-14MA	Nasal Forceps, 45 Degrees Up, 12cm WL

Ring Curettes

10031MA	Ring Curette, Donut Style, 1.5mm, Straight, Malleable, 15cm WL
10032MA	Ring Curette, Donut Style, 3mm, Straight, Malleable, 15cm WL
10033MA	Ring Curette, Donut Style, 5mm, Straight, Malleable, 15cm WL
10034MA	Ring Curette, Donut Style, 7mm, Straight, Malleable, 15cm WL
10035MA	Ring Curette, Donut Style, 1.5mm, 45 Degrees Up, 15cm WL
10036MA	Ring Curette, Donut Style, 3mm, 45 Degrees Up, 15cm WL
10037MA	Ring Curette, Donut Style, 5mm, 45 Degrees Up, 15cm WL
10038MA	Ring Curette, Donut Style, 7mm, 45 Degrees Up, 15cm WL
10039MA	Ring Curette, Donut Style, 1.5mm, 90 Degrees Up, 15cm WL
10040MA	Ring Curette, Donut Style, 3mm, 90 Degrees Up, 15cm WL
10041MA	Ring Curette, Donut Style, 5mm, 90 Degrees Up, 15cm WL
10197MA	Ring Curette, Donut Style, 1.5mm, Offset, 15cm WL
10198MA	Ring Curette, Donut Style, 3mm, Offset, 15cm WL
10199MA	Ring Curette, Donut Style, 5mm, Offset, 15cm WL

10216MA	Ring Curette, Donut Style, 3mm, 45 Degrees Lateral, 15cm WL
10217MA	Ring Curette, Donut Style, 5mm, 45 Degrees Lateral, 15cm WL
10218MA	Ring Curette, Donut Style, 7mm, 45 Degrees Lateral, 15cm WL
10219MA	Ring Curette, Donut Style, 5mm, Curved 45 Degrees Lateral, 15cm WL
10220MA	Ring Curette, Donut Style, 7mm, Curved 45 Degrees Lateral, 15cm WL
10223MA	Ring Curette, Donut Style, 3mm, Curved 45 Degrees Lateral, 15cm WL
10296MA	Ring Curette, Donut Style, 3mm, 125 Degrees Up, 15cm WL
10297MA	Ring Curette, Donut Style, 5mm, 125 Degrees Up, 15cm WL

Atraumatic Suctions

10042MA	Suction Tube, 4Fr, Atraumatic, 8.8cm WL
10043MA	Suction Tube, 6Fr, Atraumatic, 8.8cm WL
10044MA	Suction Tube, 8Fr, Atraumatic, 8.8cm WL
10045MA	Suction Tube, 10Fr, Atraumatic, 8.8cm WL
10046MA	Suction Tube, 4Fr, Atraumatic, 11.3cm WL
10047MA	Suction Tube, 6Fr, Atraumatic, 11.3cm WL
10048MA	Suction Tube, 8Fr, Atraumatic, 11.3cm WL
10049MA	Suction Tube, 10Fr, Atraumatic, 11.3cm WL
10050MA	Suction Tube, 4Fr, Atraumatic, 13.3cm WL
10051MA	Suction Tube, 6Fr, Atraumatic, 13.3cm WL
10052MA	Suction Tube, 8Fr, Atraumatic, 13.3cm WL
10053MA	Suction Tube, 10Fr, Atraumatic, 13.3cm WL
07-DAY-025C	Suction Tube Handle, Flat, Teardrop Control, 8.5cm TL

Malleable Tapered Suctions

10054MA	Tapered Suction Tube, 6Fr, Malleable, 8.8cm WL
10055MA	Tapered Suction Tube, 8Fr, Malleable, 8.8cm WL
10056MA	Tapered Suction Tube, 10Fr, Malleable, 8.8cm WL
10057MA	Tapered Suction Tube, 12Fr, Malleable, 8.8cm WL
10058MA	Tapered Suction Tube, 6Fr, Malleable, 11.3cm WL
10059MA	Tapered Suction Tube, 8Fr, Malleable, 11.3cm WL
10060MA	Tapered Suction Tube, 10Fr, Malleable, 11.3cm WL
10061MA	Tapered Suction Tube, 12Fr, Malleable, 11.3cm WL
10224MA	Tapered Suction Tube, 6Fr, Malleable, 13.3cm WL
10225MA	Tapered Suction Tube, 8Fr, Malleable, 13.3cm WL
10226MA	Tapered Suction Tube, 10Fr, Malleable, 13.3cm WL
10227MA	Tapered Suction Tube, 12Fr, Malleable, 13.3cm WL

Round Knives

10200MA	Round Knife 2mm, 45 Degrees Up, Straight, 15cm WL
10201MA	Round Knife 3mm, 45 Degrees Up, Straight, 15cm WL

Elevator

10202MA	Elevator, 1.5mm Tip, 15 Degrees Up, 15cm WL
---------	---

Spatulas

10221MA	Spatula Dissector, 2mm Tip, 30 Degrees Up, 15cm WL
10222MA	Spatula Dissector, 3mm Tip, 30 Degrees Up, 15cm WL

Kerrison Rongeurs

7375-41E	Kerrison Rongeur, Rotatable, 40 Degrees, 1mm Tip
7375-42E	Kerrison Rongeur, Rotatable, 40 Degrees, 2mm Tip
7375-43E	Kerrison Rongeur, Rotatable, 40 Degrees, 3mm Tip
7375-62	Kerrison Rongeur, Rotatable, 125 Degrees, 2mm Tip
7375-63	Kerrison Rongeur, Rotatable, 125 Degrees, 3mm Tip

Manufactured by:

MIZUHO AMERICA



Mizuho America Inc.
30057 Ahern Avenue Union City, CA 94587
www.mizuho.com

Distributed by:

Mizuho America Inc.
30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



EVANS ELITE ENDOSCOPIC INSTRUMENTS



The Gold Standard



EVANS ELITE ENDOSCOPIC INSTRUMENTS



- Designed as a comprehensive set for endoscopic and minimally invasive procedures
- Malleable instruments and suctions provide more versatility to surgeons
- Rotating Kerrison Rongeurs allow multiple positions for optimal orientation
- Low profile instruments improve visualization through keyhole access




Back-Biting Kerrison Rongeurs

- Rotatable Shaft
- 125° Back Cutting
- 2.0 mm & 3.0 mm Bite Width
- 7.0" Working Length




Kerrison Rongeurs

- Rotatable Shaft
- 40° Up Angle
- 1.0 mm, 2.0 mm & 3.0 mm Bite Width
- 7.0" Working Length




Round Knives

- 45° Angle
- Sharp
- 2.0 mm (Small) & 3.0 mm (Large)
- 15.0 cm Working Length



Elevator


- 15° Angle
- 1.5 mm Tip
- 15.0 cm Working Length



Spatula


- 30° Angle
- 2.0 mm (Small) & 3.0 mm (Large)
- 15.0 cm Working Length

The Gold Standard




Ring Curettes

- Donut Style (Blunt)
- 1.5 mm (Micro), 3.0 mm (Small), 5.0 mm (Medium) & 7.0 mm (Large) Diameters
- Straight Malleable
- 45°, 90° & 125° Angles
- Lateral, Curved Lateral & Offset
- 15.0 cm Working Length



Atraumatic Suctions

- Football Tip
- 4Fr, 6Fr, 8Fr & 10Fr
- 12.5 cm (Medium), 15.0 cm (Long) & 17.0 cm (Extra Long)




Malleable Suctions

- Tapered
- 6Fr, 8Fr, 10Fr & 12Fr
- 15.0 cm (Long) & 17.0 cm (Extra Long)

E3 SERIES INSTRUMENT TRAYS


10293MA

Evans Endoscopic Set Tray



02-082-00MA


Evans Rotatable Instruments Organization Tray



360° ROTATABLES




- 360° Fully Rotatable Shaft
- Axial Bayonet for Improved Visualization
- Round Handles for Easier Grip
- Safely Prevents Rotation when Jaws are Closed
- Balanced for Comfort




Microscissors

- Straight, 45° Up, Curved, 125° Back Cutting & Extended Blades
- 11.5 cm, 12.0 cm, 15.5 cm & 16.0 cm Working Lengths
- Serrated Blades




Micro Biopsy Cup Forceps

- Straight or 45° Up
- 2.0 mm & 4.0 mm Round Cup
- 2.0 x 5.0 mm Oval Cup
- 11.0 cm & 15.5 cm Working Lengths



Micro Grasping Forceps

- Straight or Angled Up
- 11.5 cm & 15.5 cm Working Lengths



Micro Needle Holders

- Straight or Curved
- 12.0 cm & 15.5 cm Working Lengths