



Kelly Endoscopic Endonasal Instrument Set

Manufactured by:
MIZUHO AMERICA

Mizuho America Inc.
 30057 Ahern Avenue Union City, CA 94587
www.mizuho.com



Mizuho Corporation
 3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN
<https://www.mizuho.co.jp>

Distributed by:
Mizuho America Inc.
 30057 Ahern Avenue Union City, CA 94587
www.mizuho.com

 Caution	•Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information. •Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations. •For after-sale support, contact your local sales representative.
---	--

Designed to support surgical needs for

ENDOSCOPIC APPROACHES

- Low- profile devices that maximizes visibility at the surgical site
- Comprehensive set saves time with less hassle
- Designed in collaboration with Dr. Daniel F. Kelly, M.D.

Product Line-up

- | | | |
|---|---|--|
| <ul style="list-style-type: none"> • Cottle Elevator • Probes • Nasal Forceps • Antrum Punch • Nasal Scissors • Suctions • Antrum Cannula • Weil Blakesley Ronguers | <ul style="list-style-type: none"> • Kerrison Rongeurs • Rhoton Dissectors • Decker Cup Forceps • Pistol Grip Micro Scissors • Spoon Forceps • Ring Curettes • Hemostat • Mayo Scissors | <ul style="list-style-type: none"> • Weitlaner Retractor • Metzenbaum Scissors • Spoon Forceps • Feather Blades and Handle |
|---|---|--|

Tray 1



Cottle Elevators

- Double-ended
- 20cm Total Length



Probes

- Double-ended
- 19cm Total Length



Rhoton #3 Dissectors

- 12cm Working Length
- 3mm Tip
- 45° Up & Down
- Bayonet



Mucosal Dissector

- Left & Right Angle
- Spear-like Tip



Nasal Forceps

- 13cm Working Length
- Straight & 45°Up



Antrum Punch

- 360° Rotatable
- Back-cutting
- 10cm Working Length



Nasal Scissors

- Pistol Grip handle
- 13cm Working Length
- Straight, Curved Right and Left



Antrum Cannulas

- 14.5cm Total Length
- 2.5mm, 3.0mm & 4.9mm Diameter



Kerrison Rongeurs

- 178mm & 180mm Working Length
- Angle Up & Down
- 40° Angle
- 1mm, 2mm & 3mm Biting



One-Piece Suctions

- 5Fr, 6Fr, 7Fr, 8Fr, 9Fr, and 10Fr
- Long (150mm) & Extra Long (170mm)



Kelly Endoscopic Endonasal Instruments Organization Tray 1
07-200-03MA

Tray 2



Decker Cup Forceps

- Oval & Round Cups
- 2mm, 2x5mm & 6mm Cup Sizes
- Straight and 45° Up
- 14.5cm & 15cmn Working Length



Micro Scissors

- Pistol Grip Handle
- 14.5cm Working Length
- Straight and Curved Right & Left



Spoon Forceps

- Cushing, Love-Gruenwald, Spurling Types
- 178mm Working Length
- 2x10mm, 3x10mm, 4x10mm Cup Sizes
- Straight Angle



Dissectors

Probe teardrop

- 2.5mm Tip
- Straight, 45°, 60° & 90°
- Straight, Up & Down Angles



Ring Curettes

- 1.5mm, 3mm, 5mm, 7mm Cup Sizes
- 45° & 90° Up
- 45° Down



Spatula Dissectors

- 2.4mm Tip
- Curved Up and Down



Hemostat

- 15cm Total Length



Mezenbaum Scissors

- 18cm Total Length
- Straight & Curved



Mayo Scissors

- 17cm Total Length
- Straight



Weitlaner Retractor

- 13cm Total Length



Adson Tissue Forcep

- 12cm Total Length
- 1x2mm Teeth



Dressing Forcep

- 16cm Total Length

Instruments Organization Trays



Kelly Endoscopic Endonasal Instruments Organization Tray 2
07-200-04MA

Ordering Information

Product #	Description	TRAY #1	Product #	Description	TRAY #2
07-200-03MA	Kelly Endoscopic Endonasal Instruments Organization Tray 1		07-200-04MA	Kelly Endoscopic Endonasal Instruments - Tray 1 Lid & Top Tray Only	
02-034-09MA	Rhoton #3, Round Bayonetted Dissector 12cm WL, 3mm 45 Degree Up		02-202-07	Cushing, Spoon Jaw 2x10mm, 178mm Straight	
02-034-10MA	Rhoton #3, Round Bayonetted Dissector 12cm WL, 3mm 45 Degrees Down		02-202-08	Love-Gruenwald, Spoon Jaw 3x10mm, 178mm Straight	
●07-181-11MA	Cottle Elevator, Double-Ended 20cm		02-202-09	Spurling, Spoon Jaw 4x10mm, 178mm Straight	
07-181-12MA	Probe, Double-Ended 19cm		07-110-10	Kelly Endonasal Yasargil Forceps, Bayonet, 23.5cm	
07-181-13MA	Nasal Forceps, 13cm WL Straight		NF-100BA	Titanium Blade Handle Rotatable Bayonet 100mm	
07-181-14MA	Nasal Forceps, 13cm WL 45 Degrees Up		02-082-09	Decker Cup Forceps, Oval 2x6mm, 15cm Straight	
07-181-15MA	Antrum Punch, 10cm WL, Backward Cutting, 360 Degrees Rotatable		02-082-20	Decker Cup Forceps, Round 2mm, 14.5cm 45 Degrees Up	
07-181-16MA	Antrum Punch, 10cm WL, Backward Cutting, 360 Degrees Rotatable		02-082-21	Decker Cup Forceps, Round 2mm, 14.5cm Straight	
07-181-17MA	Nasal Scissors, 13cm WL Straight		02-082-22	Kelly Decker Cup Forceps, Oval 2x5mm, 14.5cm 45 Degrees Up	
07-181-18MA	Kelly EndoNasal Nasal Scissors, 13cm WL Curved to Right		02-082-23	Kelly Endonasal Decker Cup Forceps, Oval 2x5mm, 14.5cm Straight	
07-181-19MA	Kelly Endonasal Nasal Scissors, 13cm WL Curved to Left		07-134-10	Kelly Endonasal Pistol Grip Micro Scissors, Straight 14.5cm WL	
07-181-20MA	Antrum Cannula w/Tear Drop Suction Control, 14.5cm Long Curved, 2.5mm Outer Dia		07-134-11	Kelly, Endonasal Pistol Grip Micro Scissors, Curved Right 14.5cm WL	
07-181-21MA	Antrum Cannula w/Tear Drop Suction Control, 14.5cm Long Curved, 3.0mm Outer Dia.		07-134-12	Kelly Endonasal Pistol Grip Micro Scissors, Curved Left 14.5cm	
07-181-22MA	Antrum Cannula w/Tear Drop Suction Control, 14.5cm Long Curved, 4.0mm Outer Dia.		06-674-00	Kelly Endoscopic Endonasal Hemostat, 15cm	
07-181-23MA	Kerrison Rongeur, 178mm 40 Degree Down, 2mm Biting		06-826-00	Mayo-Hegar Needle Holder, 16cm	
07-181-24MA	Kerrison Rongeur, 178mm 40 Degree Down, 3mm Biting		07-057-10	Weitlaner Retractor, 13cm	
07-181-25MA	Weil-Blakesley Rongeur, 13cm WL Straight 2mm		07-120-10	#3 Scalpel Handle	
07-181-26MA	Weil-Blakesley Rongeur, 13cm WL 30 Degree Up 2.6mm x 1.8mm		07-181-01	Adson Tissue Forceps, 12cm 1x2 Teeth	
07-181-27MA	Weil-Blakesley Rongeur, 2.8 x 2.5mm 13cm WL Straight		07-181-10	Dressing Forceps, 16cm	
07-181-28MA	Weil-Blakesley Rongeur, 13cm WL 30 Degree Up 3mm		09-061-10	Mayo Scissors, 17cm Straight	
02-201-01	Kerrison Rongeur, 18cm 40 Degree 1mm Up Biting		09-071-10	Metzenbaum Scissors, 18cm Straight	
02-201-02	Kerrison Rongeur, 18cm 40 Degrees 2mm Up Biting		09-071-11	Metzenbaum Scissors, 18cm Curved	
02-201-03	Kerrison Rongeur, 18cm 40 Degrees 3mm Up Biting		09-214-10	Allis Forceps, Straight 15cm	
02-034-10	Left Mucosal Dissector, Spear-Like		K-5310	Straight Blade 1mm (Box of 5)	
02-034-11	Right Mucosal Dissector, Spear-Like		K-5520	Hook Blade 30° 1mm (Box of 5)	
07-960-05L	Tapered Suction 5Fr Long		Product #	Description	LIFT OUT CURETTE TRAY
07-960-05XL	Tapered Suction 5Fr XL		02-034-12	Spatula Dissector, 2.4mm Curved Up	
07-960-06L	Tapered Suction 6Fr Long		02-034-13	Spatula Dissector, 2.4mm Curved Down	
07-960-06XL	Tapered Suction 6Fr XL		07-820-24	Probe Teardrop Tip, 2.5mm, Straight	
07-960-07L	Tapered Suction 7Fr Long		07-820-26	Probe Teardrop Tip, 2.5mm, 45 Degrees Up	
07-960-07XL	Tapered Suction 7Fr XL		07-820-27	Probe Teardrop Tip, 2.5mm, 90 Degrees Up	
07-960-08L	Tapered Suction 8Fr Long		07-820-29	Probe Teardrop Tip, 2.5mm, 60 Degrees Down	
07-960-08XL	Tapered Suction 8Fr XL		07-822-21	Ring Curette, 1.5mm, 45 Degrees Up	
07-960-09L	Tapered Suction 9Fr Long		07-822-23	Ring Curette, 3mm, 45 Degrees Up	
07-960-09XL	Tapered Suction 9Fr XL		07-822-25	Ring Curette, 5mm, 45 Degrees Up	
07-960-10L	Tapered Suction 10Fr Long		07-822-27	Ring Curette, 7mm, 45 Degrees Up	
07-960-10XL	Tapered Suction 10Fr XL		07-822-31	Ring Curette, 1.5mm, 90 Degrees Up	
			07-822-33	Ring Curette, 3mm, 90 Degrees Up	
			07-822-35	Ring Curette, 5mm, 90 Degrees Up	
			07-822-37	Ring Curette, 7mm, 90 Degrees Up	
			07-822-41	Ring Curette, 1.5mm, 45 Degrees Down	
			07-822-43	Ring Curette, 3mm, 45 Degrees Down	
			07-822-45	Ring Curette, 5mm, 45 Degrees Down	
			07-822-47	Ring Curette, 7mm, 45 Degrees Down	
			07-822-97	Ring Curette, 7mm, Left to Right Action	

Notes

This image shows a full page of blank, lined paper. It features approximately 20 evenly spaced horizontal gray lines across the entire width of the page, providing a guide for writing. The background is a solid off-white color. There are no margins, text, or other markings present.