



## Day-Bailes Tapered Locking Suction Tubes

Manufactured by:

**MIZUHO AMERICA**



**Mizuho America Inc.**  
30057 Ahern Avenue Union City, CA 94587  
www.mizuho.com

Distributed by:

**Mizuho America Inc.**  
30057 Ahern Avenue Union City, CA 94587  
www.mizuho.com



**Caution**

- Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
- Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.
- For after-sale support, contact your local sales representative.

# DAY-BAILES SUCTION TUBE

## Designed for a wide range of delicate spine and microsurgical procedures

- Ergonomic – slim, lightweight and easy on surgeons’ hands
- Cost Effective – interchangeable, quick release tube for replacement
- Gold suction tip – easier identification of the tip during surgery

SUCTION TUBE LENGHTS  
AVAILABLE FROM **SHORT TO EXTRA LONG**



ALSO AVAILABLE IN  
**ATRAUMATIC TIP**



### Handle Design

- Flat or Triangular hand piece
- Tear-drop design for suction control
-  indicates multi-position coupling, which allows the suction tube to fit handle at eight different positions



### Ordering Information

| Product #     | Description                                           |
|---------------|-------------------------------------------------------|
| 07-DAY-00     | Day-Bailes Suctions Organization Tray                 |
| 07-119-01SC   | Suction Tube Handle, Flat 8.5cm TL                    |
| 07-DAY-01SC   | Suction Tube Handle, Teardrop Control, 11cm TL        |
| 07-DAY-02SC   | Suction Tube Handle, Flat, Teardrop Control, 8.5cm TL |
| 07-DAY-03LSC  | Tapered Suction Tube, 3Fr, 11.3cm WL                  |
| 07-DAY-03MSC  | Tapered Suction Tube, 3Fr, 8.8cm WL                   |
| 07-DAY-03SSC  | Tapered Suction Tube, 3Fr, 6.8cm WL                   |
| 07-DAY-03XLSC | Tapered Suction Tube, 3Fr, 13.3cm WL                  |
| 07-DAY-04LSC  | Tapered Suction Tube, 4Fr, 11.3cm WL                  |
| 07-DAY-04MSC  | Tapered Suction Tube, 4Fr, 8.8cm WL                   |
| 07-DAY-04SSC  | Tapered Suction Tube, 4Fr, 6.8cm WL                   |
| 07-DAY-04XLSC | Tapered Suction Tube, 4Fr, 13.3cm WL                  |
| 07-DAY-05LSC  | Tapered Suction Tube, 5Fr, 11.3cm WL                  |
| 07-DAY-05MSC  | Tapered Suction Tube, 5Fr, 8.8cm WL                   |
| 07-DAY-05SSC  | Tapered Suction Tube, 5Fr, 6.8cm WL                   |
| 07-DAY-05XLSC | Tapered Suction Tube, 5Fr, 13.3cm WL                  |
| 07-DAY-06LSC  | Tapered Suction Tube, 6Fr, 11.3cm WL                  |
| 07-DAY-06MSC  | Tapered Suction Tube, 6Fr, 8.8cm WL                   |
| 07-DAY-06SSC  | Tapered Suction Tube, 6Fr, 6.8cm WL                   |
| 07-DAY-06XLSC | Tapered Suction Tube, 6Fr, 13.3cm WL                  |
| 07-DAY-07LSC  | Tapered Suction Tube, 7Fr, 11.3cm WL                  |
| 07-DAY-07MSC  | Tapered Suction Tube, 7Fr, 8.8cm WL                   |
| 07-DAY-07SSC  | Tapered Suction Tube, 7Fr, 6.8cm WL                   |
| 07-DAY-07XLSC | Tapered Suction Tube, 7Fr, 13.3cm WL                  |
| 07-DAY-08LSC  | Tapered Suction Tube, 8Fr, 11.3cm WL                  |
| 07-DAY-08MSC  | Tapered Suction Tube, 8Fr, 8.8cm WL                   |
| 07-DAY-08SSC  | Tapered Suction Tube, 8Fr, 6.8cm WL                   |
| 07-DAY-08XLSC | Tapered Suction Tube, 8Fr, 13.3cm WL                  |
| 07-DAY-09LSC  | Tapered Suction Tube, 9Fr, 11.3cm WL                  |
| 07-DAY-09MSC  | Tapered Suction Tube, 9Fr, 8.8cm WL                   |
| 07-DAY-09SSC  | Tapered Suction Tube, 9Fr, 6.8cm WL                   |
| 07-DAY-09XLSC | Tapered Suction Tube, 9Fr, 13.3cm WL                  |
| 07-DAY-10LSC  | Tapered Suction Tube, 10Fr, 11.3cm WL                 |
| 07-DAY-10MSC  | Tapered Suction Tube, 10Fr, 8.8cm WL                  |
| 07-DAY-10SSC  | Tapered Suction Tube, 10Fr, 6.8cm WL                  |
| 07-DAY-10XLSC | Tapered Suction Tube, 10Fr, 13.3cm WL                 |
| 07-DAY-11LSC  | Tapered Suction Tube, 11Fr, 11.3cm WL                 |
| 07-DAY-11MSC  | Tapered Suction Tube, 11Fr, 8.8cm WL                  |
| 07-DAY-11SSC  | Tapered Suction Tube, 11Fr, 6.8cm WL                  |
| 07-DAY-11XLSC | Tapered Suction Tube, 11Fr, 13.3cm WL                 |
| 07-DAY-12LSC  | Tapered Suction Tube, 12Fr, 11.3cm WL                 |
| 07-DAY-12MSC  | Tapered Suction Tube, 12Fr, 8.8cm WL                  |
| 07-DAY-12SSC  | Tapered Suction Tube, 12Fr, 6.8cm WL                  |
| 07-DAY-12XLSC | Tapered Suction Tube, 12Fr, 13.3cm WL                 |

| Product #   | Description                                          |
|-------------|------------------------------------------------------|
| 07-961-05L  | Multi-Hole Suction Tube, 5Fr, Atraumatic, 11.5cm WL  |
| 07-961-05XL | Multi-Hole Suction Tube, 5Fr, Atraumatic, 13.5cm WL  |
| 07-961-07L  | Multi-Hole Suction Tube, 7Fr, Atraumatic, 11.5cm WL  |
| 07-961-07XL | Multi-Hole Suction Tube, 7Fr, Atraumatic, 13.5cm WL  |
| 07-961-09L  | Multi-Hole Suction Tube, 9Fr, Atraumatic, 11.5cm WL  |
| 07-961-09XL | Multi-Hole Suction Tube, 9Fr, Atraumatic, 13.5cm WL  |
| 07-961-12L  | Multi-Hole Suction Tube, 12Fr, Atraumatic, 11.5cm WL |
| 07-961-12XL | Multi-Hole Suction Tube, 12Fr, Atraumatic, 13.5cm WL |

