




## HEART TISSUE STABILIZER

For beating heart surgery

- Strong reusable design.
- Superior positioning and articulation.
- Attachable to any sternal retractor system.
- Special interlock-joints for ultimate strength.
- Maximum versatility with interchangeable optional attachments for stabilization, manipulation or retraction.
- Ball-joint connection of attachments provides additional degrees of freedom.
- Partially detachable for ease of cleaning.



### Ordering Information

<b>04-5201</b>	Tissue stabilizer- arm, 24cm, working with flexible arm (04 - 9501) or with 7mm Shank diameter Retractor.
<b>04-7001</b>	Tissue stabilizer- arm, 24cm , working with ball head sucker.
<b>04-9501</b>	Retractor flexible arm 24cm. 7mm Shank diameter.

### OpInstruments GmbH

Tel: +49 69 9501 4824 | Email: [info@opinstru.de](mailto:info@opinstru.de) | [www.opinstru.de](http://www.opinstru.de)  
Add: Barckhausstr.11, 60325 Frankfurt am Main, Germany | Version: 2017-2.0A

Marketed In India By:-  
**INSTRUMENTED DESIGN & PRODUCTION COMPANY**  
Unit No-5, Navyug Industrial Estate,  
T J Road, Sweri (West)  
Mumbai-400015  
Tel-022 24188990  
Email: [idpc.ho@gmail.com](mailto:idpc.ho@gmail.com)

# SUCKERS

Suckers are manufactured by Recovre.

---



## 010101

- U-shaped sucker with four holes, single use and sterile, spherical joints, could connect with 04-9501 and 04-7001.



## 010201

- The head of apex sucker is stable and absorbs effectively the apex without myocardial injury..
- Apical-shaped sucker, single use and sterile, spherical joints, should be connect with 04-9501 and 04-7001.



## 010301

- U-shaped sucker with three holes, single use and sterile, spherical joints, could connect with 04-9501 and 04-7001.



## 010302

- U-shaped sucker with four holes, single use and sterile, spherical joints, could connect with 04-9501 and 04-7001.



## 010501

- U-shaped pressure plate, single use and sterile, spherical joints, could connect with 04-9501 and 04-7001.



## 010601

- O-shaped pressure plate, single use and sterile, spherical joints, could connect with 04-9501 and 04-7001.