

Needles Need EthyCalm

EthyCalm **+** PLUS



- ✓ Instant spray for fast venous access
- ✓ Reduces handling and stress for patients
- ✓ Saves you time and money

INVICTA
ANIMAL HEALTH

EthyCalm PLUS

Only £34 per vial... Provides over 100 applications

Indications

EthyCalm Plus is a vapo-coolant spray for topical use. Suitable for venepuncture, tick removal, microchip placement, venous catheterisation etc.

How To Use

Use immediately prior to intended procedure.
Hold the vial inverted over the skin at a distance of 20 cm.
Spray for 6 seconds.

Alternatively invert and spray EthyCalm Plus onto a swab for 8 seconds
and immediately wipe swab slowly over intended site three or four times.



Quantity Required:

Name:

Practice Name:

Address:

Postcode:

Email:

Wholesaler:

Please fax this order to **01403 790062**

Further information is available on request from:

Invicta, Unit 6 Shiprods Farm, Bashurst Hill, Itchingfield, Horsham, West Sussex RH13 0PD

T: 01403 791313 F: 01403 790062

E: mail@invictavet.com www.invictavet.com


ANIMAL HEALTH

ETH_BR05_06