

AED

- AUTOMATED EXTERNAL DEFIBRILLATOR -

For Training Contact:

888-407-8770

EmergencyTrainingPartners.com




**Emergency
Training
Partners, LLC**

ONLY USE WHEN UNRESPONSIVE AND NOT BREATHING



**♥ Make Sure 911 Has
Been Called**



♥ Turn On AED



♥ Follow Voice Prompts

- Expose Chest
- Apply Pads
- Do NOT Touch The Patient
- Shock (if indicated)
- Provide CPR