

Hair and Fingernail Testing - Alcohol & Other Drugs

Hair and Fingernail Drug Panels

	Panels:	19	18	17	16	15	14	12	10	9	7	5
Amphetamines amphetamine, MDA, MDEA, MDMA, methamphetamine*		●	●	●	●	●	●	●	●	●	●	●
Cannabinoids carboxy-Δ9-THC		●	●	●	●	●	●	●	●	●	●	●
Cocaine benzoylecgonine, cocaethylene, cocaine, norcocaine		●	●	●	●	●	●	●	●	●	●	●
Opiates 6-MAM, codeine, hydrocodone, hydromorphone, morphine, norhydrocodone		●	●	●	●	●	●	●	●	●	●	●
Phencyclidine phencyclidine (PCP)		●	●	●	●	●	●	●	●	●	●	●
Benzodiazepines alprazolam, diazepam, midazolam, nordiazepam, oxazepam, temazepam		●	●	●	●	●	●	●	●	●	●	●
Barbiturates amobarbital, butalbital, pentobarbital, phenobarbital, secobarbital		●	●	●	●	●	●	●	●	●	●	●
Methadone EDDP, methadone		●	●	●	●	●	●	●	●	●	●	●
Propoxyphene norpropoxyphene, propoxyphene		●	●	●	●	●	●	●	●	●	●	●
Oxycodone noroxycodone, oxycodone, oxymorphone		●	●	●	●	●	●	●	●	●	●	●
Meperidine normeperidine		●	●	●	●	●	●	●	●	●	●	●
Tramadol tramadol		●	●	●	●	●	●	●	●	●	●	●
Fentanyl acetyl fentanyl, acetyl norfentanyl, fentanyl, norfentanyl		●	●	●	●	●	●	●	●	●	●	●
Sufentanil alfentanil, norsufentanil, sufentanil		●	●	●	●	●	●	●	●	●	●	●
Ketamine ketamine, norketamine		●	●	●	●	●	●	●	●	●	●	●
Buprenorphine buprenorphine, norbuprenorphine		●	●	●	●	●	●	●	●	●	●	●
Zolpidem zolpidem (e.g. Ambien®)		●	●	●	●	●	●	●	●	●	●	●
Kratom 7-hydroxymitragynine, mitragynine		●	●	●	●	●	●	●	●	●	●	●
Gabapentin gabapentin		●	●	●	●	●	●	●	●	●	●	●

*D/L Methamphetamine Isomer testing is available upon request for specimens that test positive for methamphetamine. Turnaround time is 5-7 business days. Turnaround time is 2-4 business days for additional add-ons or for any testing of hair or fingernail 15-panel and above.